


Seating - Kinder Chairs

Clinical reports link a child's postural control with functional movements of the upper extremities.



K2

K1 (tray included)




The Kaye Kinder Chairs

Kinder Chairs provide many options to gain the best alignment and postural control possible so that a child can focus on developing and coordinating arm and hand skills. Clinical experts recommend a firm back and appropriate seat depth to align the pelvis and femurs.



K1 Kinder Chair with optional KFS Pelvic-Femoral Support



K2CB/KCFRE installed on K2

The Kaye Kinder Chairs have an adjustable seat depth that assures that the child's thighs are well supported to prevent undue internal rotation stress on the femoral heads. The firm back and adjustable seat depth helps maintain the pelvis in the vertical position with thoracic extension.

The adjustable footrest supports the feet when the chair is in its higher settings, such as when it is adjusted to fit under a standard table. This footrest is removable to allow taller children to rest their feet on the floor when the chair is in the lower settings.

The abductor pad at the distal third of the thighs maintains the alignment of the hips and knees and inhibits adduction and internal rotation. This, combined with the seat belt across the pelvis, decreases the tendency for the child to either push or slide out of the chair. Additional pelvic stability and femoral alignment can be obtained by adding the optional pelvic femoral stabilizer.



K3 with tray

Model Features	K1 (low/small)	K2 (high/small)	K3 (large)
Age	1-7 years	1-7 years	6-12 years
Size & Weight	Up to 48" / 60 lbs (27 kg)	Up to 48" / 60 lbs (27 kg)	Up to 56" / 100 lbs (45 kg)
Seat depth	Adjusts 7" to 12"	Adjusts 7" to 12"	Adjusts 11" to 16"
Seat width	11"	11"	14"
Back angle	Reclines 10° from vertical	Reclines 10° from vertical	Reclines 10° from vertical
Footrest	Adjusts 3" to 10.5" below seat	Adjusts 3" to 10.5" below seat	Adjusts 9" to 17.5" below seat
Height of arms	5.25" above seat	5.25" above seat	5.5" above seat
Chair height - floor to top of arm rests	Adjusts 14.75" to 22"	Adjusts 26" to 28.5"	Adjusts 24" to 29"
Tray dimensions	25"w x 19.75"d, included	Tray is optional. To order with tray, specify model K2+T	25"w x 23.5"d, included
Optional Caster Base *	Model K1CB, includes 4" swivel, locking casters	Model K2CB includes 4" swivel, locking casters	Standard swivel, locking casters included
Optional Footrest Extension **	Model KCFRE, increases distance from footrest to seat 7"	Model KCFRE, increases distance from footrest to seat 7"	N/A
Optional Pelvic-Femoral Stabilizer ***	Model KFS, 2 each 4" x 6" knee pads	Model KFS, 2 each 4" x 6" knee pads	N/A

* Caster Base: Model K1CB/K2CB - 4" locking, swivel casters; includes replacement legs to attach the caster base; K2CB adds 6" to chair height.

** Footrest Extension: Model KCFRE - Fits K1 and K2; adds 7" to the distance between the footrest and the seat.

*** Pelvic-Femoral Support: Model KFS - fits K1/K2; provides additional stability and alignment for the user's pelvis, legs, and lower body.