


## IMPORTANT SAFETY INSTRUCTIONS SAVE THESE INSTRUCTIONS



### AC Adapter for Omron® Digital Blood Pressure Monitor Model HEM-ADPT907

Thank you for purchasing the Omron® HEM-ADPT907 AC adapter.

-  Please read this instruction manual before using the AC adapter.  
Please keep this instruction manual for future reference.

#### -DANGER-

#### TO REDUCE THE RISK OF FIRE OR ELECTRIC SHOCK, CAREFULLY FOLLOW THESE INSTRUCTIONS.


-  For connection to a power supply not in the U.S.A, use an attachment plug adapter of the proper configuration for the power outlet.
-  The power unit (AC adapter) is intended to be correctly oriented in a vertical or floor mount position.

#### TO CONNECT THE AC ADAPTER



- Insert the AC adapter plug into the AC adapter jack on the monitor.
- Plug the AC adapter into a 120V AC outlet (60 cycles).

#### TO DISCONNECT THE AC ADAPTER

- Hold and pull the housing to remove the AC adapter from the electrical outlet. Do not remove by pulling the cord.
- Remove the AC adapter plug from the monitor.




#### WARNING

-  Do not plug or unplug the power cord into the electrical outlet with wet hands.
-  Do not overload power outlets. Plug the device into the appropriate voltage outlet.

#### NOTES

Operating temperature / humidity / air pressure:  
50 °F to 104 °F (10 °C to 40 °C) / 15 to 95% RH / 700 to 1060 hPa  
Storage temperature / humidity:  
-4 °F to 140 °F (-20 °C to 60 °C) / 10 to 95% RH

 : Class II equipment

 : For indoor use only

 : Direct current

 : Alternating current