

Aluminum Series Back Supports Loop Handle

STEP 1



With cushion removed, insert M6x10mm button head with one serrated washer through the inside of the backrest.

STEP 2



After the bolt is through the AL shell, place an additional serrated nut on bolt on the outside of the shell

STEP 3



With the cable fully inserted into the cable end piece, thread bolt into cable end.

Note: The M6x10mm bolt will crimp the cable inside the cable end piece, securing cable in place.

STEP 4



Repeat steps 1-3 for the other mounting hole.

Note: The cable should exit the cable ends on top.



Note: The black cable ends should be tight against the shell, as well as the cable crimped tightly inside the cable ends. Pull on the cable to ensure everything is secured properly.




PSIS Support Insert

When positioning an ADI back, the interior portion of the back support should be positioned to make contact at the PSIS to provide optimal pelvic support and trunk control. If additional support is required, ADI offers the support insert shown below as an optional feature.

To install, unzip the fabric cover and insert the PSIS support between the yellow foam surface and the fabric cover.

Position as desired.

Note:
For a secure hold, double sided Velcro is provided.



position support insert as desired

ZIP Privacy Flap

ADI privacy flaps are an optional feature. Our new design simply "Zips" onto the bottom of the fabric cover for a secure hold. This design allows you to add a privacy flap at the time of order or when needed.

STEP #1

attach privacy flap to the fabric back support cover with zipper



STEP #2

Velcro privacy flap to seatbase



Note:

a strip of Velcro hook & loop is sewn into flap to attach to all bases / slings



Fabric Cover

Privacy Flap Zipper

ZIP Privacy Flap

