



Support

Immedia Sling

MULTIFUNCTIONAL TRANSFER AID

The Immedia Sling is a multifunctional transfer aid that can be used in several situations.

Immedia Sling with handles acts as an extension of the arms, allowing the caregiver to offer support and assistance whilst maintaining a good ergonomic position and grip, without encroaching on the user's personal space.

Sling is a flexible tool that can be used to support many different transfers such as, moving higher up in bed, transfer between sitting positions, etc.

Immedia Sling has 4 handles on each side and a non-slip material on the inside to secure the user during the transfer. The outside has a low friction surface to aid transfer and removal.



Immedia Sling is a flexible aid offering a number of different functions.



Immedia Sling with handles acts as an extension of the arms, allowing the caregivers to offer support and assistance whilst maintaining a good ergonomic position and grip.

**RECOMMENDED
USE!**



1. The Immedia Sling provides excellent handles for caregivers. Used with Immedia MultiGlide or similar low friction devices, you ensure a safe working environment when users need to be repositioned in bed.




2. Even with one caregiver, it is possible to reposition the user higher in bed using the Immedia Sling as shown.



3. With one or two caregivers, a user can be repositioned further back in a chair by moving the user slightly from side to side.

TIP!

When two caregivers are working together one should always give a precise command when to push or pull.



INTENDED USE

- Support moving higher up in bed
- Support transfer between bed, wheelchair, chair, etc
- Support stand up / sit down
- Moving back / forward in chair / wheelchair

DESCRIPTION	WIDTH X LENGTH		ITEM NO.	SPECIFICATIONS
IMMEDIA SLING	CM	INCH		
~ Long, 2 x 4 handles	16 x 190	6.3 x 74.8	IM428	BACKSIDE: Nylon
~ Short, 2 x 4 handles	16 x 120	6.3 x 47.2	IM428K	PADDING: Polyester
				INSIDE: Nylon with polyurethane
				HANDLES: Polyester
WASHING INSTRUCTIONS				
Do not use fabric softener or bleach.				