

The Roosevelt

**“Quick-Change”
Incontinent Cover
Instructions**



Table of Contents

Warnings	3
Incontinent Cover Details	4
Installing the Incontinent Cover	6
Care and Cleaning	8
Fabric Care	8
Warranty	9

WARNING!
**FAILURE TO COMPLY WITH THE
FOLLOWING WARNINGS CAN RESULT
IN DEATH OR SERIOUS INJURY TO
YOUR CHILD:**

- The incontinent cover described in this instruction booklet is only to be used with the *Roosevelt* child safety seat manufactured by Merritt Manufacturing, Inc. **Do not use with any other child restraint system.**
- Prior to each use, verify that the installation of the incontinent cover is **not** creating any slack in the harness or buckle webbing.
- These are supplemental instructions which only describe how to use the incontinent cover. Please read and follow all instructions on this child restraint and in the written instructions located in the pocket on the side of the restraint.

Incontinent Cover Details

Top Side of Incontinent Cover

Upper Tab Set
Left

Upper Tab Set
Right



Lower Tabs

“Dam” designed to prevent spills from migrating to the main cover.

Incontinent Cover Details, Continued

Bottom Side of Incontinent Cover

Upper Tab Set Right
(Bottom side of cover)

Upper Tab Set Left
(Bottom side of cover)



Lower Tabs


Installing the Incontinent Cover

Note: Steps 1-12 may be performed without removing the restraint from the motor vehicle.

1. Disconnect the buckle and chest clip and lay each harness strap over the side of the restraint (Figure A).
2. Place the incontinent cover on the restraint as shown in Figure B.
3. Grasp the lower tabs (See arrows - Figure C) and insert them through their corresponding slots in the restraint's cover.
4. Route the lower tabs behind the cover (between the cover and the foam).

Warning: Do not route the Lower Tabs behind the plastic shell of the restraint.

5. Pull the tabs toward each other (underneath the cover) until the Velcro on each of them overlaps (Figure D).
6. Velcro the two tabs together underneath the cover (Figure D).



Installing the Incontinent Cover, Continued

7. Grasp one set of the upper tabs (Figure E) and insert one of them through the long vertical belt slot and the other one through the lowest harness slot located on the side of the restraint's cover nearest to the tab.
8. Route both of these tabs behind the foam (between the foam and the restraint's plastic shell).



Warning: Do not route the Upper Tabs behind the plastic shell of the restraint.

9. Pull the tabs toward each other until the Velcro on each of them overlaps (Figure F).
10. Velcro the two tabs together underneath the foam (Figure F).
11. Repeats Steps 7-10 for the other set of upper tabs.



Note: The fully installed incontinent cover is shown in Figure G.

Note: Proceed to Step 12 on the next page.



Installing the Incontinent Cover, Continued

12. Secure your child in the restraint according to the applicable steps on Pages 44 and 45 of the Roosevelt's "Instructions for Use" manual.



Note: To remove the incontinent cover, perform Steps 1-12 in reverse.

Care and Cleaning

Fabric Care

1. After removal, clean the incontinent cover separate from any other clothing per Steps 2-4.
2. Velcro all mating tabs together (See Pages 4 and 5) so that the Velcro does not fray the fabric during cleaning.
3. Machine wash on gentle cycle using mild soap (detergent) with the "wash" and "rinse" cycles set on "cold".
4. Machine dry in dryer on "no heat" or "fluff dry" cycle. Add a dry bath towel (an older one that is unlikely to "bleed" on the cover) to expedite the drying process.
5. Do not iron.

Limited Warranty

Merritt Manufacturing, Inc. (Merritt) provides the following limited warranty on the Roosevelt:

Merritt warrants to the original purchaser of the Roosevelt that, if properly used strictly in accordance with the instructions and specifications included in the original packaging, it shall be free from defects in materials or workmanship for a period of one (1) year from the date of original purchase. Proof of purchase is required, including, but not limited to furnishing to Merritt proof that the owner registration card was returned to Merritt within thirty days of the original purchase or furnishing to Merritt proof of purchase that identifies the claimant as the original purchaser. If any component covered by this warranty fails for reasons covered by this warranty, Merritt will, at its option, provide the original purchaser with repair or replacement components for the covered product or component. This limited warranty is available only for purchasers in the United States and Canada. Merritt reserves the right to discontinue or change fabrics, parts, models or products, or to make substitutions. Please note that the upholstery fabric is subject to some degree of wear as a result of normal use. This kind of normal wear is not covered under warranty.

To make a claim under this warranty, you must contact Merritt via telephone at (317) 409-0148 or write to the company using the address on the back cover of this booklet. PLEASE COMPLETE AND MAIL THE OWNER REGISTRATION CARD WITHIN THIRTY DAYS OF PURCHASE TO PRESERVE YOUR WARRANTY RIGHTS.

WARRANTY LIMITATIONS

This warranty does not include damages due to external sources, including, but not limited to accident, negligence, misuse, neglect, alteration, repair, improper installation, improper testing or use that is not in conformance with the instructions and specifications in the instruction booklet.

LIMITATION OF LIABILITY AND DAMAGES

THE WARRANTY AND REMEDIES SET FORTH ABOVE ARE THE SOLE AND EXCLUSIVE REMEDIES AND ARE PROVIDED IN LIEU OF ALL OTHERS, ORAL OR WRITTEN, EXPRESS OR IMPLIED. IN NO EVENT WILL MERRITT, OR THE RETAILER/DISTRIBUTOR SELLING THIS PRODUCT, BE LIABLE TO THE PURCHASER, OR ANY USER OF THIS RESTRAINT OR OTHER THIRD PARTY, FOR ANY DAMAGES; INCLUDING LOST PROFITS, LOST BUSINESS, INCIDENTAL, INDIRECT, SPECIAL, CONSEQUENTIAL DAMAGES OR PUNITIVE DAMAGES. IN NO EVENT SHALL MERRITT OR ITS AFFILIATED COMPANIES BE LIABLE OR RESPONSIBLE FOR ANY

Limited Warranty, Continued

AMOUNT OF DAMAGES ABOVE THE AGGREGATE DOLLAR AMOUNT PAID BY THE CUSTOMER FOR THE PURCHASE OF THE PRODUCT COVERED BY THIS WARRANTY.

LIMITATIONS OF WARRANTIES AND OTHER WARRANTY TERMS AND STATES' LAW

MERRITT FULLY DISCLAIMS AND MAKES NO IMPLIED WARRANTIES OF ANY KIND, INCLUDING, BUT NOT LIMITED TO IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR USE, AND MAKES NO EXPRESS WARRANTIES BEYOND THOSE STATED HERE. NOTWITHSTANDING THE FOREGOING, IF ANY IMPLIED WARRANTIES ARE DEEMED TO BE GIVEN BECAUSE THE APPLICABLE LAW PROHIBITS THEIR EXCLUSION, THEN ANY SUCH IMPLIED WARRANTIES, INCLUDING IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, SHALL BE LIMITED TO THE DURATION AND TERMS OF THE EXPRESS WRITTEN WARRANTY. SOME STATES DO NOT ALLOW THE LIMITATION OF IMPLIED WARRANTIES, LIMITATIONS ON HOW LONG AN IMPLIED WARRANTY LASTS OR THE EXCLUSION OR LIMITATION OF INCIDENTAL OR CONSEQUENTIAL DAMAGES. THEREFORE, IN THOSE STATES, THE LIMITATIONS LISTED ABOVE MAY NOT APPLY TO YOU. THIS WARRANTY AFFORDS YOU WITH SPECIFIC LEGAL RIGHTS AND YOU MAY HAVE ADDITIONAL RIGHTS, DEPENDING ON THE STATE IN WHICH YOU RESIDE.

BINDING ARBITRATION

Any dispute, controversy or claim (whether in contract, tort or otherwise, whether pre-existing, present or future, and including statutory, common law, intentional torts and equitable claims) against Merritt, its agents, employees, successors, assigns or affiliates (collectively for purposes of this warranty, "Merritt"), arising from or relating to this warranty and agreement, its interpretation, or the breach, termination or validity thereof, the relationships which result from this warranty and agreement (including, to the full extent permitted by applicable law, relationships with third parties who are not parties to this warranty and agreement), SHALL BE RESOLVED EXCLUSIVELY AND FINALLY BY BINDING ARBITRATION ADMINISTERED BY THE AMERICAN ARBITRATION ASSOCIATION under its code of procedure then in effect. Any award of the arbitrator shall be final and binding on each of the parties and may be entered as a judgment in any court of competent jurisdiction.