



# 1<sup>ST</sup> DEFENSE

Dynamic Hybrid Pressure  
Prevention Mattress\*

The "New" Foam Mattress  
Replacement Alternative



- The #1 issue with a foam prevention mattress is that it is static, not dynamic.
- The #2 issue with a foam prevention mattress is that it takes a "set" or "sags" in the middle section over time.
- The main reason why every resident is not on a Dynamic Pressure Prevention System is due to high cost.
- The primary complaint from residents on a dynamic air mattress is that it is uncomfortable.



**The 1st Defense Pressure Ulcer Prevention System solves each of these problems. The dynamic mode continuously redistributes pressure every 10 minutes, the air bladder chamber will not take a set over time, and the foam insert system provides added comfort.**



## FEATURES:

- Dynamic pressure redistribution mode and low air loss therapy.
- High density/resilient foam insert creates a comfortable pressure redistribution support surface.
- Foam insert slopes from 2" to 1" offering added pressure redistribution in the vulnerable heel area.
- Top cover insert pocket protects the foam from contamination.
- Air flow dial offers custom patient pressure levels.
- Dynamic pressure pump is preset to alternate at a 10-minute cycle time.
- Static mode provides a firm surface for egress and ingress from the mattress.
- Nylon top cover is waterproof, vapor permeable, low shear, and anti-microbial.
- Optional: Foam raised side rails available.

Weight Capacity: 350 lbs. | Warranty: 2 Year Non-Prorated | Item#: 82030 | Size: 36"x80"x7"

### AVERAGE INTERFACE PRESSURES, mmHg - SUBJECT: MALE / 5' 9" / 170 LBS.

Scapula Max	Scapula Min	Sacral Max	Sacral Min	Heel Max	Heel Min	Trochanter Max	Trochanter Min
22.7	17.8	17.2	13.7	18.3	11.6	33	20.8
delta	-4.9	delta	-3.5	delta	-6.7	delta	-12.2

An Xsensor pressure mapping system was employed for conducting this evaluation. The pump was set to the median position of the air flow range. Two positions were employed: back lying and ninety-degree side lying. A 5-minute acclimation period was observed prior to the actual measurement period. Once the initial acclimation time was over, a 10-minute measurement period commenced. Pressure scans were obtained every 2 seconds for the duration of the measurement period. A 4" by 4" area representing 64 individual sensors was used to isolate and average the pressure points. The maximum average and the minimum average values were obtained from each individual pressure point.

*\* Independently tested by Element Materials Technology.*