

# Mobile Standers

**Standing skills, social skills** and the freedom to experience the world upright. The Mobile Stander provides benefits like these and so many more. This mobile standing device gives new independence for wheelchair-bound clients. It helps strengthen critical motor skills so a client can progress toward independent standing. In addition, it enables easy, upright interaction with peers and activities – providing enrichment opportunities and just plain fun.

## Freedom to explore

With the removable large wheels in place, clients can self-propel in a standing position—free to explore their surroundings.

## Freedom to participate

Remove its large wheels, and the Mobile Stander easily rolls right up to a table or counter enabling participation in crafts and other stationary activities.

## Skillbuilding in action

With regular use, the Mobile Stander provides a way for clients to strengthen motor skills and work toward independent standing.

## Why is standing so important?

### Research indicates it can:

- Improve bone density
- Prolong stretch for tight muscles at hips, knees and ankles
- Improve pulmonary and digestive function
- Enable eye-level social interaction
- Create access to more activities and experiences
- Support self-confidence and cognitive growth



“Nate has gained great freedom by using the Mobile Stander for his upright mobility. Previously wheelchair dependent, he now moves freely in the classroom and hallway. His arms easily propel the stander, allowing him to interact with his peers during classroom activities—even his adaptive PE class.”

Stephanie Heiser, PT  
West Central Joint Services  
Indianapolis, IN

DESIGNED FOR USE WITH THE  
**MOVE**<sup>®</sup>  
PROGRAM



With the big wheels on, Sarah can “walk” with her teacher to the next class.



Supported in the Mobile Stander, Jordan participates in activities he enjoys.



James can set or release the wheel lock when he is ready.



The wide, accessible deck makes transfers easy.

# Have fun while developing new skills.

Because the Mobile Stander encourages vertical posturing with weight-bearing and balance challenges similar to those experienced with independent standing, it represents a step up therapeutically from the Prone Stander. This is especially true if the device is used in reverse. With the large wheels in place, the Mobile Stander provides self-propelled exploration and freedom. When the large wheels are removed, the stander makes accessing a table or countertop a cinch. With this range of possibilities, users can busily interact with their peers and their environments, not even realizing they're working on weight-shifting, balance and postural control at the same time – all skills needed for independent standing.

## Why use a Mobile Stander?

- **For social interaction:** Users are placed in an upright and mobile position, facilitating eye-level interaction with their peers and participation in class activities.
- **For independent mobility:** With the large wheels, the Mobile Stander can be self-propelled providing the user with opportunities to explore independently.
- **To improve weight-bearing:** With the use of the seat pad, the Mobile Stander can support users with partial to full weight-bearing abilities and those with lower extremity contractures to help improve musculoskeletal development.
- **To improve weight-shifting and postural control:** The self-propelled mobility offered by the Mobile Stander stimulates postural righting reactions and weight-shifting.
- **To improve trunk control:** The adjustable body support on the Mobile Stander can be lowered to promote independent trunk control. Additionally, the Mobile Stander can be used in the reverse position to further challenge trunk control.



With the big wheels off, Amanda gets close to the table. The seat pad supports her body weight as she learns to stand.



When Carolyn can actively extend her hips and knees to bear her own weight, she is placed in the reverse position. This allows free arm movement as Carolyn's standing ability develops further.



An additional seat pad and fleece-covered strap add extra support to the back.



The conversion kit turns your K110 Mobile Stander into a K130 model.



Sandals and wedges encourage standing by positioning the client's feet correctly. Wedges lift toes or heels.



Arm prompts help stabilize clients with poor head and shoulder control.



To alter the height of the communication tray or arm prompts, press the button on the clamp and adjust the post up or down.



The communication tray holds pictures of Sarah's activity choices under protective clear plastic.

# Accessories



**Arm prompts** provide fully adjustable support for users with poor head or shoulder control.

(Require mounting bracket.)



**Communication tray** is fully adjustable and holds small objects or activities.

(Requires mounting bracket.)

8½" x 11" inside dimensions  
10 lb. max. working load.



**Mounting bracket** secures either arm prompts or hand loops as well as a communication tray.



**Sandals and wedges** secure feet in the right position. Adjust sandals to rotate or slide, and add wedges for additional heel or toe elevation of up to 1¾" to align legs if one is shorter.

Small: 4" x 6", med: 4½" x 8½",  
large: 5⅛" x 10½"



**Large wheels** with brakes pop on and off with the push of a button. Wheel locks provide security, spoke guards protect fingers, and hand rims keep hands clean.



**Additional seat pads** are available in two sizes.

Small: 6" x 7"  
Large: 8" x 11"



The **conversion kit** converts the smallest stander (K110) to the next size up (K130) for less than one quarter the cost of a new stander.



Secure and highly adjustable, **hand loops** boost a client's sense of stability and confidence while moving. They attach as a pair or singly on the mounting bracket.

(Require mounting bracket.)



A pair of hand loops attached to the mounting bracket of a K110 Mobile Stander.

# Features of the Mobile Stander

**Available in three colors:**  

 raspberry, blue and lime.

**Healthcare grade vinyl** covers padded body support.

**Body support adjusts** according to user's size and requirements.

**Soft fleecy strap** provides back or neck support.

**Our unique** adjustable seat pad allows gradual increases in weight-bearing and develops standing skills.

**We've made** our deck extra wide for simpler transfers and more user participation in the sit-to-stand process.

**Large casters** allow the stander to move easily during activities, yet lock securely during transfer.

**The fully adjustable** communication tray holds small objects or activities. (accessory)

**With our Mobile Stander**, you can add arm prompts to help clients develop greater head and shoulder control. Simply remove them as control improves. (accessory)

**Large wheels** pop on and off with the push of a button. Wheel locks provide security, spoke guards protect fingers and hand rims keep hands clean. (accessory)

**Leg straps** provide abduction when desired. (Extras are available on request.)



Mobile Stander with accessories (K150)



K170



K150



K130



K110

All four sizes are available in three colors.