

Supine Stander

E420 & E430 Product Manual



E430 Large Supine stander



E420 Small Supine stander




Contents


Important information	3
Safety messages	4
Recommended use and item dimensions	5
Check your order	6
Basic components	
Headboard	6
Footboard	7-8
Small Supine Stander main board (E420)	9-10
Large Supine Stander main board (E430)	11-12
Accessories	
Tray	13
Hand anchors	14
Arm anchors	15
Abduction wedge and block	16
Operation	17
Maintenance, cleaning and warranty	18
Materials and user modifications	19

Key for users

Use this key to determine which sections of this product manual apply to you.

 **Technical Users** For professionals who order and set up Rifton products.

 **Home Users** For care-givers who use Rifton products on a regular basis.

 **Maintenance Personnel** For anyone who is responsible for service or re-ordering of Rifton products and parts.

WARNING



- Thoroughly read and understand the information in this product manual before attempting to use this product. If the procedures and instructions in this product manual are not followed, serious injury or death could occur.
- A qualified professional must assess the appropriateness and safety of all equipment for each user.
- This product is intended for use by clients of unreliable judgment. Adult supervision is required at all times.
- To prevent falls and injuries:
 - Do not use this product on rough and uneven terrain, around swimming pools, or near stairways.
 - Ensure the appropriate use of straps and supports at all times. Straps and supports are provided for the safety of the user and must be carefully adjusted for comfort and security.
 - Tighten all adjustment knobs before use and immediately after making any adjustments.
- Do not use this product for clients outside the height and weight limits specified in this manual.
- To prevent structural failure, which may result in serious injury or death:
 - Inspect this product and accessories regularly for loose or missing screws, metal fatigue, cracks, broken welds, missing attachments, general instability or other signs of excessive wear.
 - Immediately remove this product from use when any condition develops that might make operation unsafe.
 - Do not use Rifton components or products for any purpose other than their intended use.

Recommended use

⚠️ WARNING To prevent falls and injury, do not use the Supine Stander as a transportation device.


The Rifton Supine Stander is a Class 1 medical device. It supports the user in the standing position. The Supine Stander differs from a tilt table in that it allows much greater control of the head, trunk, pelvis, knees and feet while allowing variable weight-bearing. It is especially suitable for users who lack sufficient anti-gravity control of the head and shoulder girdle.

For users who have sufficient head and upper trunk control, the headboard folds down to encourage independence.

The E430 large Supine Stander is recommended for children/adults between 46"– 72" (117 cm – 183 cm) tall and up to 250 lbs (114 kg).

The E420 small Supine Stander is recommended for children/adults between 30"– 50" (75 cm – 125 cm) tall and up to 100 lbs (45 kg).

User and item dimensions

User dimensions – inches (cm)	small	large
Height	30–50 (76-127)	46–72 (117-183)
 <p>Key user dimension: height Select the appropriate Stander by the user's overall height. Choose the model that allows for growth.</p> <p>Important: User's weight must not exceed the maximum working load.</p>		

Item dimensions – inches (cm)	small	large
Board length and width	51 x 16½ (130 x 42)	71 x 19 (180 x 48)
Distance between lateral supports	7–13 (18–33)	9–16 (23–41)
Height when horizontal	28½ (72)	21–30 (53–76)
Height when vertical	54½ (138)	74 (188)
Base length and width	37½ x 27 (95 x 69)	46½ x 29½ (118 x 75)
Max. working load in lbs. (kg)	100 (45)	250 (114)

Check your order

This product comes in a single carton, containing:

1. Supine Stander (includes padding, footboard, adjustable armrests, hip strap, chest strap, leg straps, one pair of lateral blocks, one pair of headblocks and double-locking casters).
2. Accessories for the Supine Stander.

You may not have ordered all of the available accessories, but use the diagrams that follow to confirm that your order is complete.

Basic components

Headboard (see Figures 6a and 6b)

⚠️ WARNING To prevent spinal injury, ensure the headboard is in place before tilting the client to a horizontal position.

To prevent pinching or crushing, ensure that all body parts are clear of the joint between headboard and main board during adjustment.

Attaching headboard:

Lift the headboard, slide the cleat (A) into frame slot, and tighten the cleat knob (B) (see Figure 7a).

Attaching headblocks:

1. Loosen two knobs (C) on the headblock adjusting clamp.
2. Remove the lateral block with curved bracket (D), turn it upright, and insert bracket into clamp.
3. Slide curved bracket into clamp past the safety dimples to desired spacing. Brackets must overlap if spacing of less than 8" is required.
4. Tighten both knobs (C) to secure headblock.
5. Repeat this procedure on the second headblock.




Figure 6a



Figure 6b. Headboard with headblocks attached.

If you are using the Supine Stander for a client whose head does not extend above the main board, remove the headblocks and the two clamps. Remount the clamps in the top slots of the main board. Attach the headblocks as described on page 6.

Adjusting headboard: (see Figure 6a)

Position headboard by raising it and sliding cleat (A) into the frame slot. Tighten cleat knob (B). If headboard is not needed, slide cleat up, tighten knob, and fold the headboard down.

Adjusting headblocks:

Loosen the two knobs on the clamp, slide the block to the desired position, and tighten knobs securely.

Footboard

⚠ WARNING To prevent falls, ensure the footboard is in place before tilting the client to a vertical position.

To prevent pinching or crushing, correctly position the footboard at the appropriate height before placing the client on the footboard.

To prevent injury, lock all four casters at all times, regardless of whether or not the product is occupied, except when moving the stander from one location to another.

Attaching sandals: (see Figure 7a)

Either large or medium sandals can be attached to the footboard.

1. Insert carriage bolts with bushings through the bottom of the footboard.
2. Insert sandal base onto bolts and adjust sandal before tightening knobs (A) securely.
3. Place sandal onto base and lock into place with latch (B).




Figure 7a

Adjusting sandals:

To adjust sandals, loosen knobs (A), move sandals to the desired position, and tighten knobs. Sandal height can be built up or a toe/heel angle achieved by using one or more wedges (see Figures 8a and 8b).

Attaching sandal wedges:

To add a wedge, unlatch and remove sandal from the base, place wedge over the base and latch it into place; any number of wedges can be added. Place and latch the sandal on top of the wedges.

Adjusting footboard:

The footboard can be adjusted to accommodate users of varying heights. Loosen the clamp knobs (C), move board to the desired height, and tighten knobs securely. For the user's security, footboard is best left as close as possible to the floor (see Figure 8b).

Adjusting the casters:

Swivel casters allow the Supine Stander to be moved easily. The foot-operated brake, locks both the wheel and the swivel simultaneously.




Figure 8a




Figure 8b

Small Supine Stander main board (E420)

Adjusting the main board:

The hand crank (A) controls the transfer height adjustment and the tilt angle of the main board. By turning the crank, the board will begin to tilt. Continue to turn the crank until the desired angle is reached (see Figure 9a).

Tilt angle is indicated on the board angle label. This can be used to record the correct angle setting for each user.



Figure 9a

Attaching lateral blocks:

Lateral blocks (B) ship unassembled. Follow instructions below to attach them.

1. Remove knob and plastic washer, leaving rubber washer on the lateral block.
2. Slide bolt of lateral block through slot in board.
3. Replace knob (C) and plastic washer on bolt on back of board, and tighten knob (C) securely (see Figure 9b).




Figure 9b

Adjusting lateral blocks and straps:

Any lateral block or strap can be adjusted in and out or up and down. To adjust in or out, loosen the appropriate knobs, slide the block or strap to the desired position, and tighten the knobs securely.

To adjust up and down or to completely remove a lateral block or strap, loosen the knob(s) completely and remove the block or strap. The block or strap can then be moved to a different set of slots, and the knobs tightened securely.

Adjusting hip and knee straps:

⚠️ WARNING To prevent falls and injury, secure the strap safety closures (C) when in use.

To adjust height of hip and knee straps, pull back of strap away from hook-and-loop fabric at back of the main board, slide strap to the desired height, and resecure hook-and-loop fabric (see Figure 10a).

Adjusting arm rests:

⚠️ WARNING To prevent pinching or crushing, ensure body parts are clear of potential crush points while adjusting armrest.

To adjust height of armrest, loosen both knobs (D), slide armrest to desired position, and tighten knobs securely. To adjust angle of armrest, loosen the bottom knob only, lift or lower front of armrest to desired position, and tighten knob securely (see Figure 10a).




Figure 10a

Large Supine Stander main board (E430)

Adjusting the main board:

The hand crank (A) controls transfer height adjustment and the tilt angle of the main board. By turning the crank, the board in its horizontal position will move vertically from 21" to 30" (53cm to 76cm) and then begin to tilt. Continue to turn hand crank until desired angle is reached (see Figure 11a).



Figure 11a

Tilt angle is indicated on the board angle label. This can be used to record the correct angle setting for each user.

Attaching lateral blocks:

Lateral blocks (B) ship in reverse position.

To reposition lateral blocks:

- Loosen knobs (C) and slide blocks out from under the clamp.
- Rotate blocks into correct position and replace them under the clamp.
- Retighten knobs (see Figure 11c).

Use knobs (C) to adjust the blocks laterally and vertically.




Figure 11b

Adjusting blocks and straps:

Any lateral block or strap can be adjusted in and out or up and down. Loosen the two knobs on the appropriate clamp, slide block or strap to the desired position, and tighten knobs securely.

To completely remove a lateral block or strap, loosen clamp knobs enough to slide bracket out past the safety dimples (D). The block or strap can then be moved to a different set of clamps and the knobs tightened securely (see Figure 11c).



Figure 11c

Adjusting hip strap:

⚠️ WARNING To prevent falls and injury, secure the strap safety closures (C) when in use (see Figure 12a).

It is possible to place the hip strap (A) in line with the lateral blocks. Loosen the clamp knobs enough to allow the hip strap metal bracket to lie on top of the lateral block. Tighten clamp knobs securely.

Adjusting knee straps:

To adjust height of knee strap, pull back of strap away from hook-and-loop fabric at back of the main board, slide strap to the desired height, and resecure hook-and-loop fabric.

Adjusting arm rests:

⚠️ WARNING To prevent pinching or crushing, ensure body parts are clear of potential crush points while adjusting armrest.

To adjust height of armrest, loosen both knobs (B), slide armrest to desired position, and tighten knobs securely. To adjust angle of armrest, loosen the bottom knob only, lift or lower front of armrest to desired position, and tighten knob securely (see Figure 12a).




Figure 12a

Accessories

Tray

Attaching the tray:

A clear or opaque tray may be fitted to the armrests. Position both armrests at the same height and angle. Attach tray by inserting the tray latches into appropriate armrest holes (see Figure 13a and 13b).

Adjusting the tray:

The tray can be adjusted in and out by opening the tray latches and moving them to a different set of holes in the armrests (see Figure 13b).



Figure 13a



Figure 13b

Hand anchors

⚠ WARNING To prevent falls and resulting injury, do not use Rifton Anchors as a primary support or as a standing support. Suction cups may release.

To prevent falls and cuts, do not use Rifton Anchors on thin glass.

Recommended use

Hand Anchors can be used on any of our available trays. They give the client an anchor to control excessive movement of the arms and encourage stabilization of the shoulder girdle. The client may hold on with both hands for stabilization training or with one hand, keeping the other hand free.

Operation

On a clean smooth table or tray, press down on suction cup and turn handhold clockwise until secure. To remove, turn handhold counterclockwise.

Tip: Surface must be very smooth and completely clean to get a strong seal. If suction cup does not hold well, wipe with a damp cloth or apply a thin layer of petroleum jelly around the rim of the rubber pad to enhance the vacuum.



Figure 14a

Figure 20d

Arm Anchors

⚠️ WARNING

To prevent falls and resulting injury, do not use Rifton Anchors as a primary support or as a standing support. Suction cups may release.

To prevent falls and cuts, do not use Rifton Anchors on thin glass.

Recommended use

The arm Anchor is intended to stabilize one arm on a table or tray, or as a pair to stabilize both arms. Providing increased stability to the forearm, arm Anchors may be beneficial for head and trunk control and to allow functional use of an unanchored arm.

Adjustment

Arm and wrist straps can be adjusted and secured with hook-and-loop fasteners. Move handhold backward and forward or rotate by loosening the knob.

Operation

On a clean smooth table or tray, press down on suction cup, and lower lever. Place arm on arm anchor pad and secure it with arm and wrist straps. If suction cup does not hold well, wipe with a damp cloth or apply a thin layer of petroleum jelly around the rim of the rubber pad to enhance the vacuum.



Figure 15a

Adjustable abduction wedge

(see Figure 16a)

Attaching

Remove the knobs (A) and one washer off each threaded stud.

1. Place abduction wedge in appropriate slot with one plastic washer on each threaded stud between wedge and Stander.
2. Put the second washer on each of the threaded studs between knobs and Stander.
3. Thread and tighten knobs securely.

Adjusting

To adjust the amount of abduction, raise one wing of the wedge to expose the adjustment screw (B). Turn screw out to increase abduction. Raise the second wing, and adjust it in the same way.




Figure 16a



Figure 16b

Round abduction block

(see Figure 16b)

Attaching

To attach the round abduction block, remove the knob and place where appropriate. Re-tighten knobs securely.

Collar

(see Figure 16c)

Attaching

To attach the collar to the abduction block, wrap it around the block, and fasten snaps.



Figure 16c

Operation

1. To prevent falls and injury, lock all four casters.
2. Crank the Supine Stander to desired horizontal height for ease of transfer.
3. Slide the armrest and the lateral block away from the transfer area. Open the chest, hip and knee straps.
4. Make sure that both the headboard and footboard are in place before transferring the user.
5. Transfer the user onto the Supine Stander in a horizontal position, and make final adjustments to all straps, ensuring that the safety closures on the straps are secure.
6. Adjust the headblock and lateral blocks.
7. Secure the user's feet in the sandals and adjust positioning.
8. Adjust the armrests to the user's elbow height, and set to the appropriate angle.
9. Make any final adjustments, ensuring that everything is in place, and that all knobs are tight and secure.
10. Crank the Supine Stander to the desired angle as indicated on the board angle label.
11. Attach tray on armrests if desired.

Maintenance 🗑️

This product is designed and tested for an expected life of 5 years when used and maintained in accordance with this manual. At all times, users must ensure that the product remains in a safe and useable condition, including regular maintenance and inspections as specified in the manual.

To prevent structural failure, which may result in serious injury or death:

- Inspect this product and accessories regularly for loose or missing screws, metal fatigue, cracks, broken welds, missing attachments, general instability or other signs of excessive wear.
- Immediately remove this product from use when any condition develops that might make operation unsafe.
- Do not use Rifton components or products for any purpose other than their intended use.
- Replace or repair components or products that are damaged or appear to be unstable.
- Use only Rifton authorized replacement parts. Order information for replacement parts is provided on the back of this product manual.

Twice a year lubricate all four casters. Apply a drop of light oil to the ball bearings and to the axle on each side of the caster wheel.

Once a year, check the lubrication of the threaded rod. If it seems dry, apply a general purpose grease to the length of the rod.

Cleaning 🧑🏠 🗑️

Clean with disinfectant wipes or a solution of up to 10% bleach. Do not use excessive water on pads and wooden surfaces.

Clean tray with mild soap or detergent. Avoid abrasive cleaners and aromatic spirits. Disinfectants can be used, but no solvents. Never scrape with razor blades or other sharp instruments. Do not clean in hot sun or at elevated temperatures. Do not send through a dishwasher.

The components of the arm Anchor can be removed for cleaning. Wipe clean with disinfectant wipes or a solution of up to 10% bleach.

Warranty Statement 🧑🏠 🗑️


If a Rifton product breaks or fails in service during the first year, we will replace it free of charge.

Materials

- Steel hardware items (nuts, bolts, screws, etc) are typically zinc or nickel plated, or stainless steel.
- Upholstery items (pads, support blocks, padded prompts, etc) are typically polyurethane foam with a fire-retardant cover made from expanded vinyl.
- Frames are typically steel or aluminum tubing, welded together, and coated with a baked-on paint finish. Some frame components may also be stainless steel.
- Straps are typically made of polypropylene or nylon webbing.
- Wooden components are typically birch plywood, solid maple, or laminated hardwood veneers, finished with a clear polyurethane lacquer.
- Plastic components are typically injection molded from a variety of industrial resins.

All materials are latex, lead and phthalates free.

User modifications

 WARNING To prevent serious injury or death, do not modify or alter Rifton products or components, or use Rifton products or components in conjunction with products from other manufacturers. Rifton does not accept responsibility for any modifications or alterations made to our components or products after they leave our premises. Customers modifying or altering our components or products, or using them in conjunction with products from other manufacturers, do so at their own risk.