

*AMSure*<sup>®</sup>  
UROLOGICAL CARE



Comprehensive and Cost-Effective Urological  
Solutions Engineered to Achieve Superior Results

# TABLE OF CONTENTS

AMSure® Urological Care product portfolio provides solutions to help improve patient safety and infection prevention. Our product selection includes:

- Self-Catheterization Kits ..... 1
- Urethral Catheters ..... 2
- Urethral Trays and Kits ..... 3
- Foley Catheters ..... 4
- Foley Insertion Trays ..... 5
- Syringes ..... 7
- Irrigation Trays and Kits ..... 7
- Urinary Collection Devices ..... 8
- Specimen Containers ..... 9



HCPCS Disclaimer:  
HCPCS codes and Home Health Consolidated Billing codes provided by Amsino are intended as general guidelines only. Amsino does not guarantee coverage or reimbursement of any products. You must address all coverage and reimbursement issues (including the correctness and accuracy of codes) with your individual payers. It is your responsibility to ensure the accuracy and appropriateness of each claim you submit, in accordance with all applicable payer requirements.

# SELF-CATHETERIZATION KITS

## URETHRAL SELF-CATHETERIZATION KIT WITH R-POLISHED® EYES

AMSure® R-Polished technology uses radio frequency to create ultra-smooth eyelets for optimal patient comfort. The AMSure Self-Catheterization kit allows the patient, even with limited dexterity, to perform incontinence care with confidence.

- Advanced R-Polished eyelets enhance patient comfort during catheter insertion and removal
- Introducer tip reduces the potential for urinary tract infection (UTI)
- Collection bag with built-in thumb holes allows for ease of use
- Sterile, latex-free and BPA-free

Kits include:

- Pre-lubricated 16 in. (40.64 cm) long PVC urethral catheter
- 1100 mL collection bag with insertion tip and protector cap
- (3) Povidone-iodine swab sticks and benzalkonium chloride prep pad
- Vinyl powder-free gloves
- Waterproof underpad 13 in. x 18 in. (33.02 cm x 45.72 cm)
- Waste disposal bag
- Gauze pad 4 in. x 4 in. (10.16 cm x 10.16 cm)

HPCS code A4353

| ITEM CODE | FRENCH SIZE | QTY/UNIT |
|-----------|-------------|----------|
| AS85008   | 8 Fr        | 50/Case  |
| AS85010   | 10 Fr       | 50/Case  |
| AS85012   | 12 Fr       | 50/Case  |
| AS85014   | 14 Fr       | 50/Case  |
| AS85016   | 16 Fr       | 50/Case  |



# URETHRAL CATHETERS

## RED RUBBER LATEX URETHRAL CATHETER

- Smooth tip and quality rubber latex promote patient comfort during insertion
- Sterile, DEHP-free and BPA-free



HPCS code A4351

| ITEM CODE      | DESCRIPTION              | FEATURES | QTY/UNIT |
|----------------|--------------------------|----------|----------|
| <b>AS44008</b> | 8 Fr, 11 in. (27.94 cm)  | Female   | 100/Case |
| <b>AS44010</b> | 10 Fr, 11 in. (27.94 cm) | Female   | 100/Case |
| <b>AS44012</b> | 12 Fr, 16 in. (40.64 cm) | Male     | 100/Case |
| <b>AS44014</b> | 14 Fr, 16 in. (40.64 cm) | Male     | 100/Case |
| <b>AS44016</b> | 16 Fr, 16 in. (40.64 cm) | Male     | 100/Case |
| <b>AS44018</b> | 18 Fr, 16 in. (40.64 cm) | Male     | 100/Case |
| <b>AS44020</b> | 20 Fr, 16 in. (40.64 cm) | Male     | 100/Case |
| <b>AS44022</b> | 22 Fr, 16 in. (40.64 cm) | Male     | 100/Case |
| <b>AS44024</b> | 24 Fr, 16 in. (40.64 cm) | Male     | 100/Case |
| <b>AS44028</b> | 28 Fr, 16 in. (40.64 cm) | Male     | 100/Case |

## PVC INTERMITTENT URETHRAL CATHETERS

- Staggered eyes help maintain tip strength during insertion and promote efficient drainage
- Sterile, latex-free and BPA-free



HPCS code A4351

| ITEM CODE       | DESCRIPTION              | FEATURES  | QTY/UNIT |
|-----------------|--------------------------|-----------|----------|
| <b>AS860608</b> | 8 Fr, 6 in. (15.24 cm)   | Female    | 50/Case  |
| <b>AS860610</b> | 10 Fr, 6 in. (15.24 cm)  | Female    | 50/Case  |
| <b>AS860612</b> | 12 Fr, 6 in. (15.24 cm)  | Female    | 50/Case  |
| <b>AS860614</b> | 14 Fr, 6 in. (15.24 cm)  | Female    | 50/Case  |
| <b>AS860616</b> | 16 Fr, 6 in. (15.24 cm)  | Female    | 50/Case  |
| <b>AS861006</b> | 6 Fr, 10 in. (25.40 cm)  | Pediatric | 50/Case  |
| <b>AS861008</b> | 8 Fr, 10 in. (25.40 cm)  | Pediatric | 50/Case  |
| <b>AS861010</b> | 10 Fr, 10 in. (25.40 cm) | Pediatric | 50/Case  |
| <b>AS861608</b> | 8 Fr, 16 in. (40.64 cm)  | Male      | 50/Case  |
| <b>AS861610</b> | 10 Fr, 16 in. (40.64 cm) | Male      | 50/Case  |
| <b>AS861612</b> | 12 Fr, 16 in. (40.64 cm) | Male      | 50/Case  |
| <b>AS861614</b> | 14 Fr, 16 in. (40.64 cm) | Male      | 50/Case  |
| <b>AS861616</b> | 16 Fr, 16 in. (40.64 cm) | Male      | 50/Case  |
| <b>AS861618</b> | 18 Fr, 16 in. (40.64 cm) | Male      | 50/Case  |



## PVC INTERMITTENT URETHRAL CATHETER WITH R-POLISHED EYES

- Advanced R-Polished technology uses radio frequency to create ultra-smooth eyelets for optimal patient comfort
- Staggered eyes help maintain tip strength during insertion and promote efficient drainage
- Designed to minimize trauma and reduce risk of infection
- Sterile, latex-free and BPA-free

HPCS code A4351



| ITEM CODE       | DESCRIPTION              | FEATURES  | QTY/UNIT |
|-----------------|--------------------------|-----------|----------|
| <b>AS960610</b> | 10 Fr, 6 in. (15.24 cm)  | Female    | 50/Case  |
| <b>AS960612</b> | 12 Fr, 6 in. (15.24 cm)  | Female    | 50/Case  |
| <b>AS960614</b> | 14 Fr, 6 in. (15.24 cm)  | Female    | 50/Case  |
| <b>AS960616</b> | 16 Fr, 6 in. (15.24 cm)  | Female    | 50/Case  |
| <b>AS961010</b> | 10 Fr, 10 in. (25.40 cm) | Pediatric | 50/Case  |
| <b>AS961610</b> | 10 Fr, 16 in. (40.64 cm) | Male      | 50/Case  |
| <b>AS961612</b> | 12 Fr, 16 in. (40.64 cm) | Male      | 50/Case  |
| <b>AS961614</b> | 14 Fr, 16 in. (40.64 cm) | Male      | 50/Case  |
| <b>AS961616</b> | 16 Fr, 16 in. (40.64 cm) | Male      | 50/Case  |

# URETHRAL INSERTION TRAYS AND KITS

## URETHRAL INSERTION TRAY

- Trays provide clinician flexibility to choose their preferred drainage or leg bag

### PVC CATHETER

- Staggered eyes help maintain tip strength during insertion and promote efficient drainage
- Sterile, latex-free and BPA-free

### RED RUBBER CATHETER

- Smooth tip and quality rubber latex promote patient comfort during insertion
- Sterile, DEHP-free and BPA-free

Trays include:

- 1000 mL outer basin tray
- 16 in. urethral catheter
- Waterproof and fenestrated drapes
- Vinyl powder-free gloves
- (3) 10% povidone-iodine swab sticks
- Lubricating jelly
- 4 oz. specimen container

HPCS code A4353



| ITEM CODE      | URETHRAL CATHETER TYPE | QTY/UNIT |
|----------------|------------------------|----------|
| <b>AS870</b>   | 14 Fr PVC (Latex-free) | 20/Case  |
| <b>AS87114</b> | 14 Fr Red Rubber Latex | 20/Case  |
| <b>AS87118</b> | 18 Fr Red Rubber Latex | 20/Case  |



## PRE-CONNECTED DRAINAGE BAG URETHRAL CATHETER KIT

- Sterile kits provide clinician flexibility to choose their preferred drainage or leg bag
- Gloves and povidone-iodine swab sticks help maintain sterility and reduce the risk of infection

Kits include:

- 1200 mL low-profile urinary drainage bag
- 16 in. (40.64 cm) urethral catheter
- CSR wrap, waterproof and fenestrated drapes
- (2) Nitrile powder-free gloves
- (3) Povidone-iodine swab sticks
- Lubricating jelly
- 2 oz. specimen container with ID label

HPCS code A4353



| ITEM CODE        | 14 FR RED RUBBER CATHETER | 14 FR PVC CATHETER | QTY/UNIT |
|------------------|---------------------------|--------------------|----------|
| <b>AS89214R*</b> | •                         |                    | 20/Case  |
| <b>AS89214</b>   |                           | •                  | 20/Case  |

\*Contains latex

# FOLEY CATHETERS

## 2-WAY SILICONIZED LATEX FOLEY CATHETERS

- Silicone coating helps reduce urethral irritation
- Smooth beveled symmetrical eyes minimize tissue damage and aid in maximum drainage
- Firm catheter tip designed to improve ease of insertion
- Symmetrical balloon provides a tight seal at the neck of the urethra
- Sterile, DEHP-free and BPA-free



HPCS code A4338



| ITEM CODE | DESCRIPTION              | FEATURES | QTY/UNIT |
|-----------|--------------------------|----------|----------|
| AS41012   | 12 Fr, 16 in. (40.64 cm) | 5 cc     | 10/Box   |
| AS41014   | 14 Fr, 16 in. (40.64 cm) | 5 cc     | 10/Box   |
| AS41016   | 16 Fr, 16 in. (40.64 cm) | 5 cc     | 10/Box   |
| AS41018   | 18 Fr, 16 in. (40.64 cm) | 5 cc     | 10/Box   |
| AS41020   | 20 Fr, 16 in. (40.64 cm) | 5 cc     | 10/Box   |
| AS41022   | 22 Fr, 16 in. (40.64 cm) | 5 cc     | 10/Box   |
| AS41024   | 24 Fr, 16 in. (40.64 cm) | 5 cc     | 10/Box   |
| AS41026   | 26 Fr, 16 in. (40.64 cm) | 5 cc     | 10/Box   |
| AS41028   | 28 Fr, 16 in. (40.64 cm) | 5 cc     | 10/Box   |
| AS41030   | 30 Fr, 16 in. (40.64 cm) | 5 cc     | 10/Box   |
| AS42012   | 12 Fr, 16 in. (40.64 cm) | 30 cc    | 10/Box   |
| AS42014   | 14 Fr, 16 in. (40.64 cm) | 30 cc    | 10/Box   |
| AS42016   | 16 Fr, 16 in. (40.64 cm) | 30 cc    | 10/Box   |
| AS42018   | 18 Fr, 16 in. (40.64 cm) | 30 cc    | 10/Box   |
| AS42020   | 20 Fr, 16 in. (40.64 cm) | 30 cc    | 10/Box   |
| AS42022   | 22 Fr, 16 in. (40.64 cm) | 30 cc    | 10/Box   |
| AS42024   | 24 Fr, 16 in. (40.64 cm) | 30 cc    | 10/Box   |
| AS42026   | 26 Fr, 16 in. (40.64 cm) | 30 cc    | 10/Box   |
| AS42028   | 28 Fr, 16 in. (40.64 cm) | 30 cc    | 10/Box   |
| AS42030   | 30 Fr, 16 in. (40.64 cm) | 30 cc    | 10/Box   |

## 2-WAY 100% SILICONE FOLEY CATHETERS

- Smooth eyes improve patient comfort
- Reinforced tip provides a stylet for optimized insertion efficacy
- 100% silicone material reduces encrustation and calcifications, making it ideal for longer-term catheterization
- Symmetrical balloon provides a tight seal at the neck of the urethra
- Sterile, latex-free, DEHP-free and BPA-free



HPCS code A4344

| ITEM CODE | DESCRIPTION              | FEATURES | QTY/UNIT |
|-----------|--------------------------|----------|----------|
| AS41012S  | 12 Fr, 16 in. (40.64 cm) | 5 cc     | 10/Box   |
| AS41014S  | 14 Fr, 16 in. (40.64 cm) | 5 cc     | 10/Box   |
| AS41016S  | 16 Fr, 16 in. (40.64 cm) | 5 cc     | 10/Box   |
| AS41018S  | 18 Fr, 16 in. (40.64 cm) | 5 cc     | 10/Box   |
| AS41020S  | 20 Fr, 16 in. (40.64 cm) | 5 cc     | 10/Box   |
| AS41022S  | 22 Fr, 16 in. (40.64 cm) | 5 cc     | 10/Box   |
| AS41024S  | 24 Fr, 16 in. (40.64 cm) | 5 cc     | 10/Box   |
| AS41026S  | 26 Fr, 16 in. (40.64 cm) | 5 cc     | 10/Box   |
| AS41028S  | 28 Fr, 16 in. (40.64 cm) | 5 cc     | 10/Box   |
| AS42016S  | 16 Fr, 16 in. (40.64 cm) | 30 cc    | 10/Box   |
| AS42018S  | 18 Fr, 16 in. (40.64 cm) | 30 cc    | 10/Box   |
| AS42020S  | 20 Fr, 16 in. (40.64 cm) | 30 cc    | 10/Box   |
| AS42022S  | 22 Fr, 16 in. (40.64 cm) | 30 cc    | 10/Box   |
| AS42024S  | 24 Fr, 16 in. (40.64 cm) | 30 cc    | 10/Box   |

# FOLEY INSERTION TRAYS

## FOLEY INSERTION TRAY

- Deep tray basin designed to prevent contamination
- Gloves and povidone-iodine swab sticks help maintain sterility and minimized risk of infection
- Sterile, latex-free, DEHP-free and BPA-free

Foley Insertion Trays include:

- 1000 mL collection tray
- (3) Povidone-iodine swab sticks and benzalkonium chloride prep pad
- Vinyl powder-free gloves
- Waterproof underpad 13 in. x 18 in. (33.02 cm x 45.72 cm)
- Lubrication jelly
- Fenestrated drape
- Pre-filled syringe with sterile water

HPCS code A4310



| ITEM CODE | 10 cc SYRINGE | 30 cc SYRINGE | QTY/UNIT |
|-----------|---------------|---------------|----------|
| AS880     | ●             |               | 20/Case  |
| AS890     |               | ●             | 20/Case  |

## ADD-A-FOLEY CATHETER TRAY

- Needle-free sampling port on urinary bag reduces the risk of needle stick injuries
- Gloves and povidone-iodine solution help maintain sterility and reduce the risk of infection
- Sterile and latex-free

Add-A-Foley Trays include:

- 2000 mL low-profile urinary drainage bag
- Waterproof and fenestrated drapes
- Nitrile powder-free gloves
- (5) Rayon balls
- 10% povidone-iodine solution
- Forceps
- Lubricating jelly
- 2 oz. specimen container with patient label
- Central supply room wrap and Tyvek header bag
- Pre-filled syringe with sterile water

HPCS code A4354



| ITEM CODE | 10 cc SYRINGE | 30 cc SYRINGE | QTY/UNIT |
|-----------|---------------|---------------|----------|
| AS89300   | ●             |               | 10/Case  |
| AS89330   |               | ●             | 10/Case  |

# FOLEY INSERTION TRAYS

## PRE-CONNECTED DRAINAGE BAG FOLEY CATHETER TRAYS

- Sterile tray includes all necessary components for aseptic catheter insertion
- Tray improves clinical efficacy through cost-effective inventory management
- Needle-free urinary sampling port reduces the risk of needle stick injuries
- Tray is available with both siliconized latex and 100% silicone catheters to suit patient needs

Pre-connected drainage bag Foley tray contents include:

- 2000 mL low-profile urinary drainage bag
- Pre-connected 16 in. Foley catheter with 5 cc balloon
- Waterproof and fenestrated drapes
- Nitrile powder-free gloves
- (5) Rayon balls
- 10% povidone-iodine solution
- Forceps
- Lubricating jelly
- Syringe prefilled with 10 cc sterile water
- 2 oz. specimen container with a patient label
- Central supply room wrap and Tyvek header bag



AS89316

HPCS code A4314

| ITEM CODE | FOLEY CATHETER FR SIZE |       |       | CATHETER MATERIAL |                            | QTY/UNIT |
|-----------|------------------------|-------|-------|-------------------|----------------------------|----------|
|           | 14 Fr                  | 16 Fr | 18 Fr | Siliconized Latex | 100% Silicone (Latex-free) |          |
| AS89314   | ●                      |       |       | ●                 |                            | 10/Case  |
| AS89316   |                        | ●     |       | ●                 |                            | 10/Case  |
| AS89318   |                        |       | ●     | ●                 |                            | 10/Case  |
| AS89314S  | ●                      |       |       |                   | ●                          | 10/Case  |
| AS89316S  |                        | ●     |       |                   | ●                          | 10/Case  |
| AS89318S  |                        |       | ●     |                   | ●                          | 10/Case  |



# SYRINGES, IRRIGATION TRAYS AND KITS

## BULB SYRINGE

- Designed to stand on the bulb end to help minimize contamination
- Tip protector included


## THUMB CONTROL RING SYRINGE

- Re-sealable pole bag designed for safe and convenient storage
- Includes syringe luer tip adapter and protector cap

## FLAT TOP SYRINGE

- Visible dosing marks designed for ease of use
- Includes syringe luer tip adapter and protector cap

HPCS code A4322



| ITEM CODE | 60 mL SYRINGE TYPE |                    |          | PACKAGE TYPE |                |                     |      | STERILE | QTY/UNIT |
|-----------|--------------------|--------------------|----------|--------------|----------------|---------------------|------|---------|----------|
|           | Bulb               | Thumb Control Ring | Flat Top | Poly Pouch   | Form Fill Pack | Resealable Pole Bag | Bulk |         |          |
| AS116     |                    |                    | •        |              |                | •                   |      |         | 30/Case  |
| AS016     |                    | •                  |          |              |                | •                   |      |         | 30/Case  |
| AS010     | •                  |                    |          |              |                |                     | •    |         | 50/Case  |
| AS011     | •                  |                    |          | •            |                |                     |      | •       | 50/Case  |
| AS011P    | •                  |                    |          |              | •              |                     |      | •       | 50/Case  |
| AS015     |                    | •                  |          | •            |                |                     |      | •       | 50/Case  |
| AS115     |                    |                    | •        | •            |                |                     |      | •       | 50/Case  |

Note: All syringes are latex-free, DEHP-free and BPA-free

## IRRIGATION TRAYS AND KITS

Non-sterile Kits include:

- 500 mL graduated container
- Patient ID label and tip protector

HPCS code A4320



| ITEM CODE | THUMB CONTROL RING SYRINGE | BULB SYRINGE | RESEALABLE POLE BAG | QTY/UNIT |
|-----------|----------------------------|--------------|---------------------|----------|
| AS121     |                            | •            |                     | 30/Case  |
| AS126     | •                          |              |                     | 30/Case  |
| AS127     | •                          |              | •                   | 30/Case  |

Sterile Trays include:

- 1000 mL outer tray
- 500 mL graduated container
- Waterproof underpad
- Alcohol prep pad
- Patient ID label and tip protector

HPCS code A4320



| ITEM CODE | THUMB CONTROL RING SYRINGE | BULB SYRINGE | QTY/UNIT |
|-----------|----------------------------|--------------|----------|
| AS130     |                            | •            | 20/Case  |
| AS136     | •                          |              | 20/Case  |

# URINARY COLLECTION DEVICES

## URINARY DRAINAGE BAG

- Anti-reflux mechanism helps improve patient safety by restricting urine backflow
- Leak-proof seams ensure fluid containment
- Needle-free sampling port reduces the risk of needle stick injuries
- Air vents help to maintain consistent urine flow, delivering safer patient outcomes
- Sterile fluid pathway, latex-free and BPA-free



AS322 Teardrop



AS32600 Low-profile

HPCS code A4358

| ITEM CODE | BAG SHAPE   |          | ANTI-REFLUX TYPE |               |        | NEEDLE-FREE SAMPLING PORT |                             | DRAIN TYPE |             | HANGER TYPE |             |      | QTY/UNIT |
|-----------|-------------|----------|------------------|---------------|--------|---------------------------|-----------------------------|------------|-------------|-------------|-------------|------|----------|
|           | Low Profile | Teardrop | Chamber          | Flutter Valve | Device | Pre-pierced Blunt Cannula | Male Luer Lock or Luer Slip | T-tap      | Metal Clamp | Single Hook | Double Hook | Rope |          |
| AS312     |             | •        | •                |               |        | •                         |                             | •          |             |             | •           | •    | 20/Case  |
| AS322     |             | •        | •                |               |        | •                         |                             | •          |             |             | •           | •    | 20/Case  |
| AS32200   | •           |          | •                |               |        | •                         |                             | •          |             | •           |             | •    | 20/Case  |
| AS32600*  | •           |          | •                |               |        |                           | •                           | •          |             | •           |             | •    | 20/Case  |
| AS332     |             | •        |                  | •             |        | •                         |                             | •          |             | •           |             |      | 20/Case  |
| AS326     | •           |          |                  |               | •      | •                         |                             |            | •           |             | •           | •    | 20/Case  |
| AS302     |             | •        |                  |               | •      | •                         |                             | •          |             |             | •           | •    | 20/Case  |

\*100% sterile

## PEDIATRIC URINE COLLECTOR

- Adhesive surface folds to form a spill-proof specimen transport container
- Hypoallergenic feature makes it safe to use with most pediatric patients
- Visible graduation marks designed for ease of use
- Single-use disposable
- Sterile, latex-free, DEHP-free and BPA-free

HPCS code A4328

| ITEM CODE | DESCRIPTION                 | QTY/UNIT |
|-----------|-----------------------------|----------|
| AS409     | 200 mL, Pediatric Urine Bag | 50/Box   |



AS409

# URINARY COLLECTION DEVICES

## URINARY LEG BAGS

- Soft and breathable pre-attached leg straps provide patient comfort
- Anti-reflux device helps improve patient safety by reducing urine backflow
- Leak-proof seams ensure fluid containment
- Sterile fluid pathway, latex-free and BPA-free




AS307N

HPCS code A4358

| ITEM CODE | 600 mL LEG BAG | 900 mL LEG BAG | PUSH-PULL DRAIN PORT | TWIST-TURN DRAIN PORT | 18 in. (45.72 cm) EXTENSION TUBE WITH ADAPTER AND PROTECTIVE CAP | QTY/UNIT |
|-----------|----------------|----------------|----------------------|-----------------------|--|----------|
| AS306N    | •              |                | •                    |                       |  | 48/Case  |
| AS307N    |                | •              | •                    |                       |  | 48/Case  |
| AS308N    | •              |                |                      | •                     |  | 48/Case  |
| AS309N    |                | •              |                      | •                     |  | 48/Case  |
| LBE306N   | •              |                | •                    |                       | •  | 48/Case  |
| LBE307N   |                | •              | •                    |                       | •  | 48/Case  |
| LBE308N   | •              |                |                      | •                     | •  | 48/Case  |
| LBE309N   |                | •              |                      | •                     | •  | 48/Case  |

## EXTENSION TUBE FOR LEG BAGS

- Kink-resistant tubing safeguards against fluid buildup
- Flexible tubing permits leg bag to be positioned on the thigh or calf
- Sterile, latex-free and BPA-free




AS311

HPCS code A4331

# SPECIMEN CONTAINERS

- 4 oz. leak resistant specimen container
- Latex-free, DEHP-free, BPA-free


| ITEM CODE | FORM FILL SEAL | POLY POUCH | BULK | STERILE INTERIOR | TAMPER EVIDENT LABEL | PATIENT LABELS | QTY/UNIT |
|-----------|----------------|------------|------|------------------|----------------------|----------------|----------|
| AS340     |                | •          |      | •                |                      | •              | 100/Case |
| AS341     |                |            | •    | •                | •                    |                | 100/Case |
| AS343     |                |            | •    |                  |                      |                | 500/Case |
| AS345     | •              |            |      | •*               |                      |                | 100/Case |



AS340

\*100% sterile

## Your Local Distributor



### **Amsino International, Inc.**

708 Corporate Center Drive  
Pomona, CA 91768

P: 800.632.6746

F: 909.626.3888

[info@amsino.com](mailto:info@amsino.com)

[www.amsino.com](http://www.amsino.com)

