

Delta-Cast® Elite

Fibreglass-Free Cast Tape

Delta-Cast Elite



Fibreglass-Free Cast Tape

Indications of Use

- Primary & secondary casting
- Rigid support

Unique Features

- Fibreglass-free
- Clear radiograph viewing
- Excellent lamination
- Smooth finish

Delta-Cast® Elite

Fibreglass-Free Cast Tape

Used for rigid primary and secondary casting applications where enhanced radiograph viewing is important.

Delta-Cast Elite has been designed with ease of application and patient comfort in mind. Our patented fibreglass-free casting material combines strength and rigidity with soft edges for enhanced patient comfort.

Patented Fibreglass-Free Substrate

Delta-Cast Elite's patented polyester substrate features extensible yarns that easily shape to body contours, providing enhanced molding capabilities for both primary and secondary casting applications and produces less dust during cast removal.

Patented Resin System

Ensures maximum lamination and end laydown for strong, long-lasting casts. The patented glove system provides a tack-free application.



Application



Ordering Information

Order Number	Size	Colour	Unit of Measure
6712	5cm x 3.6m	WHITE	10 rolls/box
6713	7.5cm x 3.6m		10 rolls/box
6714	10cm x 3.6m		10 rolls/box
6715	12.5cm x 3.6m		10 rolls/box
4812	5cm x 3.6m	DARK BLUE	10 rolls/box
4813	7.5cm x 3.6m		10 rolls/box
4814	10cm x 3.6m		10 rolls/box
4862	5cm x 3.6m	PURPLE	10 rolls/box
4863	7.5cm x 3.6m		10 rolls/box
4864	10cm x 3.6m		10 rolls/box
4802	5cm x 3.6m	RED	10 rolls/box
4803	7.5cm x 3.6m		10 rolls/box
4804	10cm x 3.6m		10 rolls/box
4872	5cm x 3.6m	BLACK	10 rolls/box
4873	7.5cm x 3.6m		10 rolls/box
4874	10cm x 3.6m		10 rolls/box