



Philips InnoSpire Essence
Compressor Nebulizer
System

Compact home-use device
Fast, effective treatment
SideStream Disposable included

HH1338

Reliable nebulizer therapy you can count on

For consistent treatment

InnoSpire Essence provides fast, efficient and reliable aerosol therapy, designed for home use.

Proven SideStream technology

- Proven nebulizer technology
- Robust design meets reliability

Short treatment times

- Treatment delivered in 6-8 minutes*

For use in the home

- Designed for intermittent use
- Compatible with commonly prescribed nebulizer medications

PHILIPS

Highlights

For home use



Designed for intermittent use and suitable for the home, InnoSpire Essence provides the consistent delivery high-quality aerosol treatment that you can rely on.

Power and style all in one



InnoSpire Essence is designed to both fit your lifestyle and cope with frequent use. You can be confident that you will receive high quality aerosol for the life of the nebulizer.

Proven nebulizer technology



Reduced Treatment Times - Our Sidestream nebulizer is specially designed to keep your treatment times as short as possible. As the InnoSpire Essence pushes air through the tube, additional air is drawn in through the top of the nebulizer to create a fast flow of aerosol medication for you to inhale. Treatments can take as little as 6 minutes* Consistent delivery of high quality aerosol therapy - Sidestream is designed to withstand frequent and repeated use and will continue to efficiently deliver your nebulizer medication until it is time for replacement.

Versatility



Capable of nebulizing a wide variety of medications that are routinely with nebulizers used to treat various respiratory conditions.

Treatment time

InnoSpire Essence, combined with our SideStream nebulizer, is designed to provide fast drug delivery and short treatment times.

Specifications

Service

- Warranty: 3 years

Design

- Compact: Weighs only 1.5 kg

Package contents

- Include: SideStream Disposable, InnoSpire Essence compressor

Operating parameters

- Electrical rating/power: 230 V / 50 Hz, 1 Amp
- Maximum flow rate: 8 LPM
- Maximum pressure: 45 psi
- Mode of operation: Intermittent (30 min. on/off)
- Operating flow rate: 6 LPM @ 10 psi

Product details

- Weight: 1.5 kg
- Size: 165x165x108 mm
- Noise level: 58 ±3 dBA
- Treatment time (6-8 minutes*): * 2.5ml Salbutamol dose

Performance data

- MMAD: 2.9 µm
- Fine particle fraction: 77% below 5 µm
- Nebulization capacity: 8 ml



Issue date 2018-07-24

Version: 3.0.1

12 NC: 8670 001 46077
EAN: 03 83730 26225 3

© 2018 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

** Using a 2.5ml Salbutamol dose