



## IV Administration Set With GVS EasyDrop<sup>®</sup> Flow Regulator

*Our IV Administration Set with Flow Regulator consists of a universal spike, male swivel luer lock, roller clamp, flow regulator, and drip chamber. Each set is individually blister-packaged in quantities of 40 per box.*

- 20 drops/mL
- One (1) Y-injection site 6" from distal end
- Male swivel luer lock connector
- Universal vented/non-vented spike
- Roller clamp
- 15µm filter located in drip chamber
- GVS EasyDrop<sup>®</sup> Flow Regulator
- 92" length
- DEHP-Free
- Latex-Free
- Sterile



Product # TCBINF033G

Visit [www.tcbiomedix.com](http://www.tcbiomedix.com) and click on the Category icon to read detailed product descriptions, interactively examine each product and download PDFs of group or individual Product Spec Sheets