

SLING COLOR/SIZE GUIDE

	Size	Lift	Color	Weight Range
Pediatrics	XXS		Teal	0-55 lbs (0-25 kg)
	XS		Brown	55-77 lbs (25-35 kg)
	S		Red	77-132 lbs (35-60 kg)
Adults	XS		Brown	55-77 lbs (25-35 kg)
	S		Red	77-132 lbs (35-60 kg)
	M		Yellow	121-165 lbs (55-75 kg)
	L		Green	154-264 lbs (70-120 kg)
	LL		Purple	220-350 lbs (100-160 kg)
	XL		Blue	308-440 lbs (140-200 kg)
	XXL		Terracotta	440-500 lbs (200-228 kg)
	Bariatric	Bariatric M		Yellow/Dark Gray
Bariatric L		Green/Dark Gray		450-555 lbs (210-250 kg)
Bariatric XL		Blue/Dark Gray		500-750 lbs (230-340 kg)
Bariatric XXL		Terracotta/Dark Gray		700-1000 lbs (318-455 kg)

ARJO makes sizing and sling selection easy with a color-coded sling system. The sling sizes covering pediatric, adult and bariatric needs each have a distinctive color code that can be clearly seen on the edge of ARJO slings. Please display this chart according to HIPAA regulations.

ARJO ...with people in mind.

Lift Type:

Maxi Move* Sara Plus Sara 3000 Tenor Maxi Sky



* Note: The *Maxi Move* has various size hanger bars. It is important that the user select the correct sling size prior to using the *Maxi Move* to assist with transfer of a patient.