


# ATLANTA™

## TECHNICAL SPECIFICATION



### Basic dimensions

<b>A</b>	Overall height	1135mm/45in
<b>B</b>	Back height	750mm/30in
<b>C</b>	Seat height	540mm/21in
<b>D</b>	Inside back width	505mm/20in
<b>E</b>	Overall length	1350mm/53in
<b>F</b>	Reclined length	1700mm/67in

### Available seat and overall widths

<b>G (Seat width)</b>	<b>H (Overall width)</b>
380mm/15in	635mm/25in*
450mm/18in	690mm/27in*
500mm/20in	690mm/27in*
550mm/22in	740mm/29in*
600mm/24in	790mm/31in*

\* Dimension includes legrest belt. Deduct 35mm/1.5in if belt is not used.

### Arm height options

<b>I</b>	Low	630mm/25in
<b>J</b>	High	660mm/26in

### Seat depth

<b>K</b>	480mm/19in
----------	------------

### Leg rest length

<b>L</b>	420mm/16.5in
----------	--------------

### Clearance Height

<b>M</b>	75mm/3in
----------	----------

### Recline and tilt angles

#### ii Back angle recline

110° - 140°

30° range

#### iii Leg rest angle



75° - 107°

32° range

#### iv Integrated tilt and recline

54° - 74°

20° range



Chair weight unloaded is 67 Kg/147.71 lb/10.56 Stone.  
Chair user weight limit is 159 Kg/364 lb/25 Stone