


Sizes of Sara 3000 active slings , figure 1



Sizes of Sara 3000 active slings in centimeters
--

SIZE	COMMENTS	TSS.502	TSS.501	TSS.500	TSS.511
A	width anti slip backsupport	16,9	13,8	10,6	15,0
B	width polstered back part	49,6	39,4	34,6	38,6
C	length chest fixation right	24,4	18,9	14,6	18,9
D	length chest fixation left	20,5	14,2	10,6	14,6
E*	<i>circumference chest fixation</i>	56,7	46,1	37,4	46,1
F	length from clips to clips 1	55,9	43,3	38,6	42,5
G	length from clips to clips 2	64,6	47,2	42,5	46,5
H	height anti slip backsupport	15,4	14,6	11,8	13,4
I	length adjustable strap	39,4	35,4	28,3	35,4

E* These numbers are averages as the chest fixation belt can be adjusted and a number of lengths are possible. E is not shown on the above figure 1.

For the TSS.511, the length of the supporting and polstered leg part is 28,3 inches, measured from the top seam of the backsupport. The length including the strap with clips is 55,1 inches.