


# Free standing system for portable ceiling lifts

## Easytrack FS

The Easytrack FS is a freestanding, versatile solution for all types of point-to-point transfers - wheelchair to bed, bed to chair, even in the bathroom, anywhere there is a lift and transfer need.

The Easytrack FS has a safe working capacity of 440 lb (200 kg). That's exceptional for a portable product making it suitable for use in almost any home. This equipment was designed to allow those with a mobility impairment to continue to live at home and enjoy normal activities.



# Product	Description
9720010	Easytrack FS - for portable ceiling lifts

√ = Included / o = Optional

# Part	Description	
1a	700.17250 Extendable track with trolley - 1 piece	√
1b	700.17100 Base - 2 pieces	√
1c	700.17000 Post - 2 pieces	√

Optional Protective Bag kit		
2a	002.17010 Protective Bag Kit (Including : #002.17100, #002.17250 & 2 X #002.17000)	o
2b	002.17250 1 Rail bag	o
2c	002.17100 1 Base bag	o
2d	002.17000 1 Post bag	o

Accessories		
3a	700.08310 Reacher - 61 cm (24 in)	o
3b	700.08320 Reacher - 91 cm (36 in)	o

