

A three-section treatment table with a white metal frame and red padded sections. The backrest is tilted upwards, and the footrest is extended. A red line points to the backrest with the text 'ELEVATES 73°'. Another red line points to the footrest with the text 'RAISES 45° LOWERS 90°'.

ELEVATES 73°

RAISES 45°
LOWERS 90°

AM-350 Treatment Table

- 3-Section Treatment Table
- 76"L x 27"W
- Min. Height: 18"
- Max Height: 37"
- 1 1/2" Firm Density Foam
- Head: 15" x 27" with contoured face opening
- Center: 22" x 27"
- Foot: 39" x 27"
- 400lbs. Capacity
- Power-assisted Caster System