



Machine Stand Accomodates
All Traction Machines

ELEVATES 80°

RAISES 45°
LOWERS 90°

- Friction Free Lumbar Separates 5" And Rolls On 4 Steel Ball Roller Bearings.
- The Rolling Section Can Be Locked In Any Position With A Positive Gear Rack System.

AM-400 Traction Table

- | | | | |
|-----------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------|
| <ul style="list-style-type: none"> • 4-Section Traction Table • 76"L x 27"W • 99" Overall Length | <ul style="list-style-type: none"> • Min. Height: 18" • Max Height: 37" • 1 1/2" Firm Density Foam | <ul style="list-style-type: none"> • Head: 15" x 27" with contoured face opening • Thoracic: 17.5" x 27" • Center: 10" x 27" • Foot: 33" x 27" | <ul style="list-style-type: none"> • 400lbs. Capacity • Power-assisted Caster System |
|-----------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------|