

Name:	
School/Organization:	
Dealer/Rep:	

Phone:	
E-mail:	
Sales Rep:	Quote:

# AM-BA SERIES

## 34" BARIATRIC HI-LO TREATMENT TABLES

### STANDARD FEATURES:

- Load capacity: 500 lbs
- 34"W x 76"L x 18-37"H
- 34 oz. Heavy-duty vinyl with Permablok 3® bacterial protection
- 1.5" High density foam
- Rounded corners
- Steel base with dual foot bar control
- White powder coat base



AM-BA234

### TABLE CONFIGURATIONS:

AM-BA134 - 1 Section	AM-BA234 - 2 Section	AM-BA334 - 3 Section
QTY	QTY	QTY



### TABLE FEATURES:

<p>Flat Cushion</p>	<p>Lower Extremity Angle Ranges from 0 to 75°</p>	<p>Lower Extremity Angle Ranges from 0 to 78°</p> <p>Head Section Angle Ranges from +45° to -90°</p>
AM-BA134	AM-BA234	AM-BA334














800.701.5122  
armedicamfg.com

**OPTIONS:**

<b>AM-804</b> - Paper Holder with Cutter (Accommodates 18-21" Paper Rolls)
<b>AM-18</b> - Face Cutout
<b>AM-17</b> - Face Plug (Must order AM-18 also if ordering AM-BA134)
<b>AM-91</b> - Add 1" Extra Soft Foam
<b>19855</b> - Add-on Belt & Brackets
<b>SPCO</b> - Special Vinyl Color (Provided by Customer)
<b>SPBASE</b> - Special Base Color (Select Custom Base Color)
<b>12088</b> - Side Panel 9 x 33.25 Perforated
<b>12075</b> - Side Panel 9 x 33.25 Flat




**STANDARD UPHOLSTERY CHOICES:**

 Imperial Blue	 Taupe
 Blue Ridge	 Burgundy
 Forest Green	 Dove Gray
 Black	 Greystone
 Cappuccino	 Merlot
 Patina	

**STANDARD BASE COLOR:**

 Office White
--

**CUSTOM BASE COLORS:**

 Black
 Steel Gray
 Metallic Silver