

Name:	
School/Organization:	
Dealer/Rep:	

Phone:			
E-mail:			
Sales Rep:		Quote:	

# AM SERIES

## 40" BOBATH HI-LO TREATMENT TABLES


### STANDARD FEATURES:

- Load capacity: 500 lbs
- 40"W x 76"L x 18-37"H
- 34 oz. Heavy-duty vinyl with Permablok 3® bacterial protection
- 1.5" High density foam
- Rounded corners
- Single foot control pedal
- White powder coat base



### TABLE CONFIGURATIONS:

AM-140 - 1 Section	AM-240 - 2 Section	AM-340 - 3 Section
QTY	QTY	QTY



### TABLE FEATURES:

<p>Flat Cushion</p>	<p>Lower Extremity Angle Ranges from 0 to 74°</p>	<p>Lower Extremity Angle Ranges from 0 to 78°</p> <p>Head Section Angle Ranges from +45° to -90°</p>
AM-140	AM-240	AM-340














800.701.5122  
armedicamfg.com

**OPTIONS:**

<b>AM-1036</b> - Add Base with 4 Total Lock Swivel Casters
<b>AM-1062</b> - Add Base with 4 Casters and 6" Clearance (Adds 4")
<b>AM-804</b> - Paper Holder with Cutter (Accommodates 18-21" Paper Rolls)
<b>AM-18</b> - Face Cutout
<b>AM-17</b> - Face Plug (Must order AM-18 also if ordering AM-140, & AM-240)
<b>AM-91</b> - Add 1" Extra Soft Foam
<b>AM-807</b> - Rail Side Fold Down (Pair)
<b>19855</b> - Add-on Belt & Brackets
<b>SPCO</b> - Special Vinyl Color (Provided by Customer)
<b>SPBASE</b> - Special Base Color (Select Custom Base Color)
<b>12088</b> - Side Panel 9 x 33.25 Perforated
<b>12075</b> - Side Panel 9 x 33.25 Flat




**STANDARD UPHOLSTERY CHOICES:**

	Imperial Blue		Taupe
	Blue Ridge		Burgundy
	Forest Green		Dove Gray
	Black		Greystone
	Cappuccino		Merlot
	Patina		

**STANDARD BASE COLOR:**

	Office White
---	--------------

**CUSTOM BASE COLORS:**

	Black
	Steel Gray
	Metallic Silver