

Name:	
School/Organization:	
Dealer/Rep:	

Phone:			
E-mail:			
Sales Rep:		Quote:	

# FOOTSTOOL

## WITH NON-SKID RUBBER TREAD TOP & HANDRAIL

### STANDARD FEATURES:

- 11"W x 14"L x 34"H
- Welded steel construction
- Durable chrome finish
- Surface platform: 11" x 14"
- Non-skid rubber treaded mat: 10" x 13"

### STANDARD STOOL:

#### AM-843 - Footstool with Non-Skid Rubber Tread Top & Handrail

- 11"W x 14"L x 34"H
- Load capacity: 400 lbs.
- Non-skid rubber treaded mat
- Weighs 12 lbs

QTY:



### BARIATRIC FOOT STOOL:

#### AM-849 - Footstool with Non-Skid Rubber Tread Top & Handrail

- 11"W x 14"L x 34"H
- Load capacity: 700 lbs.
- Non-skid rubber treaded mat
- Reinforced welded steel base
- Weighs 13.4 lbs

QTY:



**Armedica**  
BY PIVOTAL HEALTH SOLUTIONS

800.743.7738  
armedicamfg.com