

# Mobility - Posture Control for Adults | Posture Rest Walkers



## 'H' Frame Walker

This walker style, which includes a seat that easily folds up and down, is helpful to users who have low endurance or users who are in environments where they transition between walking and sitting often.

### Available Accessories

- Extensor Assist Pad
- Forearm Supports
- Guide Handles
- Pelvic Stabilizer
- Rotational Handgrips
- Swivel Limiters

Model No.	Height to top of handle (floor to mid-buttocks)	Distance inside handles (depth/width)	Weight Limits (lbs/kg)
-----------	--	--	---------------------------

## Posture Rest 'H' Frame Walkers

W1/2BH	14.5" to 20"	9" / 13.5"	60 / 27
W1BH	16.5" to 22"	9" / 13.5"	60 / 27
W2H	20" to 26.5"	11" / 16.25"	85 / 39
W3H	25.5" to 32"	12" / 19"	130 / 60
W4H	31" to 37"	14" / 21"	180 / 82
W5H	36" to 41"	13.5" / 20"	250 / 113
W5CH	35.5" to 41"	15" / 18.5"	250 / 113