

# Seating - Benches

## Kaye Adjustable Benches

Therapists use benches in physical, occupational and speech therapy, particularly when the child needs a stable sitting surface to achieve maximum arm and hand movements for play, learning, and daily living activities. The Kaye Adjustable Benches are ideal for home, therapy, or school programs and for children who do not need back support in sitting but need a firm seat that can be adjusted so that the child's feet rest on the floor.



**Models S3A and S4A fold for easy storage and transportation!**



S1A



S5A



S2A

Model No.	User weight limits	Height adjusts	Tilts	Folds	Seat size	Posture system
S1A	180 lbs (82 kg)	7" to 11"	Yes	No	10.5" x 26"	S1AO
S3A	180 lbs (82 kg)	7" to 11"	Yes	Yes	10.5" x 26"	N/A
S2A	180 lbs (82 kg)	10.5" to 17"	Yes	No	14" x 31"	S2AO, H2AO
S4A	180 lbs (82 kg)	10.5" to 17"	Yes	Yes	14" x 31"	N/A
S5A	180 lbs (82 kg)	16" to 23"	Yes	No	14" x 31"	S2AO, H2AO
TSS	180 lbs (82 kg)	7" to 11"	Yes	No	10.5" x 12"	S1AO
TSL	180 lbs (82 kg)	10.5" to 17"	Yes	No	14" x 17.5"	S2AO, H2AO