

## Kaye Posture System

The Kaye Posture System can be added to Kaye Benches and T-Seats to provide additional stability and improved alignment in sitting. The Posture System consists of two pieces, a pelvic support and a femoral stabilizer, which can be used together or separately. The pelvic support component comes with lateral supports and can be adjusted in height. It is also available with a high back for users who require more postural support.

The femoral stabilizer has a built-in abductor to prevent both posterior pelvic tilt and rotation or obliquity of the pelvis while stabilizing the hips and thighs in neutral alignment.



S1A with S1AO



TSL with H2AO

Model Features	S1AO	S2AO	H2AO
Benches / Seats	S1A, TSS	S2A, TSL, S5A	S2A, S5A, TSL
Adjustable seat depth	5" to 10.5"	5" to 10.5"	7" to 10"
Back dimensions	9"w x 5"h	9"w x 5"h	14"w x 14"h
Side cushions	4" x 4"	4" x 4"	4" x 4"
Width between side supports	5" to 11"	5" to 11"	7" to 15"
Knee Pads	4"w x 7"h	4"w x 7"h	4"w x 7"h
Abductor	3"w	3"w	3"w