

Kaye Play Frame

The Kaye Play Frame is designed to encourage children to experiment with various strategies for reaching, batting, and swiping or kicking, skills that develop very early in infancy. Toys can be suspended and strategically positioned to encourage the child to develop variety in shoulder and arm or leg motions.



Scan code to see Evidence Reports.

33-inch (blue) ball on model BH33

24-inch (green) ball on model BH16



Ball & Roll Holders

Holders will accommodate rolls and balls. Holders are easily mounted to the wall and all mounting hardware is included.

Model No.	Width x Depth	Holds Ball /Roll Sizes
BH16	16.25" to 19.25"	20" to 24"
BH33	28" to 30"	33" to 49"