

My Red Walker | Walker Wagon - Mobility

Kaye 'Red Walker'

The Kaye Red Walker is the newest Posture Control Walker with a family-friendly design to appeal to young walkers and their families.

Both sizes have swivel front wheels, which can be locked for straight ahead tracking, and one way ratchet rear wheels.

Accessories for this walker include forearm supports, pelvic stabilizer, swivel limiters, Kaye Extensor Assist Belt, and soft sling support.

Folds for transporting and storage.



Model Features	RFA	RFB
Height to top of handgrip	13.5" to 18.5"	17" to 25"
Distance Inside Handles (depth / width)	9.5" / 13.5"	9.5 / 13.5"
Weight Limits (lbs/kgs)	60/27	85/39
Swivel Limiters	WSL-1R	WSL-1R
Forearm Supports - FSA	Adjusts from 4" to 7.75"	Adjusts from 4" to 7.75"
Pelvic Stabilizer - PS2H	Adjusts from 8.5" to 11.5"	Adjusts from 8.5" to 11.5"
Extensor Assist Pad - EA2	6"w x 3.5"h	6"w x 3.5"h
Soft Sling Support	SSS1	SSS1