

Scot-About

The Scot-About has numerous clinical applications ranging from a seating surface for therapists as they facilitate gait to a mobile surface for children to use in multiple positions. The Scot-About has a padded seat and handgrips. It can be used in conjunction with an adjustable posterior trunk support.



SB1

Model Features	SB1
Padded seat	16.5" w x 12.5" d
Adjustable seat height	8.5" to 16.5"
Handgrips	5" above seat
Swivel casters	5"
Overall width	24"



SB1 with SBAO

The new, adjustable, removable back support for the Kaye Scot-About (Model SB1) provides support at the pelvis for children when seated on the Scot-About. The back and side supports adjust in height and width to give optimal vertical alignment for the pelvis and trunk.

This accessory makes the Scot-About a useful mobile seat for children as they practice propelling themselves with their feet, yet is removable so that the Scot-About can still be used by parents and therapists.

Model Features	SBAO
Adjustable seat depth	5" to 12"
Back dimensions	9" w x 5" h
Side cushions	4" x 4"
Width between side supports	5" to 11"