

Therapy Bolsters & Stands

The Kaye Bolsters are designed to be used in therapy by infants and children who need to develop balance and posture control. They are available in two sizes, as detailed in the table below.

The versatility of the bolsters increases when they are combined with the CB20 Kaye Bolster Stands. Bolsters then can be angled to add complexity or support without losing the mobility of the bolster. Bolster stands are adjustable in height (5.5" to 11.25" from floor) and are 7" wide and 12.5" long.

The bolsters are covered with washable vinyl.



Model No.	Width x Length	Stand / Bolster
CB15	8.5" x 43"	Bolster
CB16	10" x 50"	Bolster
CB20	7" x 12.5"	Stand