

# Seating - Tilting Tables

Occupational and physical therapists emphasize the importance of an adjustable-height work surface.

## Tilting Tables

This product was designed to provide an adjustable-height, tilting surface to be used in a therapy, home, or school setting. The T-Table allows the child to be positioned at the optimum distance from the work surface. Two sizes are available, large and small.



Model Features	TT1 (small)	TT2 (large)
Table Surface	14.75" d x 20.75"	20" d x 29.75" w
Height from floor	12" to 17"	17.75" to 25.5"
Width under tables	14.5"	21.5"
Span between legs	14.5"	21.5"
Surface tilts	0 to 45°	0 to 45°
Adjustable clamps	Yes	Yes
Removable front rail	Yes	Yes
Interfaces with bolsters	Yes	Yes
Interfaces with Kaye Therapy Benches	Yes	Yes
Interfaces with T-Seats	Yes	Yes