

Vertical Stander

Designed as positioning aids for users who do not have adequate balance to stand independently and who need lower body support to attain alignment while standing. Some of the reported benefits of standing include a decrease in the incidence of contractures and fractures. The Kaye Vertical Stander is designed to accommodate children 18 months through 9 years.



VS1+C and VS1-T



VS1+C and VS1-T

Model Features	VS1
Age	18 months to 9 years
Base width / length (fully extended)	20" / 29.25" (retractable base to fit against counters)
Height from footrest to top of pelvic support	8" to 30.5"
Pelvic stabilizer inside depth	3.5" x 8.5"
Pelvic stabilizer inside width	5" to 11.75"
Height-adjustable knee cushion	One included, 5"h x 14"w (removable)
Thoracic support	Removable, height-adjustable
Shoe clamps	Adjust to angle to user's toes in or out