



Power You Can Depend On.

# Charger Reference Guide

MK Battery has a complete

line of premium quality

battery chargers for all

of your charging needs.

They are available for

immediate shipment

or delivery from our

warehouses located

throughout North America.



**MK Model: M-D-EN0801**  
(Specify chair type or plug style) (Pictured Above)

- **DESCRIPTION:** 24 volt, 8 amp automatic charger for all lead acid battery types - GEL/AGM/Wet.
- **USE:** Wheelchairs.
- **FEATURES:** Monitors charger output. Automatically shuts off when batteries are fully charged.
- **CORDS AVAIL:** 05 (Round/Square), 06 (Square), 08 (Round), 09 (Square).

**MK Model: M-L-13770 World Charger - Dual Mode/Dual Voltage**  
(Specify chair type or plug style) (Not Pictured)

- **DESCRIPTION:** 24 volt, 8 amp dual mode for both sealed type and wet lead acid batteries. Can be used with 115 or 230 AC voltage supply (switchable). AC plug not included.
- **USE:** Wheelchairs.
- **FEATURES:** Monitors charger output. Automatically shuts off when batteries are fully charged.
- **CORDS AVAIL:** 05 (Round/Square), 06 (Square), 08 (Round), 09 (Square).



**MK Model: LS 24/4**  
(Pictured Above)

- **DESCRIPTION:** 24 volt, 4 amp for sealed type or wet lead acid batteries. Can be used with 115 or 230AC voltage supply (switchable).
- **USE:** Scooters or power wheelchairs with lower capacity batteries such as U-1 size.
- **FEATURES:** Fully automatic with L.E.D.'s that signal charging and ready stages.
- **CORDS AVAIL:** 08 (Round Only).



**MK Model: M-L-12090**

(Specify chair type or plug style) (Pictured Above)


- **DESCRIPTION:** 24 volt, 8 amp, dual mode for both sealed type and wet lead acid batteries.
- **USE:** Wheelchairs.
- **FEATURES:** Monitors charger output. Automatically shuts off when batteries are fully charged.
- **CORDS AVAIL:** 05 (Round/Square), 06 (Square), 08 (Round), 09 (Square).



**MK Model: LS2606**

- **DESCRIPTION:** 12 volt, 6 amp for sealed type or wet lead acid batteries.
- **USE:** For use in single battery systems such as ventilators and other single 12 volt battery applications.
- **FEATURES:** Charge light and ready light, final stage float mode. Light weight and compact.
- **CORDS AVAIL:** Alligator clips only.

## Cord / Plug Reference Guide



MANUFACTURER	DESCRIPTION	NUMBER	STOCKED AT
E&J	White/Yellow Square Round/Silver	05 05/Round	All Warehouses Anaheim & Baltimore
INVACARE	Red/Black Square Round Silver	06 08	All Warehouses All Warehouses
QUICKIE/SUNRISE	Round/Silver Square (Black 6 Prong + Center)	08 09	All Warehouses Anaheim Only
SHOPRIDER	Round/Silver	08	All Warehouses

**Note:** All round plugs look identical. Make sure you know what type of chair prior to ordering. All chairs are wired differently, so it does make a difference.

