



MULTICHAIR 3000

The MULTICHAIR name says it all: An adjustable-height raised toilet seat, shower chair and tub transfer bench all-in-one. Weighing in at just 14 lbs., it's engineered to be light, strong, and portable, and built to last a lifetime, likely making it the first and last shower / commode chair you'll buy.

Features

- Seat height adjusts from 20-1/2" to 24-1/2" (easily fits over the highest toilets required by the ADA)
- Self-skinning padded urethane back and seat cushion with 1/2" waterfall-style foam ring around the commode opening and slot help protect the ischial and coccyx areas to reduce risk of pressure sores
- Cantilever fold-back / locking padded chair arms
- Overall depth just 19" to fit in compact areas
- Won't rust or corrode—constructed of aircraft-grade aluminum, brass, stainless steel and durable plastics

Specifications

Approximate weight: 14 pounds

- Weight capacity: 300 pounds
- Passed ANSI / BIFMA X5.1 test performed by Gaynes Labs, Inc. (a national testing laboratory). Weights of 250 lbs. and 300 lbs. were dropped on the chair, with no structural breakage or loss of serviceability

Dimensions and Notes

A = 35-1/2" to 39-1/2" (adjustable)

B = 15-1/2" (top edge of back cushion to top of cushion or cloth)

C = not used

D = 18" to 24" (adjustable, floor to top of seat cushion)

F = 14-1/2" (inside horizontal clearance between legs)

G = 17" (overall width at the crutch tips of legs)

