



MULTICHAIR 4000
MULTICHAIR 9000

*Safe and stable for home use, yet lightweight and portable for travel, the MULTICHAIR 4000 / 9000 series does it all. Its leg clamping design, unique to **nuprod**, provides an open and accessible seat frame, yet is as solid as a one-piece welded chair. The MULTICHAIR 9000 features a larger seat base for added room plus 20" rear wheels. Both easily break down to pack into the optional carrying case.*

Features

- Swing-away / removable push handles
- Cantilever fold-back / locking padded chair arms
- Adjustable swing-away footrests
- Self-skinning padded urethane back and seat cushion with 1/2" waterfall-style foam ring around the commode opening and slot help protect the ischial and coccyx areas to reduce risk of pressure sores
- Seat height adjusts from 20-1/2" to 24-1/2" (easily fits over the highest toilets required by the ADA)
- The 4000 features four 5" diameter casters with total-lock brakes, while the 9000 uses 20" wheels rear with brakes and 5" casters in front.
- The 9000, with a larger 20" x 20" seat cushion offers 111 additional sq. in. of seat space
- Won't rust or corrode—constructed of aircraft-grade aluminum, brass, stainless steel and durable plastics

Specifications

- Approximate weight: 28 pounds
- Weight capacity: 300 pounds
- Passed ANSI / BIFMA X5.1 test performed by Gaynes Labs, Inc. (a national testing laboratory). Weights of 250 lbs. and 300 lbs. were dropped on the chair, with no structural breakage or loss of serviceability

Dimensions and Notes

- A = 35-1/2" to 39-1/2" (adjustable)
- B = 15-1/2" (top edge of back cushion to top of cushion or cloth)
- C = 16-1/2" (inside minimum clearance)
- D = 20-1/2" to 24-1/2" (adjustable, floor to top of seat cushion)
- E = 12" to 20" (adjustable)
- Overall depth, 4000 = 32." Overall depth, 9000 = 39" (includes 6" push handles and 12" footrests)

