

INNOVATIVE PRODUCTS FOR SENIORS AND PEOPLE WITH DISABILITIES



MULTICHAIR 4000Tx
 (carrying case included)

Safe and stable for home use, the 4000Tx is lightweight and portable enough to be towed behind most standard wheel chairs to facilitate self-transfers for independent users. It conveniently packs in a black nylon carrying case that complies with FAA requirements for carry-on luggage, and is small enough to go through most airport x-ray machines.

Features

- Adjustable waterproof nylon cloth back with swing-away / removable push handles
- Cantilever fold-back / locking padded chair arms
- Self-skinning padded cushions with 1/2" waterfall-style foam ring around the commode opening and slot help protect the ischial and coccyx areas to reduce risk of pressure sores.
- Seat height adjusts from 20-1/2" to 24-1/2" (easily fits over the highest toilets required by the ADA)
- Adjustable swing-away footrests
- Four 4" diameter casters with total-lock brakes
- Won't rust or corrode—constructed of aircraft-grade aluminum, brass, stainless steel and durable plastics

Specifications

- Approximate weight: 28 lbs.
- Weight capacity: 300 lbs.
- Passed ANSI / BIFMA X5.1 test performed by Gaynes Labs, Inc. (a national testing laboratory). Weights of 250 lbs. and 300 lbs. were dropped on the chair, with no structural breakage or loss of serviceability

