

INNOVATIVE PRODUCTS FOR SENIORS AND PEOPLE WITH DISABILITIES AND



**MULTICHAIR 4020**  
**MULTICHAIR 4024**

*Perfect for the more independent users, the user-propelled MULTICHAIR 4020 and 4024 feature all of the **nuprod** hallmarks: Light, strong, portable, durable, and with easy adjustability, they are adaptable to nearly any living situation. Their modular design allows them to easily pack down into a compact package, making them an ideal choice for those who are not in the same living quarters every day.*

### Features

- Swing-away / removable push handles
- Cantilever fold-back / locking padded chair arms
- Self-skinning padded urethane cushions with 1/2" waterfall-style foam ring around the commode opening and slot help protect the ischial and coccyx areas to reduce risk of pressure sores.
- Seat height adjusts from 20-1/2" to 24-1/2" (easily fits over the highest toilets required by the ADA), and can be raised or lowered independent of wheel diameter
- Adjustable swing-away footrests
- Two 20" or 24" diameter wheels feature toggle wheel locks and vinyl-coated hand rims (20" wheels come with 19" diameter hand rims for easier pushing)
- Solid rubber tires, stainless steel quick-release axles and aluminum axle blocks for durability
- 5" diameter casters with total-lock brakes
- Won't rust or corrode—constructed of aircraft-grade aluminum, brass, stainless steel and durable plastics

### Specifications

- Approximate weight: 38 pounds
- Weight capacity: 300 pounds
- Passed ANSI / BIFMA X5.1 test performed by Gaynes Labs, Inc. (a national testing laboratory). Weights of 250 lbs. and 300 lbs. were dropped on the chair, with no structural breakage or loss of serviceability

### Dimensions and Notes

- A = 35-1/2" to 39-1/2" (adjustable)
- B = 15-1/2" (top edge of seat cushion to top of back cushion)
- C = 16-1/2" (inside minimum clearance)
- D = 20-1/2" to 24-1/2" (adjustable, floor to top of seat cushion)
- E = 12" to 20" (adjustable)
- Overall depth is 38"-42" (includes 6" push handles and 12" footrests)

