

INNOVATIVE PRODUCTS FOR SENIORS AND PEOPLE WITH DISABILITIES



MULTICHAIR 6200

Add more side-to-side and front-to-back seat room plus increased load capacity, and you have the MULTICHAIR 6200. The modular design enables an amazing range of adjustability, so it easily fits into nearly any bathroom. The smooth-operating connector bridge rides on corrosion-proof self-aligning Delrin® wheels, and easily removes for shower curtain closure and improved caregiver access.

Features

- Swing-away / removable push handles
- Cantilever fold-back / locking padded chair arms
- Seamless padded self-skinning urethane back cushion and 20" x 20" x 2" thick naugahyde seat cushion with 1/2" waterfall-style foam ring around the commode opening and slot help protect the ischial and coccyx areas to reduce risk of pressure sores
- Seat height adjusts from 20-1/2" to 24-1/2" (easily fits over the highest toilets required by the ADA)
- Adjustable swing-away footrests
- Four 5" casters with total-lock brakes
- (4) adjustable-length arms with suction cups help securely locate the in-tub chair base
- Constructed of materials that won't rust or corrode: Aircraft-grade aluminum, brass, stainless steel and plastic

Specifications

- Approximate weight: 52 pounds
- Weight capacity: 400 pounds
- Passed ANSI / BIFMA X5.1 test performed by Gaynes Labs, Inc. (a national testing laboratory). Weights of 250 lbs. and 300 lbs. were dropped on the chair, with no structural breakage or loss of serviceability

Dimensions and Notes

- A = 35-1/2" to 39-1/2" (adjustable)
- B = 15-1/2" (top edge of seat cushion to top of back cushion)
- C = 16-1/2" (inside minimum clearance)
- D = 20-1/2" to 24-1/2" (adjustable, floor to top of seat cushion)
- H = 12" to 20" (adjustable)
- 20" x 20" seat cushion, 8" x 10" opening with a 4" slot
- 8" x 17" back cushion
- Overall depth 39" (includes 6" push handles and 12" footrests)
- Measures 30" between center of roll-in section and center of tub section
- Different length connector bridges available

