



## MULTICHAIR 6200Tilt

*This larger-size Tilt-in-Space system utilizes a Center-of-Mass design to centralize the user's weight through the entire tilt range, allowing a more compact base to fit in small tubs other systems can't. The side-mounted tilt crank operates with low-effort, even with heavier users\*. Best-of-all, the 6200Tilt eliminates risky bathroom transfers, and is a cost-effective alternative to installing a roll-in shower.*

### Features

- Unique Tilt-in-Space system controlled by hand crank for easy caregiver / user adjustment of tilt angle of up to 30°. Crank lever can be located for left- and right-hand bathroom installations
- The seat cushion measures 20" x 20" for an additional 111 sq. in. of seating space
- Locking connector bridge self-aligns and is easily removed for shower curtain closure and improved caregiver access.
- Self-skinning padded urethane cushion with 1/2" waterfall-style foam ring around the commode opening / slot help protect the ischial / coccyx areas to reduce risk of pressure sores.
- Adjustable padded headrest
- Cantilever fold-back / locking padded chair arms
- Seat height adjusts from 22-1/2" to 26-1/2" (easily fits over the highest toilets required by the ADA).
- Swing-away / removable push handles
- Adjustable padded calf support and swing-away footrests
- Four 5" casters with total-lock brakes
- Won't rust or corrode—constructed of aircraft grade aluminum, brass, stainless steel and plastic.

### Specifications

- Approximate weight: 66 pounds
- Weight capacity: 400 pounds
- Passed ANSI/BIFMA X5.1 test performed by Gaynes Labs, Inc.(a national testing laboratory). Weights of 250 lbs. and 300 lbs.were dropped on the chair, with no structural breakage or loss of serviceability

### Dimensions and notes

- A = 37-1/2" to 40-1/2" (adjustable)
- B = 15-1/2" (top edge of seat cushion to top of back cushion)
- C = 16-1/2" (inside minimum clearance)
- D = 22-1/2" to 26-1/2" (adjustable, floor to top of seat cushion)
- E = 12" to 20" (adjustable)
- 46" overall depth fits the smallest bathroom
- 33" from center of roll-in section to center of tub section. Custom length bridges available to fit most any bathroom.

\* Up to rated capacity

