

permobil

THE POWER OF MOBILITY™



ADVENTURE AWAITS... INDOORS OR OUT!

The Permobil M300 PS Jr. combines the capabilities of our M300 base with the versatility of the PS Jr. seating system for the perfect pediatric chair. With a tight turning radius, you'll be able to maneuver through your home with ease, while our innovative CATLink (Climbing and Traction Linkage) System increases capability to get over obstacles and helps prevent high-centering. The expandable seat pan allows for a custom fit, so the chair can grow with the child. Plus, the PS Jr. is designed with convenient flip-back armrests, making transfers a breeze. And thanks to the UniTrack, standard on the PS Jr., you can have a variety of adjustments and accessories for added comfort.

M300
PS JR.



M300PSJR012213



CATLINK SUSPENSION
Casters are connected for added climbing stability.



AERON LEGRESTS
Adjustable, individual legrests, allow for better positioning.



TIGHT TURNING RADIUS
20" / 510 mm turning radius for increased maneuverability.




FLIP-BACK ARMRESTS
Armrests flip back completely out of the way for easier transfers.

M300 PS JR.

AVAILABLE SHROUD COLORS



DIGITAL ARMY CAMO



MULTICAM CAMO



REALTREE AP CAMO



ONYX BLACK



VOLCANO RED



SUBLIME GREEN



PLATINUM SILVER



COBALT BLUE



PURPLE RAIN



PHANTOM GREY



POP STAR PINK

SPECIFICATIONS

Maximum User Weight	165 lbs	75 kg
Maximum Speed Options	6 MPH	10 km/h
Driving Range ¹	13 miles	21 km
Base Width	24"	615 mm
Base Length	36.5"	925 mm
Minimum Turning Radius	20"	510 mm
Ground Clearance	3"	75 mm
Optional Seat Elevator	8"	200 mm
Battery Type	Group 24, Group 34	
Drive Electronics	120A R-net	
Maximum Slope ²	9° Up / 9° Down	

Seat to Floor Height	(Elevator)	17"	430 mm
	(Seat Tube)	17", 17.75", 18.5"	430, 450, 470 mm
Tilt Options	(Power)	0° - 45°	
	(Manual)	0° - 45°	
Recline Options	(Power)	No Power Option	
	(Manual)	85° - 125°	
Legrest Elevation Options	(Fixed)	70°	
	(Center Mount)	90° - 150°	
Armrest Pad Lengths	10", 13"	260, 335mm	
Seat Widths	11" - 16"	280 - 410 mm	
Seat Depths	10" - 18"	260 - 460 mm	
Backrest Heights	16.5" - 22.5"	420 - 570 mm	

¹ ACTUAL DRIVING RANGE WILL VARY BASED ON DRIVING AND BATTERY CONDITIONS.
² REFER TO OWNER MANUAL FOR SAFE DRIVING PRACTICES