


permobil

THE POWER OF MOBILITY™



Swing Away Mount

Allows tablet to be cleared from in front of user without the need for additional tools



Easily Configurable

Rotates to portrait or landscape orientation, while spring loaded arms can securely hold and easily release tablet.



Adjustable Supports

Can be used to provide additional support or can be removed to accommodate wider tablets.

DEVICE COMPATIBILITY

ON CHAIR - 1106972-99-0
& 1103476-99-0

SPARE PARTS - 1830709
& 1826839

DEVICE COMPATIBILITY

Supports most 10" tablets with or without a protective sleeve, that fit the following dimensions:

Minimum Width = 6.25"
Maximum Width = 8.1"
Depth = 0.875"
Maximum Length
(w/adjustable supports) = 10.25"

RAM® X-GRIP® 10" TABLET HOLDER

Your tablet can help you with information, communication, entertainment and more. But NOT if you can't get to it easily. Now you can have your favorite mobile device right in front of you, whenever you need it, with the RAM® X-Grip® 10" Tablet Holder. The X-Grip® can securely hold most 10" screen tablets with or without a protective sleeve. Using a durable, swing-away mounting system, the holder easily attaches to the UniTrack rail of your Permobil, providing the stability and convenience you need.

AVAILABLE ON ALL CORPUS 3G AND CORPUS HD CONFIGURATIONS.

*Appearance and specifications are subject to change without notice.

