



### Versatile Mounting

Available with fixed and adjustable mounting hardware to accommodate any setup.



### Quick and Easy to Use

Spring-loaded arms for securely holding and quickly removing a wide variety of smartphones with or without protective cases.



### Infinitely Adjustable

Ball joints can position the phone at the perfect angle for optimal convenience.

PART NUMBERS
ON CHAIR - 1106806-99-0
SPARE PARTS - 1830348

  

DEVICE COMPATIBILITY
Supports a wide variety of devices that fit the following dimensions:
Minimum Width = 1.875"
Maximum Width = 3.25"
Depth = 0.875"

## RAM® X-GRIP® PHONE HOLDER

You're on the go and you're connected... but why should you have to carry your mobile devices where you can't get to them easily? With the RAM® X-Grip® Phone Holder, connectivity is at your fingertips. The RAM® X-Grip® Phone Holder provides a secure mount for a wide variety of smartphones, with or without a protective case. The specialized mounting bracket allows the holder to be easily attached to the UniTrack rail on your Permobil and can be reconfigured to accommodate a variety of mounting positions.

AVAILABLE ON ALL CORPUS 3G AND CORPUS HD CONFIGURATIONS.

\*Appearance and specifications are subject to change without notice.

