

# SmartDrive

Featuring  
SwitchControl

permobil





# SwitchControl

SwitchControl gives users another programmable option for controlling their SmartDrive. With the push of a button SmartDrive can be engaged for a momentary burst of power or activated in latched mode for consistent power over

extended distances. This empowers a more diverse group of people to experience SmartDrive with increased confidence and when paired with the PushTracker E2 offers added flexibility, convenience and control options.

## SmartDrive Specifications

<b>Dimensions</b>	Length: 15.3" (389 mm) Width: 5.7" (141 mm) Height: 9.5" (242 mm)
<b>Product weight</b>	12.5 lbs (5.7 kg) Drive Unit
<b>User weight</b>	30 to 331 lbs (14 to 150 kg)
<b>Driving range*</b>	Up to 12.3 mi (19.8 km)
<b>Motor power</b>	250 W Brushless DC
<b>Maximum speed</b>	5.5 mph (8.8 km/h)
<b>Operating temp</b>	-13°F to 122°F (-25°C to 50°C)
<b>Chair type compatibility</b>	Rigid, Folding, Tilt-In Space, One-Arm Drive and others
<b>Drive wheel diameter</b>	22", 24", 25", 26" (501, 540, 559, 590 mm)
<b>Warranty</b>	2-Year Limited Warranty

<b>Bluetooth® LE</b>	FCC ID QOQBT113 IC ID 5123A-BGTBLE113
<b>HCPCS</b>	E0986

*\*Performance measurements will vary based on user and chair characteristics, driving and battery conditions.*

*The SmartDrive has been tested and conforms to all applicable requirements of ANSI/RESNA Standards for Wheelchairs - Volume 1 and Volume 2 and EN 12184.*

*NOTE: A wheelchair's specifications provided by its manufacturer may be slightly affected by the addition of the SmartDrive.*

*Please refer to your owner's manual for more detailed operation and care instructions.*

*Since U.S. Medicare coding and coverage criteria is subject to change, the provider should always confirm the HCPCS code and coverage criteria as part of the client assessment process. ©2019 Permobil. Permobil has a policy of continual product improvement and reserves the right to amend this document.*