



Order form

# M300 PS Jr



Dealer Name\*:

Account #\*:

City, State, Zip\*:

Phone:

Email:

ATP:

Rep Name:

P.O. Number:

Quote Number:

Client Reference\*:

Client Weight\*:

*\*Required field*

Special considerations that need to be addressed here (e.g., diagnosis):

## Chair Model and Wheels

Mid-wheel drive M300 power base includes: PS Jr seat frame, adjustable suspension, 4-point tie-down hardware, VoltPro charger, and R-net 120 amp electronics. Pre-selected items are required for this configuration.

### Choose Base

- 10844 **M300 Base PS Jr - M300 NPO** \$8,616.00 K0848  
*MWD M300 power base set up for no power seat functions, max speed: 6 mph, weight capacity: 165 lb.*
- 10856 **M300 Base PS Jr - M300 SPO** \$8,616.00 K0856  
*MWD M300 power base set up for one power seat function, max speed: 6 mph, weight capacity: 165 lb.*
- 10857 **M300 Base PS Jr - M300 MPO** \$8,616.00 K0861  
*MWD M300 power base set up for two power seat functions, max speed: 6 mph, weight capacity: 165 lb.*

### Choose Drive Tires

- FLATFREE **Flat Free Tires 3"W x 14"T** No charge
- PNEUMATIC **Pneumatic Tires 3"W x 14"T** No charge
- 11054 **Gray Wheel Package - MWD** No charge  
 6" casters  8" casters
- 11055 **Black Wheel Package - MWD** No charge  
 6" casters  8" casters
- 11056 **Blackout Wheel Package - MWD** \$371.00  
*Custom Option.*  
 6" casters  8" casters

## Chair Model and Wheels

### Choose Drive Tires

- 10782 **Tri-Spoke Split Rim w/ Snow Tire Pneumatic 3x8"** \$514.00  
*Not recommended for year-long use and must be ordered along with standard drive wheels. Does not include caster wheels.*

## Base Options

### Choose Base Color

- Standard Shroud Colors**
- |                                       |                                      |  |  |           |
|---------------------------------------|--------------------------------------|--|--|-----------|
| <input type="checkbox"/> Cobalt Blue  | <input type="checkbox"/> Purple Rain | <input type="checkbox"/> Onyx Black      | <input type="checkbox"/> Phantom Gray  | No charge |
| <input type="checkbox"/> Popstar Pink | <input type="checkbox"/> Volcano Red | <input type="checkbox"/> Platinum Silver | <input type="checkbox"/> Sublime Green |           |
- Standard Rim Colors**
- |   |                                      |  |  |           |
|---|--------------------------------------|--|--|-----------|
| <input type="checkbox"/> Cobalt Blue    | <input type="checkbox"/> Purple Rain | <input type="checkbox"/> Onyx Black      | <input type="checkbox"/> Phantom Gray  | No charge |
| <input type="checkbox"/> Popstar Pink   | <input type="checkbox"/> Volcano Red | <input type="checkbox"/> Platinum Silver | <input type="checkbox"/> Sublime Green |           |
| <input type="checkbox"/> Extreme Yellow |                                      |  |  |           |

### Batteries

- 107032-99-0 **Batteries, Grp 34 Gel - Pair** \$1,048.00 E2359

### Optional Items

- 105428-99-0 **Enhanced Steering Perform Unit (ESP), R-net** \$1,659.00  
*Required when ordering any ventilator tray. Recommended for switched input devices.*
- 105297-99-0 **LED Light Kit for M300** \$1,347.00
- 106223-99-0 **Stealth BPAP Vent Tray - MWD Base** \$752.00 E1029  
*12"W x 7"D surface that can support up to 40 lb.*
- 109620-99-0 **Programming Key for R-net - Permobil Version** \$575.00

## Joystick Options

Choose joystick type, handle and mount. **NOTE:** For additional input devices, please refer to the Alternative Drive Systems section.

### Joystick Type

- 105246-99-0 **R-net LED Remote Joystick - No Display** No charge  
*No display, no monojacks, limited to 5 profiles.*
- 108319-99-0 **Permobil Joystick Module (PJSM)** \$1,182.00  
*2.8" color display.*
- 105700-99-0 **R-net Remote Color Joystick w/ Lights, Heavy Duty Gimbal** \$1,346.00  
*Color LCD joystick, includes 3" knob with 3/4" diameter. Not compatible with alternative joystick handles.*

## Joystick Options

### Joystick Handle

<input type="checkbox"/>	10045	<b>Cone Shaped Handle - Standard</b>	No charge
<input type="checkbox"/>	100960-99-0	<b>Chin Cup Handle</b>	\$139.00 E2324
<input type="checkbox"/>	100961-99-0	<b>T Handle</b>	\$196.00 E2323
<input type="checkbox"/>	100962-99-0	<b>Large Ball</b>	\$139.00 E2323
<input type="checkbox"/>	100963-99-0	<b>Softball</b>	\$183.00 E2323
<input type="checkbox"/>	100964-99-0	<b>Mushroom Handle</b>	\$183.00 E2323
<input type="checkbox"/>	100965-99-0	<b>Stick Handle</b>	\$183.00 E2323
<input type="checkbox"/>	105442-99-0	<b>1.75" Black Foam Ball</b>	\$139.00 E2323
<input type="checkbox"/>	105846-99-0	<b>Bodypoint 3" U-Shaped Handle w/ Flex-Shaft</b>	\$183.00 E2323
<input type="checkbox"/>	105847-99-0	<b>Bodypoint 4" U-Shaped Handle w/ Flex-Shaft</b>	\$183.00 E2323
<input type="checkbox"/>	105848-99-0	<b>Bodypoint Rubber Dome Handle</b>	\$183.00 E2323

### Joystick Mount

<b>Fixed</b>	<input type="checkbox"/> Mount on Left	<input type="checkbox"/> Mount on Right	No charge
<b>Swing-Away</b>	<input type="checkbox"/> Mount on Left	<input type="checkbox"/> Mount on Right	\$409.00 E1028
<b>Slimline Retractable</b>	<input type="checkbox"/> Mount on Left	<input type="checkbox"/> Mount on Right	\$550.00 E1028
<input type="checkbox"/>	107960-99-0	<b>Adjustable Height Panel Bracket</b> <i>Provides 1.77" of vertical adjustment and 25° of angle adjustment.</i>	\$299.00 E1028
<input type="checkbox"/>	103472-99-0	<b>Mount Standard Joystick On Backrest, PS, PSJR, KO</b> <i>Can be chosen instead of attendant control.</i>	\$340.00
	<input type="checkbox"/> Mount on Left	<input type="checkbox"/> Mount on Right	
<input type="checkbox"/>	108268-99-0	<b>Bodypoint Tri-Lock Midline Mount Assembly</b> <i>Must select specific joystick attachment separately.</i>	\$702.00
	<input type="checkbox"/> Mount on Left	<input type="checkbox"/> Mount on Right	
<input type="checkbox"/>	108280-99-0	<b>Bodypoint Tri-Lock Mount for R-net Joystick</b>	\$292.00
<input type="checkbox"/>	108279-99-0	<b>Bodypoint Tri-Lock Mount for Compact Joystick</b>	\$292.00
<input type="checkbox"/>	10842	<b>PERMOFix Mount for mo-Vis/Compact Joystick</b> <i>PermoFix mount for attaching Compact and mo-Vis joysticks to Tri-Lock midline mount.</i>	\$292.00

### Other Options

<input type="checkbox"/>	10411	<b>Expandable Controller - R-net</b> <i>This item is required when ordering an Alternative Drive Control.</i>	\$764.00 E2377
<input type="checkbox"/>	10513	<b>Harness for Expandable Controller</b> <i>Must be ordered with any expandable controller.</i>	\$678.00 E2313
<input type="checkbox"/>	109779-99-0	<b>R-net Attendant Control</b>	\$1,499.00 E2331
<input type="checkbox"/>	103811-99-0	<b>Stroller Handle for PS Jr Seating System</b> <i>Angle adjustable.</i>	\$264.00
<input type="checkbox"/>	10582	<b>Output Module (IOM) for R-net Control System</b> <i>Enables user to control an external device (e.g., communication device, door opener, etc.) through the drive input (e.g., JSM, head array, etc.). May require specialty cable or wiring to be ordered separately.</i>	\$1,412.00

## Power Seating

### Adjustable Seat Height

- |                          |             |   |                  |
|--------------------------|-------------|---|------------------|
| <input type="checkbox"/> | 105116-99-0 | <b>Seat Tube Fixed</b>                        | No charge        |
| <input type="checkbox"/> | 106955-99-0 | <b>Power Adjustable Seat Height-8" Travel</b> | \$3,538.00 E2300 |

### Power Seat Functions

- |                          |             |                                   |                  |
|--------------------------|-------------|-----------------------------------|------------------|
| <input type="checkbox"/> | 103712-99-0 | <b>PS Jr Seat Frame - No Tilt</b> | No charge        |
| <input type="checkbox"/> | 103716-99-0 | <b>45° Power Tilt</b>             | \$7,297.00 E1002 |

### Seat Function Control Options | ICS switchbox availability may be limited to certain configurations.

- |                          |             |   |                  |
|--------------------------|-------------|---|------------------|
| <input type="checkbox"/> | 10418       | <b>Seat Function Control Switch for R-net</b><br><i>Does not include operation of power seat functions through drive control.</i>           | No charge        |
| <input type="checkbox"/> | 10419       | <b>Single Seat Function Control Kit for R-net</b><br><i>Allows operation of a single power seat function through drive control.</i>         | \$2,110.00 E2310 |
| <input type="checkbox"/> | 10585       | <b>Multiple Seat Function Control Kit for R-net</b><br><i>Allows operation of a two or more power seat functions through drive control.</i> | \$3,994.00 E2311 |
| <input type="checkbox"/> | 104135-99-0 | <b>Push Buttons on ICS Switchbox</b>  | No charge        |
| <input type="checkbox"/> | 104136-99-0 | <b>Toggle Switches on ICS Switchbox</b>   | \$357.00         |

### Switchbox Mounting

- |                          |             |  |           |
|--------------------------|-------------|--|-----------|
| <input type="checkbox"/> | 10626       | <b>Mount ICS Switchbox onto R-net Joystick Bracket</b>   | No charge |
| <input type="checkbox"/> | 104651-99-0 | <b>Mount ICS Switchbox onto Armrest</b><br><input type="checkbox"/> Mount on Left <input type="checkbox"/> Mount on Right  | No charge |
| <input type="checkbox"/> | 106256-99-0 | <b>Mount ICS Switchbox on Back</b><br><i>No charge when ordering Omni2 Control Module.</i><br><input type="checkbox"/> Mount on Left <input type="checkbox"/> Mount on Right | \$172.00  |
| <input type="checkbox"/> | 105409-99-0 | <b>Do Not Install ICS Switchbox, Send as Part</b>  | No charge |

### Power Seating Options

- |                          |             |  |                  |
|--------------------------|-------------|--|------------------|
| <input type="checkbox"/> | 104196-99-0 | <b>ICS Alternative Switchbox for R-net Electronics</b><br><i>Accommodates up to 8 momentary switches.</i>  | \$1,396.00       |
| <input type="checkbox"/> | 106743-99-0 | <b>Remote Stop System - All Electronics</b><br><i>Consists of a 315MHz radio frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair.</i> | \$1,381.00 E2321 |

## Backrest Options

### PS Jr Backrest Options

- |   |             |   |           |
|---|-------------|---|-----------|
| <input checked="" type="checkbox"/>   | 103717-99-0 | <b>PS Jr Back Canes, Height &amp; Angle Adjustable</b><br><i>Must select PS Jr style armrests with this option. Backrest width should match seat width.</i> | No charge |
| <b>PS Jr Solid Backrest</b> <input type="checkbox"/> 11" W <input type="checkbox"/> 12" W <input type="checkbox"/> 13" W <input type="checkbox"/> 14" W <input type="checkbox"/> 15" W <input type="checkbox"/> 16" W |             |   | No charge |
| <input type="checkbox"/>  | 10449       | <b>Omit Backrest</b>  | No charge |
| <input type="checkbox"/>  | 104616-99-0 | <b>Backrest Adjustment Quick Lock, PS / PS Jr</b><br><i>Provides manual back angle adjustment.</i>  | \$217.00  |

## Backrest Options

**Comfort Backrests** | *Backrest heights represent actual product dimensions. Must select tubular back canes.*  
**Available for quote and order March 2023.**

<b>Acta-Embrace Comfort-Tek</b>	<input type="checkbox"/> 10" W x 14" T	<input type="checkbox"/> 12" W x 14" T	<input type="checkbox"/> 14" W x 16" T	<input type="checkbox"/> 14" W x 18" T	\$642.00	E2613
<b>Acta-Back Comfort-Tek</b>	<input type="checkbox"/> 14" W x 16" T	<input type="checkbox"/> 14" W x 18" T	<input type="checkbox"/> 16" W x 16" T	<input type="checkbox"/> 16" W x 18" T	<input type="checkbox"/> 16" W x 20" T	\$618.00 E2613

### Comfort Backrest Options

<input type="checkbox"/> 109029-99-0	<b>BodiLink Premium 3" x 5" Lateral Supports - Pair</b> <i>Must also select hardware.</i>	\$446.00	E0956
<input type="checkbox"/> 109030-99-0	<b>BodiLink Premium 4" x 6" Lateral Supports - Pair</b> <i>Must also select hardware.</i>	\$446.00	E0956
<input type="checkbox"/> 109903-99-0	<b>BodiLink GT2 Swing-Away Lateral Hardware - Pair</b> <i>Must select with BodiLink lateral trunk supports on Comfort backrests.</i>	\$625.00	E1028

## Seat Sizing

**Seat Pan** | *Choose width and depth*

<b>Seat Width</b>	<input type="checkbox"/> 11" W	<input type="checkbox"/> 12" W	<input type="checkbox"/> 13" W	<input type="checkbox"/> 14" W	<input type="checkbox"/> 15" W	<input type="checkbox"/> 16" W	No charge		
<b>Seat Depth</b>	<input type="checkbox"/> 5" D	<input type="checkbox"/> 6" D	<input type="checkbox"/> 7" D	<input type="checkbox"/> 8" D	<input type="checkbox"/> 9" D	<input type="checkbox"/> 10" D	<input type="checkbox"/> 11" D	<input type="checkbox"/> 12" D	No charge
	<input type="checkbox"/> 13" D	<input type="checkbox"/> 14" D	<input type="checkbox"/> 15" D	<input type="checkbox"/> 16" D	<input type="checkbox"/> 17" D	<input type="checkbox"/> 18" D			

## Seat Cushion Options

### Mini-Flex Seat Cushions

<b>Leatherette</b>	<input type="checkbox"/> 10" W x 10" D	<input type="checkbox"/> 12" W x 12" D	<input type="checkbox"/> 14" W x 14" D	<input type="checkbox"/> 16" W x 16" D	\$653.00
<b>Black Mesh</b>	<input type="checkbox"/> 10" W x 10" D	<input type="checkbox"/> 12" W x 12" D	<input type="checkbox"/> 14" W x 14" D	<input type="checkbox"/> 16" W x 16" D	\$657.00

**Comfort Seat Cushions** | **Available for quote and order March 2023.**

<b>Curve Comfort-Tek</b>	<input type="checkbox"/> 10" W x 10" D	<input type="checkbox"/> 12" W x 12" D	<input type="checkbox"/> 14" W x 14" D	\$134.00	E2601	
<b>Embrace Anti-Thrust Comfort-Tek</b>	<input type="checkbox"/> 10" W x 10" D	<input type="checkbox"/> 12" W x 12" D	<input type="checkbox"/> 12" W x 14" D	<input type="checkbox"/> 14" W x 14" D	\$428.00	E2607
	<input type="checkbox"/> 14" W x 16" D					
<b>M2 Anti-Thrust</b>	<input type="checkbox"/> 14" W x 14" D	<input type="checkbox"/> 14" W x 16" D	<input type="checkbox"/> 16" W x 16" D	<input type="checkbox"/> 16" W x 18" D	\$474.00	E2607

### Cushion Options

<input type="checkbox"/> 105118-99-0	<b>Omit Seat Cushion</b>	No charge
--------------------------------------	--------------------------	-----------

## Armrest Options

### Armrest Set-Up

<input type="checkbox"/> 10825	<b>Flip Back, Height Adjustable Armrests - NC</b> <i>No charge with power tilt. Set Armrest Height to:</i>	No charge
	<input type="checkbox"/> 6" <input type="checkbox"/> 7" <input type="checkbox"/> 8" <input type="checkbox"/> 9" <input type="checkbox"/> 10" <input type="checkbox"/> 11" <input type="checkbox"/> 12" <input type="checkbox"/> 13" <input type="checkbox"/> 14"	
<input type="checkbox"/> 10824	<b>Flip Back, Height Adjustable Armrests</b> <i>Set Armrest Height to:</i>	\$426.00 E0973
	<input type="checkbox"/> 6" <input type="checkbox"/> 7" <input type="checkbox"/> 8" <input type="checkbox"/> 9" <input type="checkbox"/> 10" <input type="checkbox"/> 11" <input type="checkbox"/> 12" <input type="checkbox"/> 13" <input type="checkbox"/> 14"	

## Armrest Options

### Choose Armrest Pads

- |                                  |  |           |
|----------------------------------|--|-----------|
| <input type="checkbox"/> CL Pads | <b>Corpus Leatherette</b><br><input type="checkbox"/> 10" Left <input type="checkbox"/> 13" Left <input type="checkbox"/> 10" Right <input type="checkbox"/> 13" Right     | No charge |
| <input type="checkbox"/> CF Pads | <b>Corpus Black Fabric</b><br><input type="checkbox"/> 10" Left <input type="checkbox"/> 13" Left <input type="checkbox"/> 10" Right <input type="checkbox"/> 13" Right    | No charge |
| <input type="checkbox"/> UT Pads | <b>UniTrack Armrest Pad</b><br><input type="checkbox"/> 8" Left <input type="checkbox"/> 12.5" Left <input type="checkbox"/> 8" Right <input type="checkbox"/> 12.5" Right | No charge |
| <input type="checkbox"/> CG Pads | <b>Corpus Gel</b><br><input type="checkbox"/> 10" Left <input type="checkbox"/> 13" Left <input type="checkbox"/> 10" Right <input type="checkbox"/> 13" Right             | \$192.00  |

### Optional Items

- |                         |   |          |
|-------------------------|---|----------|
| <b>Armrest Pouch</b>    | <input type="checkbox"/> Left <input type="checkbox"/> Right  | \$385.00 |
| <b>Armrest Supports</b> | <input type="checkbox"/> Left Outer <input type="checkbox"/> Left Inner <input type="checkbox"/> Right Outer <input type="checkbox"/> Right Inner | \$385.00 |

## Legrest Options

### Legrest Type

- |                             |   |           |
|-----------------------------|---|-----------|
| <b>Swing-Away</b>           | <input type="checkbox"/> 70° w/ Angle Adjustable Footplate <input type="checkbox"/> 90° w/ Angle Adjustable Footplate | No charge |
| <b>Center Mount (Aeron)</b> | <input type="checkbox"/> 2pc Footplate <input type="checkbox"/> 1pc Footplate   | \$403.00  |

### Optional Items

- |                                      |   |           |
|--------------------------------------|---|-----------|
| <input type="checkbox"/> 104684-99-0 | <b>UT Calf Support Kit (B), Mesh Cover &amp; Adj Hardware</b><br><i>3.5"W x 4.5"T Pad. Cannot be used with swing-away legrests.</i> | No charge |
| <input type="checkbox"/> 104686-99-0 | <b>UT Calf Support Kit (J), Mesh Cover &amp; Adj Hardware</b><br><i>5"W x 6.5"T Pad. Cannot be used with swing-away legrests.</i>   | No charge |
| <input type="checkbox"/> 107503-99-0 | <b>BP Ankle Huggers, Med Side Release, Adj - Pair</b>   | \$315.00  |

## Accessories

### Belt Options

- |                                      |  |           |
|--------------------------------------|--|-----------|
| <input type="checkbox"/> 104518-99-0 | <b>Bodypoint Non-Padded Hip Belt - Medium</b>              | No charge |
| <input type="checkbox"/> 105849-99-0 | <b>Bodypoint Padded Hip Belt - Extra Small</b>             | \$167.00  |
| <input type="checkbox"/> 104520-99-0 | <b>Bodypoint Padded Hip Belt - Medium</b>                  | \$167.00  |
| <input type="checkbox"/> 100519-99-0 | <b>Positioning Belt Retractable</b>                        | \$332.00  |
| <input type="checkbox"/> 108267-99-0 | <b>Bodypoint Padded Hip Belt "Airplane" Latch - Medium</b> | \$224.00  |

### Head Support Options

- |                                      |   |          |       |
|--------------------------------------|---|----------|-------|
| <input type="checkbox"/> 100749-99-0 | <b>MF Headrest for Koala/PS Jr - Leatherette</b><br><i>Must also select hardware.</i> | \$374.00 | E0955 |
| <input type="checkbox"/> 10552       | <b>Adj/Removable Headrest Hardware</b>  | \$306.00 | E1028 |
| <input type="checkbox"/> 107499-99-0 | <b>Universal Headrest Adapter</b>   | \$302.00 |       |

## Accessories

### Leg Support Options

<input type="checkbox"/>	103077-99-0	<b>UT Hip Suppt (B), Mesh Cover - Low - Pair</b> <i>4.5"W x 3.5"T pad. Must also select hardware.</i>	\$315.00	E0953
<input type="checkbox"/>	103079-99-0	<b>UT Hip Suppt (A), Mesh Cover - Low - Pair</b> <i>6"W x 5.5"T pad. Must also select hardware.</i>	\$315.00	E0953
<input type="checkbox"/>	103095-99-0	<b>UT Thigh/Hip Suppt (C), Mesh Cover - High - Pair</b> <i>8"W x 3.5"T pad. Must also select hardware.</i>	\$315.00	E0953
<input type="checkbox"/>	10550	<b>Adj Removable - Thigh Support Hardware - Pair</b>	\$602.00	E1028

### Other Options

<input type="checkbox"/>	106210-99-0	<b>RAM Self-Leveling Cup Holder</b>	\$93.00	
<input type="checkbox"/>	106806-99-0	<b>RAM X-Grip Phone Holder</b>	\$93.00	
<input type="checkbox"/>	107100-99-0	<b>RAM X-Grip 10" Tablet Holder w/ Swing-Away Mount</b>	\$467.00	
<input type="checkbox"/>	107189-99-0	<b>USB Charger</b> <i>Provides up to 1.5A (at 5V) to charge most smartphones and tablets.</i>	\$231.00	
		<input type="checkbox"/> Mount on Left Armrest <input type="checkbox"/> Mount on Right Armrest		
<input type="checkbox"/>	105805-99-0	<b>Medical Necessities Bag Hooks</b>	\$111.00	E0950
<input type="checkbox"/>	102059-99-0	<b>Medical Necessities Bag</b>	\$321.00	
<input type="checkbox"/>	108276-99-0	<b>Bodypoint PivotFit Shoulder Harness - Medium</b>	\$289.00	
<input type="checkbox"/>	100578-99-0	<b>Small Angle Adjustable Tray</b> <i>Must also select hardware.</i>	\$339.00	
<input type="checkbox"/>	10556	<b>Adj Removable - Tray Hardware - 1 Each</b>	\$306.00	
		<input type="checkbox"/> Mount on Left <input type="checkbox"/> Mount on Right		
<input type="checkbox"/>	108231-99-0	<b>Permobil Adjustable Tray - 11"W</b> <i>Must also select hardware.</i>	\$284.00	E0960
<input type="checkbox"/>	108230-99-0	<b>Permobil Adjustable Tray - 13"W</b> <i>Must also select hardware.</i>	\$284.00	E0950
		<b>Permobil Adjustable Tray Hardware</b> <input type="checkbox"/> Left Side Adjustable <input type="checkbox"/> Right Side Adjustable	\$208.00	
<input type="checkbox"/>	103798-99-0	<b>After Market Seating System Adapter</b> <i>For UniTrack seat rail.</i>	\$109.00	

## Alternative Drive Systems

Omni2 Control Module required for all alternative drive systems. Enhanced Steering Performance (ESP) recommended for switched input devices.

### Alt Proportional Joysticks

<input type="checkbox"/>	109525-99-0	<b>mo-Vis Micro Joystick for Omni</b> <i>Operating Force = 8 g (0.28 oz). Neutral to Full Deflection = 0.13" (3.3 mm).</i> <b>Available for quote and order March 2023.</b>	\$3,200.00	E2312
<input type="checkbox"/>	109527-99-0	<b>mo-Vis Multi Joystick for Omni</b> <i>Operating Force = 50 g (1.8 oz). Neutral to Full Deflection = 0.31" (8 mm).</i> <b>Available for quote and order March 2023.</b>	\$3,200.00	E2312
<input type="checkbox"/>	103894-99-0	<b>Compact Joystick Lite for Omni - (9-pin connector)</b> <i>Operating Force = 50 g (1.8 oz). Neutral to Full Deflection = 0.7" (17mm).</i>	\$2,003.00	E2373
<input type="checkbox"/>	103893-99-0	<b>Compact Joystick for Omni - (9-pin connector)</b> <i>Operating Force = 220 g (7.8 oz). Neutral to Full Deflection = 1.1" (28mm).</i>	\$1,771.00	E2373

## Alternative Drive Systems

**Alt Proportional Joystick Mounts** | Please see *Alt Drive Accessories* section for additional mounting options.

<input type="checkbox"/>	104690-99-0	<b>Alt Joystick Swing-Away Mount Kit - Armrest</b> <input type="checkbox"/> Mount on Left <input type="checkbox"/> Mount on Right	\$310.00	E1028
<input type="checkbox"/>	107135-99-0	<b>Alt Joystick / Omni Combo Swing-Away Mount Kit - Armrest</b> <input type="checkbox"/> Mount on Left <input type="checkbox"/> Mount on Right	\$586.00	E1028
<input type="checkbox"/>	104143-99-0	<b>Alt Joystick Adjustable Mount Kit - Chin Boom</b> <i>Swing-Away / Removable Chin Boom mounted on right of backrest shell. Custom charge may be applied if ordering lower backrest heights.</i>	\$966.00	E1028
<input type="checkbox"/>	10752	<b>Chin Control Harness - Large</b> <i>For neck width up to 7½".</i>	\$687.00	E1028
<input type="checkbox"/>	10753	<b>Chin Control Harness - Small</b> <i>For neck width up to 5½".</i>	\$687.00	E1028

### Alt Proportional Joystick Accessories

<input type="checkbox"/>	104904-99-0	<b>Dual Switches for Compact Joystick - Fixed Mount</b> <input type="checkbox"/> Mount Left of JS <input type="checkbox"/> Mount Right of JS <input type="checkbox"/> Mount front of JS <input type="checkbox"/> Mount back of JS	\$288.00	
<input type="checkbox"/>	104641-99-0	<b>Dual Switches for Compact Joystick - Gooseneck Mount</b>	\$518.00	

### Total Control Head Array

<input type="checkbox"/>	109558-99-0	<b>Permobil Total Control Head Array</b> <i>Features an occipital pad sensor and two lateral arms configured with proximity sensors. Must select swing away hardware and mounting hardware below to complete assembly.</i>	\$4,819.00	E2330
<input type="checkbox"/>	11575	<b>Swing-Away Hardware for Side Sensor - Left</b>	\$263.00	E1028
<input type="checkbox"/>	11576	<b>Swing-Away Hardware for Side Sensor - Right</b>	\$263.00	E1028
<input type="checkbox"/>	11072	<b>Program Head Array as a 5-Switch System</b> <i>Select additional arms/sensors below and indicate the intended use of each sensor in the notes section.</i>	No charge	
<input type="checkbox"/>	109563-99-0	<b>Adj/Removable BodiLink HA Hardware - MAX 8</b>	\$365.00	E1028
<input type="checkbox"/>	109564-99-0	<b>Adj/Removable BodiLink HA Hardware - MAX 8 Flip Down</b> <i>Flip-down joint that allows head array to swing back 120°.</i>	\$501.00	E1028
<input type="checkbox"/>	11601	<b>(Optional) Add Left Lateral Arm w/ Proximity Sensor</b>	\$688.00	E1028
<input type="checkbox"/>	11602	<b>(Optional) Add Left Lateral Arm w/ Egg Switch</b>	\$563.00	E1028
<input type="checkbox"/>	11603	<b>(Optional) Add Left Lateral Arm w/ Pad (No Sensor or Switch)</b>	\$459.00	E1028
<input type="checkbox"/>	11604	<b>(Optional) Add Right Lateral Arm w/ Proximity Sensor</b>	\$688.00	E1028
<input type="checkbox"/>	11605	<b>(Optional) Add Right Lateral Arm w/ Egg Switch</b>	\$563.00	E1028
<input type="checkbox"/>	11606	<b>(Optional) Add Right Lateral Arm w/ Pad (No Sensor or Switch)</b>	\$459.00	E1028

### ASL FUSION Head Array

<input type="checkbox"/>	109582-99-0	<b>ASL FUSION Proportional and Digital Head Array</b> <i>Includes occipital (ASL 935) and lateral (ASL 933) pads, and pulse mechanical switch (ASL 320BK) for reset/mode changes. Must select hardware listed below to complete assembly.</i>	\$6,715.00	E2328
<input type="checkbox"/>	11580	<b>Swing-Away Hardware for Lateral Sensor - Left</b>	\$342.00	E1028
<input type="checkbox"/>	11581	<b>Swing-Away Hardware for Lateral Sensor - Right</b>	\$342.00	E1028
<input type="checkbox"/>	11582	<b>Swing-Away Hardware for Reset/Mode Switch - Right</b>	\$342.00	E1028
<input type="checkbox"/>	11583	<b>Adj/Removable Headrest Mounting Hardware</b> <i>Includes ASL 606 and BodiLink receiver.</i>	\$342.00	E1028
<input type="checkbox"/>	109583-99-0	<b>ASL FUSION Display, Programmer and Attendant Control w/ Mount</b> <i>Includes ASL 165 and quick release mount.</i>	\$730.00	K0108



## Alternative Drive Systems

### Sip & Puff Systems

<input type="checkbox"/>	10573	<b>Sip &amp; Puff System for Omni</b>	\$1,937.00	E2325, E2326
<input type="checkbox"/>	10574	<b>Sip &amp; Puff System with Fiber Optic Opti-Stop</b>	\$2,577.00	E2325, E2326
<input type="checkbox"/>	10762	<b>Adjustable Mounting Hardware - Sip &amp; Puff</b> <i>Required with any sip &amp; puff system.</i>	\$310.00	

## Alt Drive System Setup

### Omni2 Options | *Required for all alternative drive systems.*

<input type="checkbox"/>	108523-99-0	<b>Omni2 Control Module for R-net - PS Jr</b> <i>3.5" VGA high resolution display, USB charger port, built-in IR control and integrated Bluetooth technology. Must select mounting hardware below. Specify alternative input device being used:</i>	\$2,061.00	
<input type="checkbox"/>	103476-99-0	<b>Adjustable Swing-Away Hardware - Omni/Tablet</b> <input type="checkbox"/> Mount on Left <input type="checkbox"/> Mount on Right	\$310.00	

## Alt Drive Accessories

### PERMOfix Mounting Solutions

<input type="checkbox"/>	106655-99-0	<b>PERMOfix UniTrack Fixed Mount Kit</b>	\$312.00	
<input type="checkbox"/>	106656-99-0	<b>PERMOfix UniTrack Removable Mount Kit</b>	\$549.00	
<input type="checkbox"/>	106657-99-0	<b>PERMOfix UniTrack Swing-Away Mount Kit</b>	\$786.00	
<input type="checkbox"/>	105907-99-0	<b>Proximity Sensor + Swivel Assembly</b> <i>Swivel assembly allows sensor to be mounted with PERMOfix. Since Proximity Sensors require power, CANNOT be used as On/Off switches.</i>	\$397.00	E1028
<input type="checkbox"/>	105919-99-0	<b>Egg Switch + Swivel Mount Assembly</b> <i>Swivel assembly allows switch to be mounted with PERMOfix.</i>	\$268.00	
<input type="checkbox"/>	105912-99-0	<b>Solid Oval Pad + Swivel Assembly</b> <i>Swivel assembly allows switch to be mounted with PERMOfix.</i>	\$162.00	
<input type="checkbox"/>	106675-99-0	<b>Oval Plate + Swivel Assembly - No Switch</b> <i>Blank plate provides mounting option for Microlight, Buddy Button or other switch options. Swivel assembly allows pad to be mounted with PERMOfix.</i>	\$142.00	

### Switches and Cables

<input type="checkbox"/>	101769-99-0	<b>Microlight Switch   Qty:</b>	\$217.00	
<input type="checkbox"/>	101767-99-0	<b>Buddy Button   Qty:</b>	\$183.00	
<input type="checkbox"/>	106293-99-0	<b>Egg Switch Black   Qty:</b>	\$183.00	
<input type="checkbox"/>	101771-99-0	<b>Mini Cup Switch   Qty:</b>	\$183.00	
<input type="checkbox"/>	107500-99-0	<b>Proximity Switch   Qty:</b>	\$723.00	
<input type="checkbox"/>	105935-99-0	<b>3.5mm Ext. Cable for Sensor/Switch, 0.3m Long</b>	\$36.00	
<input type="checkbox"/>	105931-99-0	<b>3.5mm Ext. Cable for Sensor/Switch, 1m Long</b>	\$39.00	
<input type="checkbox"/>	105932-99-0	<b>3.5mm Ext. Cable for Sensor/Switch, 2m Long</b>	\$47.00	
<input type="checkbox"/>	105933-99-0	<b>3.5mm Ext. Cable for Sensor/Switch, 3m Long</b>	\$51.00	
<input type="checkbox"/>	105934-99-0	<b>3.5mm Y-splitter Cable for Sensor/Switch</b> <i>Allows two Proximity Sensors/Mechanical Switches to connect to the same direction jack.</i>	\$39.00	

## Additional Order Instructions (for Permobil):

## Notes & Comments (not for Permobil):

## Chair Order Policy:

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- The client's name or code and the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

## Order Acknowledgement:

I, \_\_\_\_\_, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: \_\_\_\_\_

Printed Name: \_\_\_\_\_

Title: \_\_\_\_\_

Date: \_\_\_\_\_

*NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty. The HCPCS codes provided are not intended to be billing or legal advice, rather our interpretation of the code definitions. Use of the codes does not ensure coverage or payment for the item. For coverage information, verify the policy of the appropriate payer.*