



stingray

Standard Specifications

MEASUREMENTS:

INCHES

	Size 1	Size 2
Seat Width*	11.75	13.75
Seat Depth	7-11.75	9.5-15.5
Seat Back Height*	17-22.5	21.25-26.5
Seat to Footplate*	5.5-11.75	5.5-14.5
Seat to Floor*	20.5	20.5
Seat-to-Back Angle	170	170
Seat Angle	45	45
Overall Height	42	46
Overall Length	29	29
Overall Width	23	23
Folded Length	35.5	35.5
Folded Width	23	23
Folded Height	19.25	19.25
Headrest Extension	5.5	5.5
Shoulder Strap Height	11.5	16
Push Handle Adjustments	35-41.5	35-41.5
Weight Capacity Regular	88 lbs	88 lbs
Weight of Chair	29.5 lbs	29.5 lbs