

# R82 Stingray - E1233

## USD Retail Price List/ Order Form



Account #: \_\_\_\_\_ Purchase Order No.: \_\_\_\_\_ Date: \_\_\_\_\_  
 Contact: \_\_\_\_\_ Mark For: \_\_\_\_\_

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email	Phone:

### US Configuration


**Effective: 3/1/2018 - Rev 5**

- Stingray Frame
- Stingray Seat
- Cushions
- Head Support
- Side Support
- Hip Support
- Pommel
- Vests and Belts
- Tray
- Transport Fittings



- A** (Back Height) \_\_\_\_\_
- B** (Arm Height) \_\_\_\_\_
- C** (Lower Leg Length) \_\_\_\_\_
- D** (Seat Depth) \_\_\_\_\_
- E** (Chest Width) \_\_\_\_\_
- F** (Seat Width) \_\_\_\_\_
- G** (Height: Seat to Top of Head) \_\_\_\_\_
- H** (Shoulder Width) \_\_\_\_\_
- Age \_\_\_\_\_
- Weight \_\_\_\_\_

### Measurements



	<b>Size 1 Inches</b>	<b>Size 2 Inches</b>
Effective Seat Width (A)	7-11 <sup>3</sup> / <sub>4</sub> "	9-13 <sup>3</sup> / <sub>4</sub> "
Seat Width Between Knee Supports	11"	12 <sup>3</sup> / <sub>4</sub> "
Seat Width Between Hip Supports (Accessories)	7", 8 <sup>1</sup> / <sub>2</sub> ", 10 <sup>1</sup> / <sub>4</sub> "	9", 10 <sup>1</sup> / <sub>2</sub> ", 12"
Effective Seat Depth (B)	7 - 11 <sup>3</sup> / <sub>4</sub> "	9 <sup>3</sup> / <sub>4</sub> - 15 <sup>1</sup> / <sub>2</sub> "
Back Height at Shoulders (C)	10 <sup>1</sup> / <sub>2</sub> - 16"	14 - 19 <sup>3</sup> / <sub>4</sub> "
Backrest Height (C1)	16 - 21 <sup>1</sup> / <sub>2</sub> "	19 <sup>1</sup> / <sub>2</sub> - 25 <sup>1</sup> / <sub>4</sub> "
Footrest to Seat Distance (D)	5 <sup>1</sup> / <sub>2</sub> - 17 <sup>1</sup> / <sub>2</sub> "	5 <sup>1</sup> / <sub>2</sub> - 17 <sup>1</sup> / <sub>2</sub> "
Seat Surface Height at Front Edge (E)	20 <sup>1</sup> / <sub>4</sub> "	20 <sup>1</sup> / <sub>4</sub> "
Back Width, Inside (F)	11 <sup>3</sup> / <sub>4</sub> "	14"
Back Width, Outside (F1)	17 <sup>1</sup> / <sub>2</sub> "	20 <sup>1</sup> / <sub>4</sub> "
Overall Width (G)	23"	23"
Length with Turnable Front Wheels (H1)	28 <sup>3</sup> / <sub>4</sub> "	28 <sup>3</sup> / <sub>4</sub> "
Length with Large Front Wheels	32 <sup>3</sup> / <sub>4</sub> "	32 <sup>3</sup> / <sub>4</sub> "
Overall Height, Push Brace 15° (I)	39", 40", 41 <sup>1</sup> / <sub>4</sub> "	39", 40", 41 <sup>1</sup> / <sub>4</sub> "
Foot Plate, Width x Depth	12" x 7 <sup>3</sup> / <sub>4</sub> "	12" x 7 <sup>3</sup> / <sub>4</sub> "
Folded Length	35"	35"
Folded Height	19"	19"
Total Mass, Frame	29 <sup>1</sup> / <sub>2</sub> lbs.	29 <sup>1</sup> / <sub>2</sub> lbs.
Total Mass, Seat	12 lbs.	14 <sup>1</sup> / <sub>2</sub> lbs.
Max Load	88 lbs.	88 lbs.
Back Recline / Seat Plane Angle	Back: 170° / Seat: 45°	Back: 170° / Seat: 45°
Approx. Minimum User Height	29 <sup>1</sup> / <sub>2</sub> "	35 <sup>1</sup> / <sub>2</sub> "
Approx. Maximum User Height	36"	47 <sup>1</sup> / <sub>2</sub> "


**Stingray Standard Features**



- Tilt-in-Space - 45°
- Back Recline - 170°
- Turnable Seat - 180°
- One-Size Frame
- Accepts Two Seat Sizes for Maximum Growth Potential
- Adjustable Seat Depth
- Adjustable Back Height
- Suspension


Description	Part No.	USD Price	HCPCS
<input type="checkbox"/> Solid Tires, Individual Wheel Locks	8724212	\$2,445	E1233
<input type="checkbox"/> Solid Tires, Individual Wheel Locks & Hand Brakes *	8724211	\$2,820	E1233
<input type="checkbox"/> Solid Tires, Central Wheel Locks	8724222	\$2,545	E1233
<input type="checkbox"/> Solid Tires, Central Wheel Locks & Hand Brakes	8724221	\$2,920	E1233
<input type="checkbox"/> Solid Tires, Individual Wheel Locks & Quick Release Rear Axels	8722212	\$2,695	E1233

*\*Custom order; Approx. lead time: 6-8 weeks. Non-returnable item*



Description	Size	Part No.	USD Price	HCPCS
<input type="checkbox"/> Seat and Back size 1	1	8711120		
<input type="checkbox"/> Seat size 1	1	X8711000	\$410	E2292
<input type="checkbox"/> Back size 1	1	X8712000	\$505	E2291
<input type="checkbox"/> Seat and Back size 2	2	8712220		
<input type="checkbox"/> Seat size 2	2	X8721000	\$410	E2292
<input type="checkbox"/> Back size 2	2	X8722000	\$505	E2291
<input checked="" type="checkbox"/> Fully Reclining Back (Required)	One Size	X8702000	\$405	E1014

*Adjustable seat depth and height adjustable back*



Description	Size	Part No.	USD Price	HCPCS
<b>Required, must select a color</b>				
<input type="checkbox"/> Blue = 1		<input type="checkbox"/> Red = 2		
<input type="checkbox"/> Pink = 5		<input type="checkbox"/> Grey = 4		
<input type="checkbox"/> Seat and Back Cushion (Pr.)	1	875051x		N/C
<input type="checkbox"/> Seat and Back Cushion (Pr.)	2	875052x		N/C

*\*Colors Green, Orange & Yellow available as a custom order; Approx. lead time: 6-8 weeks. Non-returnable item*



Description	Measurements	Size	Part No.	USD Price	HCPCS
<input type="checkbox"/> Footrest Rod, Elevating	8 1/4"	Short	X8711100	\$110	K0108
<input type="checkbox"/> Footrest Rod, Elevating	11"	Long	X8712200	\$110	K0108










<input type="checkbox"/> Angle Adjustable Footplate		One Size	X8710010	\$95	K0040
---	--	----------	----------	------	-------

Description	Size	Part No.	USD Price	HCPCS
-------------	------	----------	-----------	-------



<input type="checkbox"/> Transit Option	One Size	8720902	\$200	
---	----------	---------	-------	--

High-Low Frames	Description	Size	Part No.	USD Price	HPCPS
	<input type="checkbox"/> High-Low:v Indoor Frame	1	8911871-20	\$1,770	
	<input type="checkbox"/> High-Low:v Indoor Frame	1, 2	8911872-20	\$2,193	
	<input type="checkbox"/> Push Brace	One Size	8910545	\$160	
	<input type="checkbox"/> Push Brace, Angle Adjustable	One Size	8910555	\$385	

Accessories	Description	Measurements	Size	Part No.	USD Price	HPCPS
	<input type="checkbox"/> Hip Supports, Black (Pr.)	6½" x 7"	1	8740230	\$225	E0956
	<input type="checkbox"/> Hip Supports, Black (Pr.)	7" x 8½"	2	8740240	\$225	E0956
	<input type="checkbox"/> Hip Supports, Pad, Black (Pr.)		1	8740250	\$225	E0956
	<input type="checkbox"/> Multi-Adjustable, Removable Hardware (Pr.) <b>(Required)</b>		1	x8740250	\$225	E1028
	<input type="checkbox"/> Hip Supports, Pad, Black (Pr.)		2	8740260	\$225	E0956
	<input type="checkbox"/> Multi-Adjustable, Removable Hardware (Pr.) <b>(Required)</b>		2	x8740260	\$225	E1028
<b>*Hip Support and Multi-Adjustable Hardware must be ordered together</b>						
	<input type="checkbox"/> Side Supports, Multi Adjustable (Pr.)	4" x 4"	1	8720301	\$195	E0956
	<input type="checkbox"/> Side Supports, Multi Adjustable (Pr.)	5" x 5¾"	2	8720302	\$195	E0956
	<input type="checkbox"/> Covers for Knobs on Side Supports (Pr.) *		One Size	8720360	\$30	
<b>*Custom order; Approx. lead time: 6-8 weeks. Non-returnable item</b>						
	<input type="checkbox"/> Medial Thigh Support (One-Piece) 3½" x 1½" x 2¾"		1	8720501	\$135	E0957
	<input type="checkbox"/> Adjustable, Removable Hardware <b>(Required)</b>		1	x8720501	\$165	E1028
	<input type="checkbox"/> Medial Thigh Support (One-Piece) 3½" x 3½" x 3½"		2	8720502	\$135	E0957
	<input type="checkbox"/> Adjustable, Removable Hardware <b>(Required)</b>		2	x8720502	\$165	E1028
	<input type="checkbox"/> Medial Thigh Support (One-Piece) 3½" x 3½" x 4¾"		3	8720503	\$135	E0957
	<input type="checkbox"/> Adjustable, Removable Hardware <b>(Required)</b>		3	x8720503	\$165	E1028
<b>*Medial Thigh Support and Adjustable, Removable Hardware must be ordered together</b>						
	<input type="checkbox"/> Anti-Thrust Wedge, Black (Pr.)	4" x 6"	1, 2, 3	8910405-3	\$46	K0108

Description	Measurements	Size	Part No.	USD Price	HPCPS
<input type="checkbox"/> Contoured Headrest, Black		1	8740810	\$200	E0955
<input type="checkbox"/> Contoured Headrest, Black		2	8740820	\$200	E0955
<input type="checkbox"/> Sculpted Head Support	9½" x 4½"	Small	8720831	\$205	E0955
<input type="checkbox"/> Sculpted Head Support	11" x 5¾"	Large	8720832	\$205	E0955
<b>Preferred option for Stingray seating</b>					
<input type="checkbox"/> Anatomic Headrest *		Small	8720862	\$255	E0955
<input type="checkbox"/> Anatomic Headrest *		Large	8720863	\$255	E0955
<b>*Requires Turtle Bar System Short or Long</b>					
<b>Suggested User Shoulder Height of 16" or Greater</b>					
<b>Swing-away Mounting Bracket - Height, Depth, Angle and Lateral Adjustable</b>					
<input type="checkbox"/> Turtle Bar System, Short	5"	Short	98619-25	\$230	E1028
<input type="checkbox"/> Turtle Bar System, Long	10"	Long	98619-26	\$230	E1028
<b>*Required with Anatomic Headrest</b>					
<input type="checkbox"/> 2-Part Headrest		One Size	8720870	\$270	E0955




Description	Measurements	Size	Part No.	USD Price	HCPCS
<input type="checkbox"/> Headrest Fittings		One Size	8720850	\$51	K0108
<input type="checkbox"/> Universal Headrest Mount *		One Size	87590	\$61	E1399

**\*Required when using other manufactures Head Support**

**Other types of headrest available, please see Head Supports Order Form**

Description	Measurements	Size	Part No.	USD Price	HCPCS
<input type="checkbox"/> Cross Vest	4½" x ½" x 11"	0	89625-060	\$190	K0108
<input type="checkbox"/> Cross Vest	5" x ½" x 12"	1	89625-160	\$190	K0108
<input type="checkbox"/> Cross Vest	5.2" x ½" x 12¾"	2	89625-260	\$190	K0108
<input type="checkbox"/> Cross Vest	6¼" x 4¾" x 15¼"	3	89625-360	\$190	K0108
<input type="checkbox"/> Cross Vest	7½" x 6¾" x 7"	4	89625-460	\$190	K0108
<input type="checkbox"/> Vest	1¼" x 2½" x 5½" x 3½" x 11"	0	89111-060	\$200	K0108
<input type="checkbox"/> Vest	1¼" x 2½" x 5¾" x 4¼" x 12"	¾	89111-07560	\$200	K0108
<input type="checkbox"/> Vest	2¼" x 2½" x 7" x 5" x 13¼"	1	89111-160	\$200	K0108
<input type="checkbox"/> Vest	2½" x 3½" x 9¼" x 5½" x 15¼"	2	89111-260	\$200	K0108
<input type="checkbox"/> Vest	2¾" x 4¾" x 11 x 5¾" x 18¼"	3	89111-360	\$200	K0108
<input type="checkbox"/> Vest	2¾" x 6¾" x 13¾" x 7½" x 19"	4	89111-460	\$200	K0108



<input type="checkbox"/> Neoprene Vest	7" x 5½" x 4" x 4¾"	1	89113-1	\$200	K0108
<input type="checkbox"/> Neoprene Vest	6½" x 9" x 6¾" x 5¾"	2	89113-2	\$200	K0108
<input type="checkbox"/> Neoprene Vest	11½" x 7¾" x 7" x 6¾"	3	89113-3	\$200	K0108



<input type="checkbox"/> 4-Point Pelvic Positioner, Black	10¼"-11¼" x 18" x 19½"	1	85427-1	\$135	E0978
<input type="checkbox"/> 4-Point Pelvic Positioner, Black	12"-13½" x 25" x 19½"	2	85427-2	\$135	E0978



<input type="checkbox"/> 2-Point Pelvic Belt, Black	10¼"-11¼" x 2½" x 18"	1	85424-1	\$97	E0978
---	-----------------------	---	---------	------	-------



<input type="checkbox"/> Foot Straps, Plastic Lock (Pr.)	23"	One Size	85096	\$56	K0040
<input type="checkbox"/> Foot Straps, Velcro (Pr.)	21"	One Size	85102	\$51	K0040

**Other types of vests and belts available, please see Vest and Belts Order Form**

Description	Measurements	Size	Part No.	USD Price	HCPCS
<input type="checkbox"/> Folding Top (Canopy), Black		1	8740110	\$330	
<input type="checkbox"/> Folding Top (Canopy), Black		2	8740120	\$330	
<input type="checkbox"/> Foam Back Plate *	5" x 12"	1	8776325	\$23	
<input type="checkbox"/> Foam Back Plate *	5½" x 13¾"	2	8776326	\$23	



**Future Growth - Fills out the gap between seat and the back**  
**\*Custom order; Approx. lead time: 6-8 weeks. Non-returnable item**



<input type="checkbox"/> Tray with Edge, Grey *	17.3" x 13.3"	One Size	8776350	\$100	E0950
<input type="checkbox"/> Adjustable Removable Tray Hardware *		One Size	x8776350	\$165	E1028

**\*Tray and Adjustable Removable Hardware must be ordered together**


<input type="checkbox"/> Rain Cover *		1	8720141	\$56	
<input type="checkbox"/> Rain Cover *		2	8720142	\$56	

**\*Requires Folding Top (Canopy)**



<input type="checkbox"/> Handle Grip	14"	1	8720401	\$200	
<input type="checkbox"/> Handle Grip	16"	2	8720402	\$200	

Description	Measurements	Size	Part No.	USD Price	HCPCS
<input type="checkbox"/> Backpack, Black		Small	87043183-1	\$56	
<i>*Custom order; Approx. lead time: 6-8 weeks. Non-returnable item</i>					
<input type="checkbox"/> Medical Essential Carrier		One Size	8720600	\$97	
<input type="checkbox"/> Vent Tray *	13" x 17"	One Size	8776290	\$460	E1029
<i>*Requires Central Wheel Locks</i>					
<input type="checkbox"/> O2 Tank Holder *			99999	\$350	E2208
<i>*Cannot be used during transit and manufacturer installation is required</i>					
<i>*Custom order; Approx. lead time: 6-8 weeks</i>					
<input type="checkbox"/> IV Pole, Height Adjustable and Collapsible		One Size	99999	\$200	K0105
<i>*Custom order; Approx. lead time: 6-8 weeks</i>					



Spare Parts	Description	Measurements	Size	Part No.	USD Price	HCPCS
-------------	-------------	--------------	------	----------	-----------	-------

**Required, must select a color**

**x = Color**

Blue = 1

Red = 2

Pink = 5

Grey = 4



<input type="checkbox"/> Seat and Back Cushion (Pr.)		1	875051x	\$443	
<input type="checkbox"/> Seat and Back Cushion (Pr.)		2	875052x	\$443	
<i>*Colors Green, Orange &amp; Yellow available as a custom order; Approx. lead time: 6-8 weeks. Non-returnable item</i>					



<input type="checkbox"/> Adductor Covers, Black (Pr.)		1	8750410	\$56	
<input type="checkbox"/> Adductor Covers, Black (Pr.)		2	8750420	\$56	



<input type="checkbox"/> Hardware Kit for Hip Supports, Fixed		One Size	872023-KIT	\$35	
---	--	----------	------------	------	--



<input type="checkbox"/> Release Strap for Back Angling		1	87102250	\$45.50	
<input type="checkbox"/> Release Strap for Back Angling		2	87102255	\$45.50	



<input type="checkbox"/> Footrest Rod, Elevating	8¼"	Short	8775317	\$110	
<input type="checkbox"/> Footrest Rod, Elevating	11"	Long	8775318	\$110	



<input type="checkbox"/> Angle Adjustable Footplate		One Size	8775010	\$260	
---	--	----------	---------	-------	--

*When being transported in a vehicle, it is R82's position that the wheelchair user transfer from the wheelchair into appropriate vehicle seating during transportation. Please note that at this time, the Department of Transportation has not officially approved any tie-down system for a user while in a wheelchair; therefore R82 cannot recommend a transit system.*

*The R82 Stingray mobility base was tested with a footrest and R82 transit anchors. The performance and safety of accessories and/or other seating options not tested by R82 is unknown, and no warranty or claim about their performance is or can be made.*