

R82 Wombat Living

USD Retail Price List/ Order Form



Account #: _____ Purchase Order No.: _____ Date: _____
 Contact: _____ Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email	Phone:

US Configuration



Effective: 04.25.19 - Rev 6

- Wombat Living
- Head Support
- Side Support
- Armrest
- Hip Support
- Pommel
- Vests and Belts
- Tray
- Helper's Handle



- A** (Back Height) _____
- B** (Arm Height) _____
- C** (Lower Leg Length) _____
- D** (Seat Depth) _____
- E** (Chest Width) _____
- F** (Seat Width) _____
- G** (Height from Seat to Top of Head) _____
- H** (Shoulder Width) _____
- Age _____
- Weight _____

Measurements



	Size 1 Inches	Size 2 Inches	Size 3 Inches
Effective seat width (A)	11¾"	14½"	16½"
Seat width between hip supports	4¼-11¾"	6¾-14½"	8½-16½"
Seat width, pelvis support w/sides	6¾-9½"	9-11¾"	11½-14¼"
Width, pelvis support without sides	9½"	9½"	11½"
Effective seat depth (B)	7-12¼"	9¾-15"	12½-19¾"
Backrest height (C)	11-13½"	12¼-15"	15¼-19¾"
Backrest height (C) (high back)	14¼-17"	15½-18"	18-22¾"
Footrest to seat distance (D)	6-10"	6-13¾"	6-17¼"
Seat surface height at front edge (4" Castors)	15¼-26¾"	15¼-26¾"	16½-27½"
Back width, flat cushion (F)	9¾"	11¾"	13¾"
Back width, contoured cushion, between side supports (F)	6¾"	9"	10½"
Arm rest height above seat cushion long rod and pad	2½-7"	2½-7"	2½-7"
Backrest angle (a1)	90° to 130°	90° to 130°	90° to 133°
Seat plane angle (a2)	-15° to 20°	-15° to 20°	-10° to 22°
Leg to seat surface angle (a3)	-30° to 30°	-30° to 30°	-30° to 30°
Overall width (G)	21½"	21½"	24"
Length (H)	23½"	23½"	28¼"
Max load	99.2 lbs.	99.2 lbs.	154.3 lbs.
Total mass	34.1 lbs.	34.1 lbs.	48½ lbs.

NOTES:

Wombat Living Standard Features



- Tilt in Space
- Back Recline and Prone Angle
- Adjustable Seat Depth
- Adjustable Back Height
- Height and Angle Adjustable Foot Support
- Without Seat/Back Cushions

Description	Pricing	
		95XXXX-00
Wombat Living, size 1	\$3,125	1
Wombat Living, size 2	\$3,125	2
Wombat Living, size 3	\$3,125	3
Wombat Living, size 1 with High Back••	\$3,125	5
Wombat Living, size 2 with High Back••	\$3,125	6
Wombat Living, size 3 with High Back••	\$3,125	7
•• Custom Order; Approx. lead time: 6-8 weeks		
Frame Color: Anthracite (Dark Grey)	\$0	2
4" Front/Rear Wheels	\$0	2
4" Front/Rear Wheels with central brake - <i>upcharge</i> ••	\$235	3
•• Custom Order; Approx. lead time: 6-8 weeks		
w/Gas Spring: Height Pedal on Frame, Tilt Lever Under Seat	\$0	1
w/Gas Spring: Height Pedal on Frame, Tilt Lever on Push Brace - <i>upcharge</i>	\$265	2
w/power actuator: Powered height setting, tilt lever under seat - <i>upcharge</i> ••	\$1,420	5
w/power actuator: Powered height setting, tilt lever on push brace - <i>upcharge</i> ••	\$1,420	6
••Custom Order for frame size 1 & 2; Apporx. lead time: 6-8 weeks		

Important!

Not recommended for users with a lot of involuntary movements

Accessories	Description	Size	Part No.	USD Price
-------------	-------------	------	----------	-----------

Required, must select a color

xx = Color

- Red = 24 Petroleum = 25 Grey = 28



<input type="checkbox"/> Back Cushion, Standard w/Cover (Ea.)	1	955031-2xx	\$140
<input type="checkbox"/> Back Cushion, Standard w/Cover (Ea.)	2	955032-2xx	\$140
<input type="checkbox"/> Back Cushion, Standard w/Cover (Ea.)	3	955033-2xx	\$140



<input type="checkbox"/> Back Cushion, Contoured w/Grey Cover (Ea.)	1	955031-238	\$240
<input type="checkbox"/> Back Cushion, Contoured w/Grey Cover (Ea.)	2	955032-238	\$240
<input type="checkbox"/> Back Cushion, Contoured w/Grey Cover (Ea.)	3	955033-238	\$240

Description	Size	Part No.	USD Price
-------------	------	----------	-----------



<input type="checkbox"/> Seat Cushion, Standard w/Cover (Ea.)	1	955021-2xx	\$270
<input type="checkbox"/> Seat Cushion, Standard w/Cover (Ea.)	2	955022-2xx	\$270
<input type="checkbox"/> Seat Cushion, Standard w/Cover (Ea.)	3	955023-2xx	\$270

Accessories	Description	Size	Part No.	USD Price
-------------	-------------	------	----------	-----------



<input type="checkbox"/> Seat Cushion, Contoured w/Grey Cover (Ea.)	1	955021-238	\$270
<input type="checkbox"/> Seat Cushion, Contoured w/Grey Cover (Ea.)	2	955022-238	\$270
<input type="checkbox"/> Seat Cushion, Contoured w/Grey Cover (Ea.)	3	955023-238	\$270



<input type="checkbox"/> Seat Cushion, Contoured w/ Incontinence Cover (Ea.)	1	955021-237	\$290
<input type="checkbox"/> Seat Cushion, Contoured w/ Incontinence Cover (Ea.)	2	955022-237	\$290
<input type="checkbox"/> Seat Cushion, Contoured w/ Incontinence Cover (Ea.)	3	955023-237	\$290

Description	Measurements	Size	Part No.	USD Price
-------------	--------------	------	----------	-----------



- Lower Back Ext, Black (Ea.) * H:2½" x W:10½" 3 95140-60 \$97
*Only available on Wombat Living sz 3



- Lumbar Foam Wedge (Ea.) 1 955030-1 \$25
- Lumbar Foam Wedge (Ea.) 2 955030-2 \$25
- Lumbar Foam Wedge (Ea.) 3 955030-3 \$25

Description	Measurements	Size	Part No.	USD Price
-------------	--------------	------	----------	-----------




- Hip Support (Pr.) H:3¼" x W:7½" 1,2 95110-1 \$200
- Covers for Hip Supports (Required) 1,2 95110-60 \$0
- Hip Support (Pr.) H:4" x W:10" 3 95111-1 \$200
- Covers for Hip Supports (Required) 3 95111-60 \$0



- Pelvis Support without Sides (Ea.) ♦ 1,2 95140-61 \$231
- Pelvis Support without Sides (Ea.) ♦ 3 95140-63 \$231
- ♦ Only available on Wombat Living with High Back



- Pelvis Support with Sides (Ea.) ♦ 1 95140-64 \$417
- Pelvis Support with Sides (Ea.) ♦ 2 95140-65 \$417
- Pelvis Support with Sides (Ea.) ♦ 3 95140-66 \$417
- ♦ Only available on Wombat Living with High Back



- Pommel (Ea.) H:3½" x W:3½" x D:3½" 1 95130-1 \$145
- Pommel (Ea.) H:3½" x W:3½" x D:4½" 2 95130-2 \$145
- Pommel (Ea.) H:3½" x W:3½" x D:6" 3 95130-3 \$145
- Pommel (Ea.) H:3½" x W:3½" x D:7" 4 95130-4 \$145



- Knee Abduction Pommel One Size 95230-10 \$170



- Padding for Knee Abduction Pommel One Size 95430-1026C \$65

Description	Measurements	Size	Part No.	USD Price
-------------	--------------	------	----------	-----------



- Armrest, Short (Pr.) 1,2 95112-1 \$350
- Armrest, Short (Pr.) 3 95114 \$350



- Armrest, Long Rod and Pad (Pr.) 1,2 95113-1 \$350
- Armrest, Long Rod and Pad (Pr.) 3 95115 \$350

Accessories	Description	Size	Part No.	USD Price
-------------	-------------	------	----------	-----------



- Fixed Side Supports w/Covers (Pr.) 1,2 95121 \$170

- Swing-Away Lateral Support Bracket, RH (Ea.) One Size 83400-R \$225
- Swing-Away Lateral Support Bracket, LH (Ea.) One Size 83400-L \$225

- Swing-Away Lateral Support Bracket, Offset, RH (Ea.) One Size 83600-R \$225
- Swing-Away Lateral Support Bracket, Offset, LH (Ea.) One Size 83600-L \$225



- Bar, Straight (Ea.) 4¾" One Size 81932 \$0
- Bar, Straight, Stainless (Ea.) 3¾" One Size 82569 \$0




- Bar, Offset, 1" (Ea.) 4¾" One Size 81968 \$0
- Bar, Offset, 1¾" (Ea.) 4¾" One Size 82000 \$0
- Bar, Offset, 2½" (Ea.) 4¾" One Size 81969 \$0




<input type="checkbox"/> Curved Cushion (Ea.)	3½" x 7½"	One Size	6175406	\$105
<input type="checkbox"/> Curved Cushion (Ea.)	2½" x 6¾"	One Size	6175405	\$105
<input type="checkbox"/> Curved Cushion (Ea.)	2½" x 6"	One Size	6175404	\$105
<input type="checkbox"/> Standard Cushion (Ea.)	6" x 4¾"	One Size	6175403	\$105
<input type="checkbox"/> Standard Cushion (Ea.)	4" x 4"	One Size	6175402	\$105
<input type="checkbox"/> Standard Cushion (Ea.)	2⅔" x 4½"	One Size	6175401	\$105

Description	Measurements	Size	Part No.	USD Price
-------------	--------------	------	----------	-----------


■ XX = 60 (Black)




<input type="checkbox"/> Neoprene Vest (Ea.)	A:7 x B:5½ x C:4 x D:4½"	1	89113-1	\$200
<input type="checkbox"/> Neoprene Vest (Ea.)	A:10½ x B:11 x C:6½ x D:5¾"	2	89113-2	\$200
<input type="checkbox"/> Neoprene Vest (Ea.)	A:11¼ x B:9 x C:7 x D:6½"	3	89113-3	\$200




<input type="checkbox"/> Cross Vest (Ea.)	A: 5 x B:½ x C:12"	1	89625-1XX	\$190
<input type="checkbox"/> Cross Vest (Ea.)	A: 5½ x B:½ x C:13"	2	89625-2XX	\$190
<input type="checkbox"/> Cross Vest (Ea.)	A: 6¼ x B:2½ x C:15¼"	3	89625-3XX	\$190



<input type="checkbox"/> Chest Belt, Upholstered (Ea.)	A:11¼xB:7½xC:9¼-12"	1	85428-1	\$270
<input type="checkbox"/> Chest Belt, Upholstered (Ea.)	A:14xB:10xC:10-12¾"	2	85428-2	\$270
<input type="checkbox"/> Chest Belt, Upholstered (Ea.)	A:16½xB:12½xC:10-13"	3	85428-3	\$270



<input type="checkbox"/> 4-Point Padded Hip Belt	L:10¼-11½"	1	85427-1	\$135
<input type="checkbox"/> 4-Point Padded Hip Belt	L: 12½-13½"	2	85427-2	\$135



<input type="checkbox"/> 2-Point Pelvic Belt, Black	L: 13-14½"	1	85424-1	\$97
---	------------	---	---------	------




<input type="checkbox"/> Fix-Locks, Black, Bag of 6 Pcs **		One Size	862109-35	\$42
--	--	----------	-----------	------

****Required for use of Vests**

<input type="checkbox"/> Y-Hip Strap, Upholstered, 1" End Girth (Ea.)		1	85427-1	\$135
<input type="checkbox"/> Y-Hip Strap, Upholstered, 1" End Girth (Ea.)		2	85427-2	\$135

Accessories	Description	Measurements	Size	Part No.	USD Price
-------------	-------------	--------------	------	----------	-----------



<input type="checkbox"/> Upholstered Ankle Strap (Pr.)	A: 5½-7 x B:8½-10"	1	85426-1	\$185
<input type="checkbox"/> Upholstered Ankle Strap (Pr.)	A: 6¼-7¾ x B:10-11¾"	2	85426-2	\$185
<input type="checkbox"/> Upholstered Ankle Strap (Pr.)	A: 7-8½ x B:11-12½"	3	85426-3	\$185



<input type="checkbox"/> Neoprene Sandals (Pr.)		One Size	85696-1	\$130
---	--	----------	---------	-------


Other types of vests and belts available, please see Vest and Belts Order Form

Description	Measurements	Size	Part No.	USD Price
-------------	--------------	------	----------	-----------

Swing-away Mounting Brackets - Height, Depth, Angle and Lateral Adjustable


<input type="checkbox"/> Turtle Bar System ♦	A: 5", B: 4-8½"	Short	98619-25	\$230
<input type="checkbox"/> Turtle Bar System ♦	A: 10", B: 4-8½"	Long	98619-26	\$230

♦ **Required with Anatomic Head Support; 5099001, 5099002 or 5099003**



<input type="checkbox"/> Anatomic Head Support ♦	A:5½ x B:5½ x C:9½"	Small	5099001	\$255
<input type="checkbox"/> Head Support Cover (Required)		Small	9923-83	\$0
<input type="checkbox"/> Anatomic Head Support ♦	A:6 x B:5½ x C:10"	Medium	5099002	\$255
<input type="checkbox"/> Head Support Cover (Required)		Medium	9945-83	\$0
<input type="checkbox"/> Anatomic Head Support ♦	A:6½ x B:5½ x C:12½"	Large	5099003	\$255
<input type="checkbox"/> Head Support Cover (Required)		Large	9967-83	\$0

♦ **Turtle Bar System Short or Long Required; 98619-25 or 98619-26**



XX = Turtle Bar Size

92 = Short Turtle Bar

93 = Long Turtle Bar



9971xxx-xx
9973xxx-xx

Adjustable Headrest without Dynamic Feature - Requires Bracket (see below)

- Adj. Width Headrest, Large Laterals, Black Cover ♦ 99737xx-85 \$305
 - Adj. Width Headrest, Large Laterals, w/Speakers, Black Cover ♦♦ 99717xx-85 \$405
- ♦♦Custom order; Approx. lead time: 6-8 weeks.
♦ Turtle Bar System Short or Long Required; x98619-25 or x98619-26



9975xxx-xx

Adjustable Headrest with Dynamic Feature - Requires Bracket (see below)

- Adj. Width Headrest, Large Laterals, Black Cover ♦ 99736xx-85 \$385
 - Adj. Width Headrest, Large Laterals, w/Speakers, Black Cover ♦♦ 99756xx-85 \$385
- ♦♦Custom order; Approx. lead time: 6-8 weeks.
♦ Turtle Bar System Short or Long Required; x98619-25 or x98619-26

Headband System for Anatomic Head Support Only One Size 99500 \$230

Universal Headrest Mount ♦ One Size 87590 \$61

♦ Required when using other manufactures Head Support

Allows most Whitmyer, Stealth, Otto Bock headrest, I2I to mount to our products (x:panda, Serval, etc.). Please contact Customer Service for confirmation

Other types of headrest available, please see Head Supports Order Form

Description	Measurements	Size	Part No.	USD Price
<input type="checkbox"/> Directional Stabilizer Kit ♦♦		One Size	89592-KIT	\$215
♦♦Custom order; Approx. lead time: 6-8 weeks.				
<input type="checkbox"/> Tray (Ea.)		1,2	95125	\$300
<input type="checkbox"/> Tray (Ea.)		3	95126	\$300



Accessories

Description	Measurements	Size	Part No.	USD Price
-------------	--------------	------	----------	-----------



Tray, Grey 1,2 8976350 \$225



Push Brace, Seat Mounted 1,2 95148-1 \$230



Push Brace, Seat Mounted 3 95148-3 \$230

Angle Adjustable Push Brace (Ea.) 3 95147 \$210



Helper's Handle, Fold Down (Pr.) One Size 95145 \$80



Helper's Handle, Fixed (Pr.) One Size 95146 \$80

Anti-Tip Device One Size 95953 \$160

Spare Parts

Description	Measurements	Size	Part No.	USD Price
-------------	--------------	------	----------	-----------



Wombat Living Step 1 to 2 Kit, Seat & Back 1 95560 \$280

Wombat Living Step 1 to 2 Kit, Seat & Back 2 95561 \$280



Wombat Living Short to High Back Kit 1 95560-1 \$157

Wombat Living Short to High Back Kit 2 95561-1 \$157

Description	Size	Part No.	USD Price
-------------	------	----------	-----------

Required, must select a color

xx = Color

Red = 24

Petroleum = 25

Grey = 28

<input type="checkbox"/> Cover for Standard Back Cushion	1	955611-2xx	\$109
<input type="checkbox"/> Cover for Standard Back Cushion	2	955612-2xx	\$109
<input type="checkbox"/> Cover for Standard Back Cushion	3	955613-2xx	\$109
<input type="checkbox"/> Cover for Contoured Back Cushion, Grey	1	955611-238	\$137
<input type="checkbox"/> Cover for Contoured Back Cushion, Grey	2	955612-238	\$137
<input type="checkbox"/> Cover for Contoured Back Cushion, Grey	3	955613-238	\$137
<input type="checkbox"/> Cover for Seat Cushions	1	955621-2xx	\$109
<input type="checkbox"/> Cover for Seat Cushions	2	955622-2xx	\$109
<input type="checkbox"/> Cover for Seat Cushions	3	955623-2xx	\$109

Description	Size	Part No.	USD Price
-------------	------	----------	-----------

<input type="checkbox"/> Cover w/Foam for Pelvis Support without Sides	1,2	95122-1	\$130
<input type="checkbox"/> Cover w/Foam for Pelvis Support without Sides	3	95122-2	\$130
<input type="checkbox"/> Cover w/Foam for Pelvis Support with Sides	1	95123-1	\$154
<input type="checkbox"/> Cover w/Foam for Pelvis Support with Sides	2	95123-2	\$154
<input type="checkbox"/> Cover w/Foam for Pelvis Support with Sides	3	95123-3	\$154
<input type="checkbox"/> Covers for Fixed Side Supports (95121)	1,2	95121-60	\$84

Spare Parts	Description	Measurements	Size	Part No.	USD Price
-------------	-------------	--------------	------	----------	-----------



<input type="checkbox"/> Footplate, Angle & Height Adj.	12" x 7½"	95145-1	\$115
<input type="checkbox"/> Footplate, Angle & Height Adj.	14" x 9"	85145-2	\$115



<input type="checkbox"/> Footplate, Angle & Height Adj. ♦	14" x 9"	95145-2	\$115
<input type="checkbox"/> Footplate, Angle & Height Adj. ♦	16" x 9½"	95145-3	\$115
<input type="checkbox"/> Footplate, Angle & Height Adj. ♦	18¾" x 9½"	95145-4	\$115
♦ Includes 4" mounting fitting			



<input type="checkbox"/> Footplate with Rod	7½" x 12"	1	95144-1	\$225
<input type="checkbox"/> Footplate with Rod	9" x 14"	2	95144-2	\$225
<input type="checkbox"/> Footplate with Rod	9¼" x 16"	3	95144-3	\$225
<input type="checkbox"/> Footplate with Rod	9½" x 19"	4	95144-4	\$225

Description	Size	Part No.	USD Price
-------------	------	----------	-----------

<input type="checkbox"/> Batteries for Wombat Living with Power (Set)	One Size	2222012	\$305
♦ Both batteries must be charged at the same time			



<input type="checkbox"/> Charger & Plug for USA/Canada (Ea.)	One Size	2222161	\$320
--	----------	---------	-------

Prices in USD & is subject to change without notice