

R82 Wombat Living

USD Retail Price List/ Order Form

Account #: _____ Purchase Order No.: _____ Date: _____
 Contact: _____ Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email:	Phone:

US Configuration


Effective: 04.27.20 - Rev 7

- Wombat Living
- Head Support
- Side Support
- Armrest
- Hip Support
- Pommel
- Vests and Belts
- Tray
- Helper's Handle




- A (Back Height) _____
- B (Arm Height) _____
- C (Lower Leg Length) _____
- D (Seat Depth) _____
- E (Chest Width) _____
- F (Seat Width) _____
- G (Height from Seat to Top of Head) _____
- H (Shoulder Width) _____
- Age _____
- Weight _____

Measurements



Size 1
Inches



Size 2
Inches

Size 3
Inches

Effective seat width (A)	11 ³ / ₄ "	14 ¹ / ₂ "	16 ¹ / ₂ "
Seat width between hip supports	4 ¹ / ₄ -11 ³ / ₄ "	6 ³ / ₄ -14 ¹ / ₂ "	8 ¹ / ₂ -16 ¹ / ₂ "
Seat width, pelvis support w/sides	6 ³ / ₄ -9 ¹ / ₂ "	9-11 ³ / ₄ "	11 ¹ / ₂ -14 ¹ / ₄ "
Width, pelvis support without sides	9 ¹ / ₂ "	9 ¹ / ₂ "	11 ¹ / ₂ "
Effective seat depth (B)	7-12 ¹ / ₄ "	9 ³ / ₄ -15"	12 ¹ / ₂ -19 ³ / ₄ "
Backrest height (C)	11-13 ¹ / ₂ "	12 ¹ / ₄ -15"	15 ¹ / ₄ -19 ³ / ₄ "
Backrest height (C) (high back)	14 ¹ / ₄ -17"	15 ¹ / ₂ -18"	18-22 ³ / ₄ "
Footrest to seat distance (D)	6-10"	6-13 ³ / ₄ "	6-17 ¹ / ₄ "
Seat surface height at front edge (4" Castors)	15 ¹ / ₄ -26 ³ / ₄ "	15 ¹ / ₄ -26 ³ / ₄ "	16 ¹ / ₂ -27 ¹ / ₂ "
Back width, flat cushion (F)	9 ³ / ₄ "	11 ³ / ₄ "	13 ³ / ₄ "
Back width, contoured cushion, between side supports (F)	6 ³ / ₄ "	9"	10 ¹ / ₂ "
Arm rest height above seat cushion long rod and pad	2 ¹ / ₂ -7"	2 ¹ / ₂ -7"	2 ¹ / ₂ -7"
Backrest angle (a1)	90° to 130°	90° to 130°	90° to 133°
Seat plane angle (a2)	-15° to 20°	-15° to 20°	-10° to 22°
Leg to seat surface angle (a3)	-30° to 30°	-30° to 30°	-30° to 30°
Overall width (G)	21 ¹ / ₂ "	21 ¹ / ₂ "	24"
Length (H)	23 ¹ / ₂ "	23 ¹ / ₂ "	28 ¹ / ₄ "
Max load	99.2 lbs.	99.2 lbs.	154.3 lbs.
Total mass	34.1 lbs.	34.1 lbs.	48 ¹ / ₂ lbs.

NOTES:

Wombat Living Standard Features



- Tilt in Space
- Back Recline and Prone Angle
- Adjustable Seat Depth
- Adjustable Back Height
- Height and Angle Adjustable Foot Support
- Without Seat/Back Cushions

Description	Pricing	
		95XXXX-00
<input type="checkbox"/> Wombat Living, size 1	\$3,125	1
<input type="checkbox"/> Wombat Living, size 2	\$3,125	2
<input type="checkbox"/> Wombat Living, size 3	\$3,125	3
<input type="checkbox"/> Wombat Living, size 1 with High Back●●	\$3,125	5
<input type="checkbox"/> Wombat Living, size 2 with High Back●●	\$3,125	6
<input type="checkbox"/> Wombat Living, size 3 with High Back●●	\$3,125	7
●●Custom Order; Approx. lead time: 6-8 weeks		
<input type="checkbox"/> Frame Color: Anthracite (Dark Grey)	\$0	2
<input type="checkbox"/> 4" Front/Rear Wheels	\$0	2
<input type="checkbox"/> 4" Front/Rear Wheels with central brake - <i>upcharge</i> ●●	\$235	3
●●Custom Order; Approx. lead time: 6-8 weeks		
<input type="checkbox"/> w/Gas Spring: Height Pedal on Frame, Tilt Lever Under Seat	\$0	1
<input type="checkbox"/> w/Gas Spring: Height Pedal on Frame, Tilt Lever on Push Brace - <i>upcharge</i>	\$265	2
<input type="checkbox"/> w/power actuator: Powered height setting, tilt lever under seat - <i>upcharge</i> ●●	\$1,420	5
<input type="checkbox"/> w/power actuator: Powered height setting, tilt lever on push brace - <i>upcharge</i> ●●	\$1,420	6
●●Custom Order for frame size 1 & 2; Approx. lead time: 6-8 weeks		

Important!


Not recommended for users with a lot of involuntary movements

Accessories	Description	Size	Part No.	USD Price
-------------	-------------	------	----------	-----------

Required, must select a color

xx = Color

- Red = 24 Petroleum = 25 Grey = 28



<input type="checkbox"/> Back Cushion, Standard w/Cover (Ea.)	1	955031-2xx	\$140
<input type="checkbox"/> Back Cushion, Standard w/Cover (Ea.)	2	955032-2xx	\$140
<input type="checkbox"/> Back Cushion, Standard w/Cover (Ea.)	3	955033-2xx	\$140



<input type="checkbox"/> Back Cushion, Contoured w/Grey Cover (Ea.)	1	955031-238	\$240
<input type="checkbox"/> Back Cushion, Contoured w/Grey Cover (Ea.)	2	955032-238	\$240
<input type="checkbox"/> Back Cushion, Contoured w/Grey Cover (Ea.)	3	955033-238	\$240



Description	Size	Part No.	USD Price
<input type="checkbox"/> Seat Cushion, Standard w/Cover (Ea.)	1	955021-2xx	\$270
<input type="checkbox"/> Seat Cushion, Standard w/Cover (Ea.)	2	955022-2xx	\$270
<input type="checkbox"/> Seat Cushion, Standard w/Cover (Ea.)	3	955023-2xx	\$270


Accessories	Description	Size	Part No.	USD Price
-------------	-------------	------	----------	-----------



<input type="checkbox"/> Seat Cushion, Contoured w/Grey Cover (Ea.)	1	955021-238	\$270
<input type="checkbox"/> Seat Cushion, Contoured w/Grey Cover (Ea.)	2	955022-238	\$270



Seat Cushion, Contoured w/Grey Cover (Ea.) 3 955023-238 \$270



Seat Cushion, Contoured w/ Incontinence Cover (Ea.) 1 955021-237 \$290

Seat Cushion, Contoured w/ Incontinence Cover (Ea.) 2 955022-237 \$290

Seat Cushion, Contoured w/ Incontinence Cover (Ea.) 3 955023-237 \$290

Description	Measurements	Size	Part No.	USD Price
-------------	--------------	------	----------	-----------



Lower Back Ext, Black (Ea.) * H:2½" x W:10½" 3 95140-60 \$97

*Only available on Wombat Living sz 3



Lumbar Foam Wedge (Ea.) 1 955030-1 \$25

Lumbar Foam Wedge (Ea.) 2 955030-2 \$25

Lumbar Foam Wedge (Ea.) 3 955030-3 \$25

Description	Measurements	Size	Part No.	USD Price
-------------	--------------	------	----------	-----------



Hip Support (Pr.) H:3¼" x W:7½" 1,2 95110-1 \$200

Covers for Hip Supports (Required) 1,2 95110-60 \$0

Hip Support (Pr.) H:4" x W:10" 3 95111-1 \$200


Covers for Hip Supports (Required) 3 95111-60 \$0



Pelvis Support without Sides (Ea.) ♦ 1,2 95140-61 \$231

Pelvis Support without Sides (Ea.) ♦ 3 95140-63 \$231

♦ Only available on Wombat Living with High Back




Pelvis Support with Sides (Ea.) ♦ 1 95140-64 \$417

Pelvis Support with Sides (Ea.) ♦ 2 95140-65 \$417

Pelvis Support with Sides (Ea.) ♦ 3 95140-66 \$417

♦ Only available on Wombat Living with High Back



Pommel (Ea.) H:3½" x W:3½" x D:3½" 1 95130-1 \$145

Pommel (Ea.) H:3½" x W:3½" x D:4½" 2 95130-2 \$145

Pommel (Ea.) H:3½" x W:3½" x D:6" 3 95130-3 \$145

Pommel (Ea.) H:3½" x W:3½" x D:7" 4 95130-4 \$145



Knee Abduction Pommel One Size 95230-10 \$170



Padding for Knee Abduction Pommel One Size 95430-10260 \$65

Description	Measurements	Size	Part No.	USD Price
-------------	--------------	------	----------	-----------



Armrest, Short (Pr.) 1,2 95112-1 \$350

Armrest, Short (Pr.) 3 95114 \$350



Armrest, Long Rod and Pad (Pr.) 1,2 95113-1 \$350




Armrest, Long Rod and Pad (Pr.) 3 95115 \$350

Accessories	Description	Size	Part No.	USD Price
-------------	-------------	------	----------	-----------




Fixed Side Supports w/Covers (Pr.) 1,2 95121 \$170

Swing-Away Lateral Support Bracket, RH (Ea.) One Size 83400-R \$225


<input type="checkbox"/> Swing-Away Lateral Support Bracket, LH (Ea.)	One Size	83400-L	\$225	
<input type="checkbox"/> Swing-Away Lateral Support Bracket, Offset, RH (Ea.)	One Size	83600-R	\$225	
<input type="checkbox"/> Swing-Away Lateral Support Bracket, Offset, LH (Ea.)	One Size	83600-L	\$225	
 <input type="checkbox"/> Bar, Straight (Ea.)	4 ³ / ₄ "	One Size	81932	\$0
<input type="checkbox"/> Bar, Straight, Stainless (Ea.)	3 ³ / ₄ "	One Size	82569	\$0
 <input type="checkbox"/> Bar, Offset, 1" (Ea.)	4 ³ / ₄ "	One Size	81968	\$0
<input type="checkbox"/> Bar, Offset, 1 ³ / ₄ " (Ea.)	4 ³ / ₄ "	One Size	82000	\$78
<input type="checkbox"/> Bar, Offset, 2 ¹ / ₂ " (Ea.)	4 ³ / ₄ "	One Size	81969	\$0
 <input type="checkbox"/> Curved Cushion (Ea.)	3 ¹ / ₂ " x 7 ¹ / ₂ "	One Size	6175406	\$105
<input type="checkbox"/> Curved Cushion (Ea.)	2 ¹ / ₂ " x 6 ³ / ₄ "	One Size	6175405	\$105
<input type="checkbox"/> Curved Cushion (Ea.)	2 ¹ / ₂ " x 6"	One Size	6175404	\$105
<input type="checkbox"/> Standard Cushion (Ea.)	6" x 4 ³ / ₄ "	One Size	6175403	\$105
<input type="checkbox"/> Standard Cushion (Ea.)	4" x 4"	One Size	6175402	\$105
<input type="checkbox"/> Standard Cushion (Ea.)	2 ² / ₃ " x 4 ¹ / ₂ "	One Size	6175401	\$105

Description	Measurements	Size	Part No.	USD Price
-------------	--------------	------	----------	-----------


■ XX = 60 (Black)




<input type="checkbox"/> Neoprene Vest (Ea.)	A:7 x B:5 ¹ / ₂ x C:4 x D:4 ¹ / ₂ "	1	89113-1	\$200
<input type="checkbox"/> Neoprene Vest (Ea.)	A:10 ¹ / ₂ x B:11 x C:6 ¹ / ₂ x D:5 ³ / ₄ "	2	89113-2	\$200
<input type="checkbox"/> Neoprene Vest (Ea.)	A:11 ¹ / ₄ x B:9 x C:7 x D:6 ¹ / ₂ "	3	89113-3	\$200




<input type="checkbox"/> Cross Vest (Ea.)	A: 5 x B:1 ¹ / ₂ x C:12"	1	89625-1XX	\$190
<input type="checkbox"/> Cross Vest (Ea.)	A: 5 ¹ / ₂ x B:1 ¹ / ₂ x C:13"	2	89625-2XX	\$190
<input type="checkbox"/> Cross Vest (Ea.)	A: 6 ¹ / ₄ x B:2 ¹ / ₂ x C:15 ¹ / ₄ "	3	89625-3XX	\$190



<input type="checkbox"/> Chest Belt, Upholstered (Ea.)	A:11 ¹ / ₄ x B:7 ¹ / ₂ x C:9 ¹ / ₄ -12"	1	85428-1	\$270
<input type="checkbox"/> Chest Belt, Upholstered (Ea.)	A:14xB:10xC:10-12 ³ / ₄ "	2	85428-2	\$270
<input type="checkbox"/> Chest Belt, Upholstered (Ea.)	A:16 ¹ / ₂ x B:12 ¹ / ₂ x C:10-13"	3	85428-3	\$270



<input type="checkbox"/> 4-Point Padded Hip Belt	L:10 ¹ / ₄ -11 ¹ / ₂ "	1	85427-1	\$135
<input type="checkbox"/> 4-Point Padded Hip Belt	L: 12 ¹ / ₂ -13 ¹ / ₂ "	2	85427-2	\$135



<input type="checkbox"/> 2-Point Pelvic Belt, Black	L: 13-14 ¹ / ₂ "	1	85424-1	\$97
---	--	---	---------	------




<input type="checkbox"/> Fix-Locks, Black, Bag of 6 Pcs **		One Size	862109-35	\$42
--	--	----------	-----------	------

****Required for use of Vests**

<input type="checkbox"/> Y-Hip Strap, Upholstered, 1" End Girth (Ea.)		1	85427-1	\$135
---	--	---	---------	-------

<input type="checkbox"/> Y-Hip Strap, Upholstered, 1" End Girth (Ea.)		2	85427-2	\$135
---	--	---	---------	-------

Accessories	Description	Measurements	Size	Part No.	USD Price
-------------	-------------	--------------	------	----------	-----------



<input type="checkbox"/> Upholstered Ankle Strap (Pr.)	A: 5 ¹ / ₂ -7 x B:8 ¹ / ₂ -10"	1	85426-1	\$185
<input type="checkbox"/> Upholstered Ankle Strap (Pr.)	A: 6 ¹ / ₄ -7 ³ / ₄ x B:10-11 ³ / ₄ "	2	85426-2	\$185
<input type="checkbox"/> Upholstered Ankle Strap (Pr.)	A: 7-8 ¹ / ₂ x B:11-12 ¹ / ₂ "	3	85426-3	\$185





Neoprene Sandals (Pr.)

One Size

85696-1


\$130

Other types of vests and belts available, please see Vest and Belts Order Form


Description	Measurements	Size	Part No.	USD Price
-------------	--------------	------	----------	-----------

Swing-away Mounting Brackets - Height, Depth, Angle and Lateral Adjustable

- Turtle Bar System ♦ A: 5", B: 4-8½" Short 98619-25 \$230
 - Turtle Bar System ♦ A: 10", B: 4-8½" Long 98619-26 \$230
- ♦Required with Anatomic Head Support; 5099001, 5099002 or 5099003



- Anatomic Head Support ♦ A:5½ x B:5½ x C:9½" Small 5099001 \$255
 - Head Support Cover (Required) Small 9923-83 \$0
 - Anatomic Head Support ♦ A:6 x B:5½ x C:10" Medium 5099002 \$255
 - Head Support Cover (Required) Medium 9945-83 \$0
 - Anatomic Head Support ♦ A:6½ x B:5½ x C:12½" Large 5099003 \$255
 - Head Support Cover (Required) Large 9967-83 \$0
- ♦Turtle Bar System Short or Long Required; 98619-25 or 98619-26



XX = Turtle Bar Size

92 = Short Turtle Bar

93 = Long Turtle Bar

Adjustable Headrest without Dynamic Feature - Requires Bracket (see below)

- Adj. Width Headrest, Large Laterals, Black Cover ♦ 99737xx-85 \$305
 - Adj. Width Headrest, Large Laterals, w/Speakers, Black Cover ♦♦ 99717xx-85 \$405
- ♦♦Custom order; Approx. lead time: 6-8 weeks.
♦Turtle Bar System Short or Long Required; x98619-25 or x98619-26



9971xxx-xx

9973xxx-xx

Adjustable Headrest with Dynamic Feature - Requires Bracket (see below)

- Adj. Width Headrest, Large Laterals, Black Cover ♦ 99736xx-85 \$385
 - Adj. Width Headrest, Large Laterals, w/Speakers, Black Cover ♦♦ 99756xx-85 \$385
- ♦♦Custom order; Approx. lead time: 6-8 weeks.
♦Turtle Bar System Short or Long Required; x98619-25 or x98619-26



9975xxx-xx

Headband System for Anatomic Head Support Only One Size 99500 \$230

Universal Headrest Mount ♦ One Size 87590 \$61

♦Required when using other manufactures Head Support

Allows most Whitmyer, Stealth, Otto Bock headrest, I2I to mount to our products (x:panda, Serval, etc.). Please contact Customer Service for confirmation

Other types of headrest available, please see Head Supports Order Form

Description	Measurements	Size	Part No.	USD Price
-------------	--------------	------	----------	-----------

Directional Stabilizer Kit ♦♦ One Size 89592-KIT \$215

♦♦Custom order; Approx. lead time: 6-8 weeks.

Tray (Ea.) 1,2 95125 \$300

Tray (Ea.) 3 95126 \$300



Accessories	Description	Measurements	Size	Part No.	USD Price
-------------	-------------	--------------	------	----------	-----------

Tray, Grey 1,2 8976350 \$225



Push Brace, Seat Mounted 1,2 95148-1 \$230

Push Brace, Seat Mounted 3 95148-3 \$230





Angle Adjustable Push Brace (Ea.) 3 95147 \$210




Helper's Handle, Fold Down (Pr.) One Size 95145 \$80

Helper's Handle, Fixed (Pr.) One Size 95146 \$80




Anti-Tip Device One Size 95953 \$160

Spare Parts	Description	Measurements	Size	Part No.	USD Price
-------------	-------------	--------------	------	----------	-----------



Wombat Living Step 1 to 2 Kit, Seat & Back 1 95560 \$280

Wombat Living Step 1 to 2 Kit, Seat & Back 2 95561 \$280



Wombat Living Short to High Back Kit 1 95560-1 \$157

Wombat Living Short to High Back Kit 2 95561-1 \$157

Description	Size	Part No.	USD Price
-------------	------	----------	-----------

Required, must select a color

xx = Color

Red = 24

Petroleum = 25

Grey = 28

Cover for Standard Back Cushion 1 955611-2xx \$109

Cover for Standard Back Cushion 2 955612-2xx \$109

Cover for Standard Back Cushion 3 955613-2xx \$109

Cover for Contoured Back Cushion, Grey 1 955611-238 \$137

Cover for Contoured Back Cushion, Grey 2 955612-238 \$137

Cover for Contoured Back Cushion, Grey 3 955613-238 \$137

Cover for Seat Cushions 1 955621-2xx \$109

Cover for Seat Cushions 2 955622-2xx \$109

Cover for Seat Cushions 3 955623-2xx \$109

Description	Size	Part No.	USD Price
-------------	------	----------	-----------

Cover w/Foam for Pelvis Support without Sides 1,2 95122-1 \$130

Cover w/Foam for Pelvis Support without Sides 3 95122-2 \$130

Cover w/Foam for Pelvis Support with Sides 1 95123-1 \$154

Cover w/Foam for Pelvis Support with Sides 2 95123-2 \$154

Cover w/Foam for Pelvis Support with Sides 3 95123-3 \$154

Covers for Fixed Side Supports (95121) 1,2 95121-60 \$84

Spare Parts	Description	Measurements	Size	Part No.	USD Price
-------------	-------------	--------------	------	----------	-----------



Footplate, Angle & Height Adj. 12" x 7½" 95145-1 \$115

Footplate, Angle & Height Adj. 14" x 9" 85145-2 \$115



Footplate, Angle & Height Adj. ♦ 14" x 9" 95145-2 \$115

Footplate, Angle & Height Adj. ♦ 16" x 9½" 95145-3 \$115

Footplate, Angle & Height Adj. ♦ 18¾" x 9½" 95145-4 \$115

♦ Includes 4" mounting fitting



Footplate with Rod 7½" x 12" 1 95144-1 \$225



<input type="checkbox"/> Footplate with Rod	9" x 14"	2	95144-2	\$225
<input type="checkbox"/> Footplate with Rod	9¼" x 16"	3	95144-3	\$225
<input type="checkbox"/> Footplate with Rod	9½" x 19"	4	95144-4	\$225

Description	Size	Part No.	USD Price
-------------	------	----------	-----------

<input type="checkbox"/> Batteries for Wombat Living with Power (Set)	One Size	2222012	\$305
◆ Both batteries <u>must</u> be charged at the same time			



<input type="checkbox"/> Charger & Plug for USA/Canada (Ea.)	One Size	2222161	\$320
--	----------	---------	-------

Prices in USD & is subject to change without notice