



**New** Models  
1124 and 1105

## Evolution Series Rehab Shower/Commode Chair Models 1124 and 1105

Models 1124 and 1105 are the next additions to our Evolution Series Chairs

The modular components and frame adjustments (**seat height and seat slope on 1124 and seat height on 1105**) will help to alleviate the need for custom design and its consequent added costs and delivery times.

These models will be an ideal chair for facility use as they can be used as a fitting chair in preparation for a clients discharge. The design also addresses the evolving needs of the rehab client once they return to a home environment.



Model 1124



Model 1105

### EVOLUTION CHAIR ACCESSORIES

- Option Codes
- A1 Arm Troughs
  - A4 \*Height Adjustable Flip Back Padded Arms
  - C7 \*Comfortuff™ Contoured Seat (Front Open)
  - C8 Ensolite U-Shaped Seat 2" Waterfall
  - D9 \*Swing Away Adjustable Laterals
  - D8 \*Sling Back with Stroller Handles
  - F Seat Belt (hook and loop)
  - F2 Bodypoint® Trunk Belt
  - G9 \*Height Adjustable Swing Away Footrests – Removable
  - G10 \*Height Adjustable Swing Away Footrests – Institutional/Captive
  - H Handrings w/projections
  - H2 PVC Dipped Handrings
  - H5 Pneumatic Tires
  - H7 Solid Treaded 24" Tire
  - M Heel Loops
  - N Wheel Lock Handle Extensions
  - P Pan Hanger and Pan

(\*See Price List order form for option photos)

EVOLUTION SERIES



## Evolution Series Rehab Shower/Commode Chair Models 1124/1105 Adjustable Seat Slope

	<b>1124</b>	<b>1105</b>
Seat width	18"	18"
Seat depth	18"	18"
Floor to top of seat height (front of seat)	$20 \frac{3}{4}" - 23 \frac{3}{4}"$	$20 \frac{3}{4}" - 23 \frac{3}{4}"$
Back height (from seat)	19"	19"
Arm height (from seat – standard arm)	$7 \frac{3}{4}"$	$7 \frac{3}{4}"$
Under seat clearance (Rear)	<i>Flat</i>	<i>1" Slope</i> <i>2" Slope</i>
Seat slope adjustments not available on 1105     Lowered	$15 \frac{1}{4}"$	$14 \frac{1}{4}"$ $13 \frac{1}{4}"$
Raised	$18 \frac{1}{4}"$	$17 \frac{1}{4}"$ $16 \frac{1}{4}"$
Under seat clearance with option P (floor to pan rails)	$13 \frac{3}{4}" - 16 \frac{3}{4}"$	$13 \frac{3}{4}" - 16 \frac{3}{4}"$
Footrest height adjustment with option G9 or G10	16" – 20"	16" – 20"
Overall width	$25 \frac{1}{4}"$	$24 \frac{1}{2}"$
Overall Length	$35 \frac{1}{2}"$	$30 \frac{1}{2}"$
Overall Length w/option G9/G10 (Ht. Adjustable Footrests)	$44 \frac{3}{4}"$	$41 \frac{3}{4}"$
Weight Capacity	350lbs.	350lbs.

1124 shown with option G9

