



REHAB SHOWER/COMMUNE CHAIR

- Complete Packages for Easier Ordering Available
- Stainless Steel Frame
- Several Frame and Seat Sizes Available
18"Wx18"D, 18"Wx20"D, 20"Wx18"D and
22"Wx18"D
- Three Mobility Choices:
Four 5" Casters, 20", or 24"
Wheels
- 350 lb. (158kg) Weight Capacity

All ActiveAid chairs are stainless steel for long-lasting and durable use. The large selection of seats allows for customizing for the individual. The 600 has many available options so the chair can be fitted to the user, or it could be set up to be used in a multi-use facility. From the least involved client to the most involved, this chair is readily able to fit their needs.

The 600 series shower/commode chairs have features that address the needs of the homebased users, long-term care, home care clients, and caregivers. This is the go-to chair because of all the options and sizes available. From 5" casters to large wheels and many choices to pick from, this series is designed to make those daily activities just a little bit easier.

FRAMES



P13675
18"Wx18"D
Stainless Steel
Frame



P13676
18"Wx20"D
Stainless Steel
Frame



P13677
20"Wx18"D
Stainless Steel
Frame



P13678
22"Wx18"D
Stainless Steel
Frame



P60800-2
5" Dual Locking
Casters
(pair)



P13860
20" Wheels
(pair)



P13863
20" Pneumatic
Street Tread
Wheels
(pair)

WHEELS/REAR CASTERS



P13862
24" Wheels
(pair)



P13864
24" Pneumatic
Street Tread
Wheels
(pair)



P13865
24" Flat Free Street
Tread Wheels
(pair)



P13855
20" PVC Handrims
(replaces standard
handrims, pair)



P13857
24" PVC Handrims
(replaces standard
handrims, pair)



P13649
Bolted Axles with
Wheel Mounts and
Tipping Levers
(pair)



P13650
Quick Release Axles
with Wheel Mounts
and Tipping Levers
(pair)

AXLES

WHEEL LOCKS



P60532
Wheel Locks
(pair)



P13135
Wheel Locks with
Extended Handle
(pair)



P60841-2
5" Casters
(pair)



P60800-2
5" Dual Locking
Casters
(pair)

SEATS

*Seat works in all 4 positions



P60993*
18"Wx18"D
Comfortuff 1.5" waterfall foam
Front/Rear Open Seat
A-5"xB-8" opening



P60974
18"Wx18"D
Comfortuff sloped foam
Front Open Seat
A-4"xB-6.5" opening



P60701*
18"Wx18"D
Ensolite® 1" foam
Front/Rear Open Seat
A-4.55"xB-7.75" opening



P61058
18"Wx18"D
Ensolite® 2" waterfall foam
Front Open Contoured Seat
A-5.5"xB-8.25" opening



P61086
18"Wx18"D
Ensolite® 2" foam
Solid Seat



P60675
18"Wx18"D
Ensolite® 1" foam
Cloverleaf Seat
A-7.5"xB-9" oval opening



P60913*
18"Wx18"D
Ensolite® 2" waterfall foam
Front/Rear Open Seat
A-5.25"xB-8" opening



P61121
18"Wx20"D
Ensolite® 1" foam
Front/Rear Open Seat
A-4.5"xB-7.5" opening



P61122
18"Wx20"D
Ensolite® 1" foam
Solid Seat



PP60903
18"Wx20"D
Ensolite® 1" foam
Side Open Seat
A-4.5"xB-7.75" opening



P60903
18"Wx18"D
Ensolite® 1" foam
Front/Rear Open Seat
A-4.5"xB-7.75" opening



PP61122
18"Wx18"D
Ensolite® 1" foam
Solid Seat



P61038
20"Wx18"D
Comfortuff 1.5" waterfall foam
Front Open Seat
A-5"xB-8" opening



P60905
22"Wx18"D
Ensolite® 1" foam
Front/Rear Open Seat
A-4.5"xB-7.75" opening

ARMRESTS



P13278
18"D Left Hand
Armrest with Pad



P13277
18"D Right Hand
Armrest with Pad



P13663
20"D Left Hand
Armrest with Pad

LEG RESTS



P13664
20"D Right Hand
Armrest with Pad



P13665
Left Leg Rest



P13666
Right Leg Rest



P13582
Left Elevating Leg
Rest



P13583
Right Elevating Leg
Rest



PL13581
Left Gooseneck
Leg Rest



PR13581
Right Gooseneck
Leg Rest

600 OPTIONS

FOOT PLATES



P13667
18"W 1-pc. Foot Plate/Leg Rest



P13668
20"W 1-pc. Foot Plate/Leg Rest



P13669
22"W 1-pc. Foot Plate/Leg Rest



PL11099
Left Foot Plate



PR11099
Right Foot Plate



PL11100
Left Foot Plate with Heel Loop



PR11100
Right Foot Plate with Heel Loop

BACKS



P13891
Left Oversized Aluminum Foot Plate



P13892
Right Oversized Aluminum Foot Plate



PL11101
Left Padded Foot Plate



PR11101
Right Padded Foot Plate



P13651
18"W Sling Back with Stroller Handles



P13652
18"W Solid Back



P13658
20"W Sling Back with Stroller Handles

BELTS



P13659
20"W Solid Back



P13661
22"W Sling Back with Stroller Handles



P13662
22"W Solid Back



PW60909
Waist Belt
(for 18"W frame only, 54" circumference with 12" overlap)



PW60884
Waist Belt
(for 20"W and 22"W frames only, 66" and 68" circumference with 12" overlap)



PT60884
Trunk Belt
(64", 66" and 68" circumference with 12" overlap)



P60879
Calf Strap



P11093
Bodypoint® Trunk Belt



P11102
Bodypoint® Trunk Belt (pair)



P11116
Bodypoint® Trunk Belt & Waist Belt



P11112
Bodypoint® XL Trunk Belt



P11113
Bodypoint® Belt w/Quick Release Knobs



P61081
Bodypoint® Calf Strap
(only available with 18" back)



P11114
Bodypoint® Knee Belt

POSITIONING AIDS



PL13574
Left Hand Swing-Away Adjustable Lateral



PR13574
Right Hand Swing-Away Adjustable Lateral



PL13618
Left Hand Arm-Mounted Hip Guide



PR13618
Right Hand Arm-Mounted Hip Guide



P61006
Abductor/Deflector Large



P13670
Front Anti-Tippers

PANS



P13840
18"Wx18"D Pan/Hanger



P13841
18"Wx20"D Pan/Hanger




P13842
20"Wx18"D Pan/Hanger




P13843
22"Wx18"D Pan/Hanger

Seat Opening Measurements



Cloverleaf Seat Opening Measurements





600 with Four
5" Locking Casters



600 with 20" Wheels

SPECIFICATIONS

	18"Wx18"D Frame	18"Wx20"D Frame	20"Wx18"D Frame	22"Wx18"D Frame
Weight Capacity	(158kg) 350 lbs.	350 lbs.	350 lbs.	350 lbs.
Seat Width	18" (45cm)	18" (45cm)	20" (50cm)	22" (55cm)
Width between Arm Rests	22.25" (56cm)	22.25" (56cm)	24.25" (61cm)	26.25" (66cm)
Seat Depth	18" (45cm)	20" (50cm)	18" (45cm)	18" (45cm)
Seat Height Range From Floor to Top of Seat (transfer height)	Front to Back - 24" Wheels - 23"-21.25" (58-53cm) 5"rear casters - 23"-22.25" (58-56cm)			
Leg Rest Range	with Padded Footplate			
Leg Rest	17.75"-18.75" (45-47cm)		16.5"-17.5" (41-44cm)	
Gooseneck	14"-20.5" (35-52cm)		12.75"-19.25" (32-48cm)	
Elevating - Up	17"-24" (43-60cm)		15.75"-22.75" (40-57cm)	
Elevating - Down	16.5"-23.5" (41-59cm)		15.25"-22.25" (48-56cm)	
One-Piece	19" (48cm)			
Usable Seat Depth w/Solid Back	18" (45cm)	20" (50cm)	18" (45cm)	18" (45cm)
Arm Height from Top of Seat	8" (20cm)			
Back Height From Top of Seat	8.25" (20cm)			
Gap between top of seat and bottom of back	16.5" (41cm)			
Seat to Back Angle (or Range) (Recline)	8.5" (21cm)			
Under Seat Height Clearance	11°			
Under Seat Width Clearance	18.25" (46cm)			
Overall Width	14" (35cm)		14" (35cm)	16" (40cm)
with 24" Wheels	15.75" (40cm)	15.75" (40cm)	17.75" (45cm)	19.75" (50cm)
with 5" Casters and Armrest w/Pad	25.25" (64cm)	25.25" (64cm)	27.25" (69cm)	29.25" (74cm)
with 5" Casters and Armrest with Arm Troughs	22.5" (57cm)	22.5" (57cm)	24.5" (62cm)	26.5" (67cm)
Overall Length	26.5" (67cm)	26.5" (67cm)	28.5" (72cm)	30.5" (77cm)
Min. with 24" Wheels	42.5" (107cm)	44.25" (57cm)	42.5" (107cm)	42.5" (107cm)
Min. with 5" Casters	41.5" (105cm)	43.5" (110cm)	41.5" (105cm)	41.5" (105cm)
Back of Unit to Start of Commode Opening			15.5" (39cm)	
with 24" Wheels			14.75" (37cm)	
with 5" Casters				
Turning Diameter Min. w/ 5" rear casters	42" (106cm)	44" (111cm)	44" (111cm)	45" (114cm)
Seat Used to Generate Specifications	#60905 22"Wx18"Dx2"H			