



Axiom Phoenix

Standard Features

- Light - Weight aluminum frame
- Aluminum footrest w/adjustable straps
- Rear wheels - 16" quick release mags
- Front wheel - 16" quick release mag
- Hand brake for slowing and stopping
- Rear footbrake for Captain transfer/parking
- 5 - point adjustable harness
- Padded seat w\recline
- Zippered storage pocket behind seat
- Behind seat storage net
- Under seat basket

FRAME Specs

assembled size w/canopy open	assemble size w/out canopy	folded w/wheels	folded w/out wheel	QR wheels	wheel size	wheel type	hand brake	rear brake	footrest	folding steps	chair weight	weight limit
70" x 26" x 50"	70" x 26" x 41"	54" x 26" x 20"	46" x 22" x 13"	Yes	16x16	mags	Yes	Yes	aluminum	2	30 lbs	200

SEATING specs

seatback height	seat depth	seat width	recline angle	5 point harness	storage basket	storage pocket	storage cargo net	sun canopy	wrist strap
32"	16"	18"	140°	Yes	Yes	Yes	Yes	Yes	Yes

