



# Axiom Improv Order Form

AXIOM IMPROV  
PUSH CHAIR

AXIOM IMPROV PUSH CHAIR  
FOLDED W/OUT WHEELS



Dealer Name: \_\_\_\_\_  
 Company Name: \_\_\_\_\_  
 Phone: \_\_\_\_\_  
 E-mail: \_\_\_\_\_  
 Ship to Address \_\_\_\_\_  
 Reference Name \_\_\_\_\_

Color :                      Size:

Red                       2

Navy                       3

**Specifications :**

Size	2	3
Weight Limit	150 LBS	200 LBS
Seatback Height	28"	31"
Seat Depth	15"	15"
Seat Width	15"	18"
Lower Leg Length	17"	17"
Assembled Size	63.5x24x46	63.5x27x45
Folded w/ Wheels	51.5x24x19	51.5x27x19
Folded w/o Wheels	43.5x19x12.25	43.5x22x12.25
Unit Weight	32.8 LBS	35.2 LBS

**Accessories :**

- |   |   |
|---|---|
| <input type="checkbox"/> Head Rest                | <input type="checkbox"/> Anti-Tippers           |
| <input type="checkbox"/> Seat Abductor            | <input type="checkbox"/> Bug Canopy             |
| <input type="checkbox"/> Hip Abductor - Single    | <input type="checkbox"/> Seat Bottom Insert Pad |
| <input type="checkbox"/> Hip Abductor - Double    | <input type="checkbox"/> Seat Back Insert Pad   |
| <input type="checkbox"/> 2pc Lateral Support      | <input type="checkbox"/> Swivel Wheel Axis Kit  |
| <input type="checkbox"/> Positioning Wedges (Set) | <input type="checkbox"/> Travel Bag             |
| <input type="checkbox"/> Harness Extension        | <input type="checkbox"/> Rain Canopy            |
| <input type="checkbox"/> Sheerling Insert         | <input type="checkbox"/> Fender                 |
| <input type="checkbox"/> Flashing Light           |   |

**Customizations:**

- Seatback Height Increase to \_\_\_\_\_ " (max height 35")
- Beach Wheels
- Seat Depth Increase to \_\_\_\_\_ " (max height 18")
- 8" Swivel Wheel Axis Kit

Always require a measurement chart.  
 ALL CUSTOMS ARE FINAL SALE

**Included with Chair**

*Seat	*16" Single Fixed Front Wheel
*5 Point Harness	*Two 16" Quick Release Rear Wheels
*Sun Canopy	*Hand Brake
*Basket underneath chair	*Rear Footbrake
*Zipper Pouch	*2 Water bottle Pocket Holder
*Adjustable Push Handle	*Saftey Wrist Strap

**\*Multi positional patient system with integrated seat, operated by the caregiver.**